


Jun 08, 2020

🌐 Targeted proteomic LC-MS/MS analysis

DOI

dx.doi.org/10.17504/protocols.io.bf9xjr7n

LC-MS/MS



Yan Chen¹, Jennifer Gin¹, Christopher J Petzold¹

¹Lawrence Berkeley National Laboratory

LBNL omics



Yan Chen

LBNL

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bf9xjr7n>

External link: <https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0211582>

Protocol Citation: Yan Chen, Jennifer Gin, Christopher J Petzold 2020. Targeted proteomic LC-MS/MS analysis. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bf9xjr7n>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 12, 2020

Last Modified: June 08, 2020

Protocol Integer ID: 36887

Keywords: LC-MS, Targeted proteomics, MRM, targeted proteomic data acquisition, targeted proteomic lc, throughput quantitative proteomic application, quantitative proteomic application, proteomic data, ms analysis this protocol details step, ms analysis this protocol, ms analysis, m, flow uhplc

Abstract

This protocol details steps in targeted proteomic data acquisition with a standard-flow UHPLC-QQQ system. It was adapted from Chen, Y. et al. "A rapid methods development workflow for high-throughput quantitative proteomic applications." PLoS ONE 14,2 e0211582. 14 Feb. 2019, doi:10.1371/journal.pone.0211582.

Materials

MATERIALS


 Acetonitrile LCMS quality **JT Baker Catalog #9829-02**

 Pierce & Warriner; Formic Acid **Thermo Fisher Catalog #TS-28905**

 LCMS grade water **VWR International (Avantor) Catalog #BJLC365-2.5**

 Isopropanol **VWR International (Avantor) Catalog #BJ650447-4L**

STEP MATERIALS

 2.7 μm

Analytical column: Ascentis Express Peptide C18 column ( 2.7 μm particle size, 160-Å pore size, 5-cm length ×  2.1 mm internal diameter).

Guard column: Ascentis guard column (5-mm ×  2.1 mm ID with  2.7 μm particle size, 160-Å pore size)

LC-MS system: Agilent 6460 QQQ mass spectrometer system coupled with an Agilent 1290 Infinity II UHPLC system (Agilent Technologies, Santa Clara, CA)

Protocol materials

 Acetonitrile LCMS quality **JT Baker Catalog #9829-02**

 Pierce™ Formic Acid **Thermo Fisher Catalog #TS-28905**


 LCMS grade water **VWR International (Avantor) Catalog #BJLC365-2.5**

 Isopropanol **VWR International (Avantor) Catalog #BJ650447-4L**

 Agilent Tune Mix: G1969-85000 **Catalog #G1969-85000**

Troubleshooting

Safety warnings

-  Wear proper PPE (gloves, safety goggle, and lab coat), and prepare solvents in a chemical hood. Store organic solvents in a flammable storage cabinet.

Before start



Prepare the following solvents:

Needle wash solvents: Add  100 mL isopropanol into  900 mL water.

Solvent A: Add  0.1 % volume formic acid into LC-MS grade water.

Solvent B: Add  0.1 % volume formic acid into LC-MS grade acetonitrile

Proteomics: HPLC and Mass Spectrometry

- 1 Thaw peptide samples  On ice , and transfer  30 µL of each sample to LC autosampler vials (Agilent, Cat.# 5182-0567, 5182-0564) or 96-well plate (Bio-Rad, Cat.# HSP9655).
- 2 Liquid chromatography tandem-mass spectrometry (LC-MS/MS) analysis is performed with an Agilent 6460 QQQ mass spectrometer system coupled to an Agilent 1290 Infinity II UHPLC system (Agilent Technologies, Santa Clara, CA).

Equipment

6460QQQ

NAME

Triple quadrupole mass spectrometer

TYPE

Agilent Technologies

BRAND

6460QQQ

SKU

<https://www.agilent.com/en/products/liquid-chromatography-mass-spectrometry-lc-ms/lc-ms-instruments/triple-quadrupole-lc-ms>

LINK

Equipment

1290 UHPLC

NAME

Ultra-high performance liquid chromatography system

TYPE

Agilent Technologies

BRAND

1290UHPLC

SKU

www.agilent.com

LINK


3 Samples were loaded on the temperature controlled autosampler at 4 °C . The Agilent 1290 Infinity II UHPLC is equipped with a Sigma–Aldrich Ascentis Peptides ES-C18 analytical column (2.1 mm ID , 50 mm length , 2.7 µm particle size , and 160-Å pore size (Sigma-Aldrich, Cat.# 53301-U)) coupled to a 2.1 mm ID , 5 mm length guard column (Sigma-Aldrich, Cat.# 53536-U) with the same particle and pore size. The column is operated at 60 °C .

4 Twenty micrograms 20 µg of peptides are loaded onto the column from each sample and separated using a linear mobile-phase gradient consisting of 0.1% formic acid in water (Solvent A) and 0.1% formic acid in acetonitrile (Solvent B) operating at a flow rate of 0.4 ml/min. A 10 minutes gradient of chromatographic separation is as follows:

	Step	%A	%B	Time (min)
	1	98	2	0.00
	2	98	2	0.50
	3	65	35	5.50
	4	20	80	6.25
	5	20	80	8.00
	6	98	2	8.50
	7	98	2	10.00

Chromatographic gradient table

%B vs. Time



Chromatographic gradient

Note

The gradient length depends on the application of interest and the number of target peptides for detection.

- 5 The eluted peptides were ionized via an Agilent Jet Stream ESI source operating in positive ion mode with the following source parameters:

Gas temperature	250 °C
Gas flow	13 liters/min
Nebulizer pressure	35 psi
Sheath gas temperature	250 °C
Sheath gas flow	11 liters/min
Capillary voltage	3500 V
Nozzle voltage	0 V

6 SRM transitions were monitored through either scanning or scheduled MRM mode with the following criteria:

- Minimum dwell time of a transition is 2 ms
- Maximum number of concurrent transitions are 200
- Total cycle time less than 1.5 sec.

7 The MS raw data were acquired using Agilent MassHunter version B.08.02

Software

MassHunter

NAME

Agilent Technologies, Inc

DEVELOPER

8 Acquired SRM data were imported and analyzed by Skyline software version 20.1 (MacCoss Lab Software).

Software

Skyline

NAME

University of Washington

DEVELOPER

9 For a small sample set, peak integration is visually checked and validated by browsing replicate peak area and retention time tabs. For a large set of more than 100 samples, peak integration is completed by using skyline advanced peak picking models (e.g., mProphet).

Citation

Reiter L, Rinner O, Picotti P, Hüttenhain R, Beck M, Brusniak MY, Hengartner MO, Aebersold R (2011)

. mProphet: automated data processing and statistical validation for large-scale SRM experiments..
Nature methods.

<https://doi.org/10.1038/nmeth.1584>

LINK

- 10 A peptide quantitative report was exported to a .csv file for further analysis and sharing. The skyline file with the data and SRM methods information was uploaded to a project folder on the PanoramaWeb server. The project information (data and SRM methods) is made public when a manuscript is accepted for journal publication.

Citation

Sharma V, Eckels J, Taylor GK, Shulman NJ, Stergachis AB, Joyner SA, Yan P, Whiteaker JR, Halusa GN, Schilling B, Gibson BW, Colangelo CM, Paulovich AG, Carr SA, Jaffe JD, MacCoss MJ, MacLean B

(2014). Panorama: a targeted proteomics knowledge base..
Journal of proteome research.

<https://doi.org/10.1021/pr5006636>

LINK

Citation

Sharma V, Eckels J, Schilling B, Ludwig C, Jaffe JD, MacCoss MJ, MacLean B (2018)

. Panorama Public: A Public Repository for Quantitative Data Sets Processed in Skyline..
Molecular & cellular proteomics : MCP.

<https://doi.org/10.1074/mcp.RA117.000543>


LINK

UHPLC-QQQ performance monitoring and QC

- 11 The mass spectrometer is subjected to a full autotune at least quarterly (and more frequently, if necessary) to optimize ion transmission and update EM voltage.

 Agilent Tune Mix: G1969-85000 **Catalog #G1969-85000**

- 12 Agilent 6460 QQQ mass spectrometer is subjected to a check tune prior to analyzing a batch of samples for resolution and mass accuracy verification.
- 13 The UHPLC-QQQ performance is tracked **daily** by injecting 50 fmol BSA tryptic digest standard and monitoring eight BSA peptides using a scheduled MRM assay. Peptide retention time, peak shape, and peak intensity of these peptides are monitored to determine if maintenance is required.
- 14 The performance tracking is automated and monitored via PanoramaWeb server through an established AutoQC pipeline.



An example of autoQC instrument performance track on the peak area of eight BSA tryptic peptide measurements



Software

AutoQC loader

NAME

University of Washington

DEVELOPER

Citation

Bereman MS, Beri J, Sharma V, Nathe C, Eckels J, MacLean B, MacCoss MJ (2016)
. An Automated Pipeline to Monitor System Performance in Liquid Chromatography-
Tandem Mass Spectrometry Proteomic Experiments.
Journal of proteome research.

<https://doi.org/10.1021/acs.jproteome.6b00744>

LINK



Citations

Step 10

Sharma V, Eckels J, Taylor GK, Shulman NJ, Stergachis AB, Joyner SA, Yan P, Whiteaker JR, Halusa GN, Schilling B, Gibson BW, Colangelo CM, Paulovich AG, Carr SA, Jaffe JD, MacCoss MJ, MacLean B. Panorama: a targeted proteomics knowledge base.

<https://doi.org/10.1021/pr5006636>

Step 10

Sharma V, Eckels J, Schilling B, Ludwig C, Jaffe JD, MacCoss MJ, MacLean B. Panorama Public: A Public Repository for Quantitative Data Sets Processed in Skyline.

<https://doi.org/10.1074/mcp.RA117.000543>

Step 14

Bereman MS, Beri J, Sharma V, Nathe C, Eckels J, MacLean B, MacCoss MJ. An Automated Pipeline to Monitor System Performance in Liquid Chromatography-Tandem Mass Spectrometry Proteomic Experiments

<https://doi.org/10.1021/acs.jproteome.6b00744>

Step 9

Reiter L, Rinner O, Picotti P, Hüttenhain R, Beck M, Brusniak MY, Hengartner MO, Aebersold R. mProphet: automated data processing and statistical validation for large-scale SRM experiments.

<https://doi.org/10.1038/nmeth.1584>