

Feb 01, 2024

Tau aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader

DOI

dx.doi.org/10.17504/protocols.io.dm6gp3nw8vzp/v1

Patricia Yuste-Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3nw8vzp/v1>

Protocol Citation: Patricia Yuste-Checa, F Ulrich Hartl 2024. Tau aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp3nw8vzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 25, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 94537

Keywords: ASAPCRN, tau aggregation, plate reader, fluorescence, thioflavin

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to efficiently monitor Tau aggregation by thioflavin T fluorescence in a plate reader.

Attachments








Tau agg monitored by...


1.8MB

Materials

Buffers:


- Heparin:  50 μL in water, Heparin sodium salt from porcine intestinal mucosa (Merck, H3393). Can be stored at  4 $^{\circ}\text{C}$. Prepare a fresh working dilution at  4.5 μL in water ([M] 250 micromolar (μM)).
- [M] 0.2 Molarity (M) MgCl_2
- [M] 1 millimolar (mM) Thioflavin T (ThT). Can be stored at  -20 $^{\circ}\text{C}$.

 Heparin sodium salt from porcine intestinal mucosa **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H3393**

 96 well half-area plate of black polystyrene with a clear bottom **Corning Catalog #3881**

Troubleshooting



Tau aggregation monitored by thioflavin-T (ThT) fluorescence in a plate reader

- 1 Mix reagents to a final concentration of [M] 10 micromolar (μM) Tau (TauRD), [M] 2.5 micromolar (μM) Heparin, [M] 2 millimolar (mM) MgCl_2 , [M] 10 micromolar (μM) ThT in 1x PBS pH 7.2 . Prepare a mix for 4.5 reactions per condition (per condition, four technical replicates in each plate). Final volume is  80 μL per well.



Note

Molecular or chemical chaperones can be included in the reaction in order to study their effect on Tau aggregation.

- 2 Dispense  80 μL of the mix into a well of a 96 well half-area plate of black polystyrene with a clear bottom (Corning 3881).
- 3 If possible, the outer wells should not be used and instead filled with water.
- 4 Seal the plate with parafilm to avoid evaporation.
- 5 Set the following parameters in a SPARK multimode microplate reader (TECAN) and start the reaction:
 - Fluorescence measurement: ThT signal, excitation 440 nm, emission 480 nm, (use gain regulation) measured every 2 minutes.
 - Temperature  37 $^{\circ}\text{C}$
 - Constant shaking (50 seconds linear shaking: amplitude 4.5 mm, frequency 420 rpm - 50 seconds orbital shaking: amplitude 1.5 mm, frequency 360 rpm).

Note

- SPARK multimode microplate reader (TECAN) is highly recommended due to its high sensitivity and the gain regulation mode that increases the fluorescence detection window. When using other plate readers, the sample ThT signal easily gets saturated even when reducing initial gain to the minimum gain capacity.
- Cysteine-free TauRD (Tau residues 244-371, C291A/P301L/C322A/V337M), rapidly aggregates under these conditions, reaching the ThT plateau after 2 h aggregation.

