

Jan 10, 2019

Tau Thioflavin T Assay

DOI

<https://dx.doi.org/10.17504/protocols.io.wygffjw>

Alexandra Netter-Glangeaud¹

¹StressMarq Biosciences Inc

StressMarq



Patricia Thomson

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.wygffjw>

External link: <https://www.stressmarq.com/support/technical-support/protocols/tau-protocol/>

Protocol Citation: Alexandra Netter-Glangeaud 2019. Tau Thioflavin T Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.wygffjw>

Manuscript citation:



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 10, 2019

Last Modified: January 10, 2019

Protocol Integer ID: 19144

Keywords: Tau, Thioflavin T, Fibrils, PFFs, assay

Materials

MATERIALS

⊗ Thioflavin T **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T3516**

⊗ Active Human Recombinant Tau441 (2N4R) P301S mutant Protein Monomer **StressMarq Biosciences Catalog #SPR-327**

⊗ Active Human Recombinant Tau (K18) P301L mutant Protein Monomer **StressMarq Biosciences Catalog #SPR-328**

⊗ Active Human Recombinant Tau441 (2N4R) P301S mutant Protein Pre-formed Fibrils **StressMarq Biosciences Catalog #SPR-329**

⊗ Active Human Recombinant Tau (K18) P301L mutant Protein Pre-formed Fibrils **StressMarq Biosciences Catalog #SPR-330**

⊗ Microplate 96 well PS F-bottom (chimney well) black non-binding **greiner bio-one Catalog #655900**

- 1 A 1mM stock solution of Thioflavin T was prepared in dH₂O (prepared fresh and filtered through a 0.2 µm syringe filter).
- 2 The thioflavin T was diluted in PBS pH 7.4 so that the final Thioflavin T concentration in each well was 25 µM (volume per well = 100 µL).
- 3 10 uM Heparin was added to each well.
- 4 Tau aliquots were thawed at room temperature just before use.
- 5 Either fibril or monomer (or both) was added to the appropriate wells. The well contents were pipetted up and down to mix.
- 6 The plate was sealed and placed in a shaking incubator (800 rpm) at  37 °C
- 7 Fluorescence was measured on a Molecular Devices Gemini XPS Microplate Reader using Softmax Pro software version 6.5.1.

XPS Microplate Reader Settings:**Temperature:** 37°C**Read Type:** Well Scan**Wavelength:** Excitation at 450 nm and Emission at 485 nm**PMT Gain:** Automatic**Flashes per read:** 6**Shake:** 20 seconds before read

- 8 The plate was re-sealed and placed into the shaking incubator at 37°C.
- 9 Readings were taking at regular intervals from 1 hour to 72 hours.