

Jul 16, 2018

tCyCIF protcol (Leica BOND) version 0.1

DOI

<https://dx.doi.org/10.17504/protocols.io.rpxd5pn>

Jia-Ren Lin¹, Benjamin Izar¹, Shaolin Mei¹, Shu Wang¹, Peter Sorger¹

¹Laboratory of Systems Pharmacology, Harvard Medical School

Human Cell Atlas Metho...

Laboratory of Systems P...

2 more workspaces



Jia Ren Lin

Harvard Medical School, Laboratory of Systems Pharmacology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rpxd5pn>

External link: www.cycif.org

Protocol Citation: Jia-Ren Lin, Benjamin Izar, Shaolin Mei, Shu Wang, Peter Sorger 2018. tCyCIF protcol (Leica BOND) version 0.1. protocols.io <https://dx.doi.org/10.17504/protocols.io.rpxd5pn>

**Manuscript citation:**CITE AS: eLife 2018;7:e31657 DOI: [10.7554/eLife.31657](https://doi.org/10.7554/eLife.31657)

Lin, Jia-Ren, et al. "Highly multiplexed immunofluorescence imaging of human tissues and tumors using t-CyCIF and conventional optical microscopes." *bioRxiv* (2018): 151738.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 16, 2018

Last Modified: July 16, 2018

Protocol Integer ID: 13783

Keywords: Multiplexed Imaging, Systems biology, FFPE, cyclic immunofluorescence, cycif, based cyclic immunofluorescence, fluorescence imaging, fluorescence imaging of formalin, immune markers in diverse tissue, plex fluorescence image, tcycif protocol, tumor antigen, cell genomic, immune marker, resolution imaging, diverse tissue, cancer

Abstract

The architecture of normal and diseased tissues strongly influences the development and progression of disease as well as responsiveness and resistance to therapy. We describe a tissue-based cyclic immunofluorescence (t-CyCIF) method for highly multiplexed immuno-fluorescence imaging of formalin-fixed, paraffin-embedded (FFPE) specimens mounted on glass slides, the most widely used specimens for histopathological diagnosis of cancer and other diseases. t-CyCIF generates up to 60-plex images using an iterative process (a cycle) in which conventional low-plex fluorescence images are repeatedly collected from the same sample and then assembled into a high dimensional representation. t-CyCIF requires no specialized instruments or reagents and is compatible with super-resolution imaging; we demonstrate its application to quantifying signal transduction cascades, tumor antigens and immune markers in diverse tissues and tumors. The simplicity and adaptability of t-CyCIF makes it an effective method for pre-clinical and clinical research and a natural complement to single-cell genomics.

Guidelines

Key resources, reagents and software used in this study are listed on-line at the HMS LINCS Center Publication Page <http://lincs.hms.harvard.edu/lin-elife-2018/> (RRID:SCR_016370). This page provides links to an OMERO image database from which individual images can be obtained; stitched and registered image panels can be obtained at www.cycif.org (RRID:SCR_016267) and a video illustrating the t-CyCIF method can be found at <https://vimeo.com/269885646>.



Materials

MATERIALS

 Hydrogen peroxide solution 50 wt. % in H₂O, stabilized **Merck MilliporeSigma (Sigma-Aldrich) Catalog #516813-500ML**

 1N NaOH **Merck MilliporeSigma (Sigma-Aldrich) Catalog #1091371000**

 Phosphate Buffered Saline: 5 L of 20X **Santa Cruz Biotechnology Catalog #sc-362299**

 Odyssey® Blocking Buffer (PBS) **LI-COR Catalog #927-40000**

 Hoechst 33342, Trihydrochloride, Trihydrate - 10 mg/mL Solution in Water **Invitrogen - Thermo Fisher Catalog #H3570**

Pre-processing and pre-staining tissues for t-CyCIF (Leica BOND RX)

- 1 Pre-processing of FFPE tissue and tumor slices mounted on slides was performed on a Leica BOND RX automated stained using the following protocol:

Step	Reagent	Supplier	Incubation (min)	Temp. (°C)
1	*No Reagent	N/D	30	60
2	BOND Dewax Solution	Leica	0	60
3	BOND Dewax Solution	Leica	0	R.T.
4	BOND Dewax Solution	Leica	0	R.T.
5	200 proof ethanol	User*	0	R.T.
6	200 proof ethanol	User*	0	R.T.
7	200 proof ethanol	User*	0	R.T.
8	Bond Wash Solution	Leica	0	R.T.
9	Bond Wash Solution	Leica	0	R.T.
10	Bond Wash Solution	Leica	0	R.T.
11	Bond ER1 solution	Leica	0	99
12	Bond ER1 solution	Leica	0	99
13	Bond ER1 solution	Leica	20	99
14	Bond ER1 solution	Leica	0	R.T.
15	Bond Wash Solution	Leica	0	R.T.
16	Bond Wash Solution	Leica	0	R.T.
17	Bond Wash Solution	Leica	0	R.T.
18	Bond Wash Solution	Leica	0	R.T.
19	Bond Wash Solution	Leica	0	R.T.
20	IF Block	User*	30	R.T.
21	Antibody Mix	User*	60	R.T.
22	Bond Wash Solution	Leica	0	R.T.
23	Bond Wash Solution	Leica	0	R.T.
24	Bond Wash Solution	Leica	0	R.T.

25	Hoechst Solution	User*	30	R.T.
26	Bond Wash Solution	Leica	0	R.T.
27	Bond Wash Solution	Leica	0	R.T.
28	Bond Wash Solution	Leica	0	R.T.

Steps 2-10: Dewaxing & Rehydration with Leica Bond Dewax Solution Cat. AR9222.

Steps 11-14: Antigen retrieval with BOND Epitope Retrieval solution 1 (ER1; Cat. AR9961)

Steps 15-19: Washing with Leica Bond Wash Solution (Cat. AR9590).

Steps 20-28 Pre-staining procedures as shown in Figure 1A:

Step 20: IF Block - Immunofluorescence blocking in Odyssey blocking buffer (LI-COR, Cat. 927401)

Step 21: Antibody Mix - Incubation with secondary antibodies diluted in Odyssey blocking buffer

Step 25: Staining with Hoechst 33342 at 2 µg/ml (w/v) in in Odyssey blocking buffer

Equipment

new equipment

NAME

Leica BOND RX

BRAND

3342171

SKU

<https://www.leicabiosystems.com/ihc-ish-fish/ihc-ish-instruments/products/leica-bond-rx/>

SPECIFICATIONS



Prestaining (manual)

- 2 Dewaxed specimens were blocked by incubation with Odyssey blocking buffer for 30 mins by applying the buffer to slides as a 250-500 µl droplet at room temperature; evaporation was minimized by using a slide moisture chamber (Scientific Device Laboratory, 197-BL). Slides were then pre-stained by incubation with diluted secondary antibodies (listed above) for 60 minutes, followed by washing 3 times with 1xPBS. Finally, slides were incubated with Hoechst 33342 (2 µg/ml) in 250-500 µl Odyssey blocking

buffer for 30 min in a moisture chamber and washed 3 times with 1xPBS in vertical staining jars. After imaging, cells were subjected to a round of fluorophore inactivation (see below). Following fluorophore inactivation, slides were washed 4 times with 1x PBS by dipping them in a series of vertical staining jars to remove residual inactivation solution.

Performing cyclic immunofluorescence

- 3 All primary antibodies (fluorophore-conjugated and unconjugated) were diluted in Odyssey blocking buffer. Slides carrying tissues that had been subjected to pre-staining, or to a previous t-CyCIF stain and bleach cycle, were incubated at 4°C for ~12 hr with diluted primary or fluorophore-conjugated antibody (250-500 μ l per slide) in a moisture chamber. Long incubation times were a matter of convenience and many antibodies only require short incubation with sample. Slides were then washed four times in 1x PBS by dipping in a series of vertical staining jars.
For indirect immunofluorescence, slides were incubated in diluted secondary antibodies in a moisture chamber for 1 hr at room temperature followed by four washes with 1xPBS. Slides were incubated in Hoechst 33342 at 2 μ g/ml in Odyssey blocking buffer for 15 min at room temperature, followed by four washes in 1xPBS. Stained slides were mounted prior to image acquisition (see the Mounting section below).

Mounting & de-coverslipping

- 4 Immediately prior to imaging, slides were mounted with 1xPBS or, if imaging was expected to take longer than 30 minutes, for example, in the case of samples larger than 2-4 cm² (corresponding to about 200 fields of view with a 10X objective) PBS was supplement with 10% Glycerol. Slides were covered using 24 × 60mm No. 1 coverslips (VWR 48393-106) to prevent evaporation while facilitating subsequent de-coverslipping via gravity. Following image acquisition, slides were placed in a vertical staining jar containing 1xPBS for at least 15 min. Coverslips were released from slides (and the tissue sample) via gravity as the slides were slowly drawn out of the staining jar

Fluorophore inactivation (bleaching)

- 5 After imaging, fluorophores were inactivated by placing slides horizontally in 4.5% H₂O₂ and 24 mM NaOH made up in PBS for 1 hour at RT in the presence of white light. Following fluorophore inactivation, slides were washed 4 times with 1x PBS by dipping them in a series of vertical staining jars to remove residual inactivation solution.

Image acquisition (RareCyte CyteFinder)

- 6 Stained slides from each round of CyCIF were imaged with a CyteFinder slide scanning fluorescence microscope (RareCyte Inc. Seattle WA) using either a 10X (NA=0.3) or 40X

long-working distance objective (NA = 0.6). Imager5 software (RareCyte Inc.) was used to sequentially scan the region of interest in 4 fluorescence channels. These channels are referred to by the manufacturer as a: (i) "DAPI channel" with an excitation filter having a peak of 390 nm and half-width of 18nm and an emission filter with a peak of 435nm and half-width of 48nm; (ii) "FITC channel" having a 475nm/28nm excitation filter and 525nm/48nm emission filter (iii); "Cy3 channel" having a 542nm/27nm excitation filter and 597nm/45nm emission filter and (iv); "Cy5 channel" having a 632nm/22nm excitation filter and 679nm/34nm emission filter. Imaging was performed with 2×2 binning to increase sensitivity, reduce exposure time and reduce photo bleaching.

Image processing (Software requirements)

- 7 Code and scripts used in this study are listed in Table 3 and also on-line at the HMS LINCS Center publication page (<http://lincs.hms.harvard.edu/lin-elif-2018/>). ImageJ is available at <https://imagej.nih.gov/ij/> BaSic is available at <https://www.helmholtz-muenchen.de/icb/research/groups/quantitative-single-cell-dynamics/software/basic/index.html>. Matlab scripts used in this paper and ASHLAR registration/stitching algorithm is available at our GitHub repositories (<https://github.com/sorgerlab/cycif> and <https://github.com/sorgerlab/ashlar>).

Software	
ImageJ/Fiji	NAME
Windows 7	OS
National Institutes of Health	DEVELOPER
http://wsr.imagej.net/distros/win/ij152-win-java8.zip	REPOSITORY
https://imagej.nih.gov/ij/	SOURCE LINK

Image processing (Background subtraction and image registration)

- 8 Background subtraction was performed using the previously established rolling ball algorithm (with a 50-pixel radius) in ImageJ. Adjacent background-subtracted images from the same sample were then registered to each using an ImageJ script as described previously²⁰. All images with 2×2 binning in acquisition were partially de-convoluted with unsharp masking. DAPI images from each cycle were used to generate reference

coordinates by Rigid-body transformation. To generate virtual hyper-stacked images, the transformed coordinates were applied to images from four channel imaging of each t-CyCIF cycle.

Software

ImageJ/Fiji

NAME

Windows 7

OS

National Institutes of Health

DEVELOPER

<http://wsr.imagej.net/distros/win/ij152-win-java8.zip>

REPOSITORY

<https://imagej.nih.gov/ij/>

SOURCE LINK

Image processing (Image Stitching, Shading and flat-field correlation)

- 9 The BaSiC algorithm²³ was used for shade and flat-field correction in the create of the multi-panel montage images shown in Figures 2B, 6B, 9A and 11A. Additional information can be found on the BaSiC website (<https://www.helmholtz-muenchen.de/icb/research/groups/quantitative-single-cell-dynamics/software/basic/index.html>). An example of the performance of BaSiC is shown in Figure 2-figure supplement 1. The ImageJ plugin of BaSiC was applied for whole image stacks using the default options. After processing with BaSiC, images stack were stitched with ImageJ/Fiji "Grid stitch" plugin with default options. ASHLAR was used to stich, register and scale images available at <http://www.cycif.org/>.

Image processing (Single-cell segmentation & quantification)

- 10 To obtain intensity values for single cells, images were segmented using a previously described²⁰ Watershed algorithm based on nuclear staining by Hoechst 33342. Images were initially thresholded using the OTSU algorithm and binarized in the Hoechst channel, which was then used to generate a nuclear mask image. The mask images were then subjected to the Watershed algorithm in ImageJ to obtain single-cell regions of interest (ROIs). From the nuclei, the cytoplasm was captured by centripetal expansion of either of 3 pixels in images obtained with a 10X objective or of 6 pixels in images obtained with a 40X objective, until cell reaching the cell boundaries (cell membrane). The cytoplasm was then defined as the region between the cell membrane and the nucleus. Following cell segmentation, these cell boundaries were used to compute mean and integrated intensity values from all channels. Because ROIs are (initially) defined



only by the nuclear signal, this approach is likely to over- or under- segment cells with irregular shapes, which can lead to nuclear, cytosolic or cell membrane “signal contamination” between neighboring and/or stacked cells. Further experimental (e.g. including membrane markers to guide whole-cell rather than nuclear-only segmentation) and analytical algorithms to more accurately segment individual cells (e.g. using deep learning methods to register and apply additional features) would help to improve segmentation. All imageJ scripts used in this manuscript can be found in our Github repository (<https://github.com/sorgerlab/cycif>).

Software

ImageJ/Fiji	NAME
Windows 7	OS
National Institutes of Health	DEVELOPER
http://wsr.imagej.net/distros/win/ij152-win-java8.zip	REPOSITORY
https://imagej.nih.gov/ij/	SOURCE LINK