

May 20, 2025

TEBV fabrication and perfusion

 [PLOS One](#)

 Peer-reviewed method

DOI

dx.doi.org/10.17504/protocols.io.3byl4wpwzv05/v1

Jingyi Zhu¹, Halie Hotchkiss¹, Kevin Shores¹, Stacey Maskarinec¹, George A. Truskey¹

¹Duke University

PLOS ONE Lab Protocols
Tech. support email: plosone@plos.org



Jingyi Zhu

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4wpwzv05/v1>

Protocol Citation: Jingyi Zhu, Halie Hotchkiss, Kevin Shores, Stacey Maskarinec, George A. Truskey 2025. TEBV fabrication and perfusion. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4wpwzv05/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 04, 2025

Last Modified: May 20, 2025

Protocol Integer ID: 123465

Keywords: Fabrication, THP-1 monocytes perfusion , Embedding, Fixation, Immunofluorescence staining , Engineered vessels, Microphysiological model systems, tebv fabrication, perfusion this protocol detail, perfusion, fabrication

Funders Acknowledgements:

American Heart Association

Grant ID: 388000147

American Heart Association

Grant ID: 957133

National Institutes of Health

Grant ID: HL159275

National Institutes of Health

Grant ID: HL138252

Abstract

This protocol details the TEBV fabrication and perfusion.

Attachments



[S2 Protocol.docx](#)

1.4MB

Guidelines

Troubleshooting:

Troubleshooting advice can be found in Table 2.

A	B	C	D
Step	Problem	Possible reason	Solution
Procedure 1 Step 12	Collagen solution does not gel after 60 minutes incubation. Air pockets form after 60 minutes incubation.	Improper pH of the collagen solution. The chamber is not sealed tightly. Insufficient mixing before injection.	Adjust the pH of the mixture. Change the O-ring seal the chamber tightly. Keep injecting while the collagen solution leaks into the space between two adjacent vessels. Briefly pipette to mix thoroughly before injection.
Procedure 2 Step 14	Micro bubbles observed in vessel wall after dehydration.	Aggressive pipetting during mixing.	Mix all the components gently in an Eppendorf tube. Briefly centrifuge before injecting collagen solution into the molds.
Procedure 2 Step 17	Perfusion chamber cannot be tightened evenly.	The PDMS clamps are too thick.	Trim the PDMS clamps. Change O-ring.
Procedure 2 Step 20	Additional vessels begin to leak following the repair process.	The extraction of media causes a meniscus to form, resulting in increased vessel leakage.	Disassemble the chamber while keeping the media in place and use absorbent Kimwipes to contain the liquid.
Procedure 3 Step 28	Bubbles inside lumen after rotation.	Medium leaks out or evaporates during rotation.	Fix the chamber orientation so that the mandrels are parallel to the rotator platform. Fabricate PDMS caps to seal the mandrels after the injection of the endothelial cell suspension.

	A	B	C	D
	Procedure 4 Step 33.6	Bubbles keep growing in the TEBV perfusion chamber during perfusion.	The chamber does not seal tightly. Bubbles in the side loop tubing. Possible leakage in vessels.	Seal the chamber tightly and evenly. Remove bubbles as much as possible in the chamber and side loop. Perform leakage test and replace leaky vessels.

Timing:

The following timeline assumes fabrication of one chamber of TEBVs. PDMS and tubing materials required for TEBV fabrication and perfusion can be prepared in advance.

- Procedure 0: PDMS clamps fabrication and sterilize TEBV chamber components 2-3 days
 - Procedure 1: Harvest hNDFs and prepare collagen mixture ~1 hour
 - Procedure 2: Plastic compression and dehydration ~30 minutes
 - Procedure 3: Endothelial cell seeding and set up rotation ~45 minutes
 - Procedure 4: Perfusion set up ~ 3.5 hours
 - Procedure 5: Daily media change 5-10 minutes
 - Procedure 6: Downstream experiments: Activation of endothelium and THP-1 monocytes perfusion
- TEBV fixation, *en face* staining of TEBVs

Note

Procedures 1 (Harvesting hNDFs and preparing the collagen mixture) and 2 (Plastic compression and dehydration) were modified from the previously published protocol described in reference 4 in the main text.

Materials

Cells

-  Human Dermal Fibroblasts, neonatal hNDF **Invitrogen Catalog #C0045C** , passages 5–11
-  RFP Expressing Human Umbilical Vein Endothelial Cells **Angio-Proteomie Catalog #cAP-0001RFP** , passages 2–11
-  HUVEC – Human Umbilical Vein Endothelial Cells **Lonza Catalog #C2517A** , passages 2–11
-  THP-1 monocytes **ATCC Catalog #TIB-202** , passages 2–25

Note

Alternative cell sources may be used for fabrication (see reference 4 in the main text).

Cells culture supplies

- hNDF growth media:
 -  DMEM, high glucose, no glutamine **Gibco - Thermo Fisher Scientific Catalog #11960044** with 4.5 g/L D-glucose supplemented with 10%
 -  Fetal Bovine Serum, certified, heat inactivated **Gibco - Thermo Fisher Scientific Catalog #10082147** , 1%
 -  Penicillin-Streptomycin **Gibco - Thermo Fisher Scientific Catalog #15140122** , 1×
 -  MEM Non-Essential Amino Acids **Gibco - Thermo Fisher Scientific Catalog #11140050** , 1×
 -  Sodium Pyruvate **Gibco - Thermo Fisher Scientific Catalog #11360070** , 1×
 -  GlutaMAX™ Supplement **Gibco - Thermo Fisher Scientific Catalog #35050061** and 0.1%
 -  2-Mercaptoethanol **Gibco - Thermo Fisher Scientific Catalog #21985023**
- HUVEC and RFP-HUVEC media:  Human Endothelial Cell Media **Cell Applications, Inc. Catalog #211-500** supplemented with 1%  Penicillin-Streptomycin **Gibco - Thermo Fisher Scientific Catalog #15140122**
- THP-1 cell media:
 -  RPMI 1640 Medium (ATCC modification) **Gibco - Thermo Fisher Scientific Catalog #A1049101** supplemented with 10%
 -  Fetal Bovine Serum, certified, heat inactivated **Gibco - Thermo Fisher Scientific Catalog #10082147** , 1%
 -  Penicillin-Streptomycin **Gibco - Thermo Fisher Scientific Catalog #15140122**
-  Nunc™ EasYFlask™ Cell Culture Flasks **Thermo Fisher Scientific Catalog #156499**
- 1X  PBS, pH 7.4 **Gibco - Thermo Fisher Scientific Catalog #10010023**

-  Trypsin-EDTA (0.25%), phenol red **Gibco - Thermo Fisher Scientific Catalog #25200056** : diluted 1:5 in sterile PBS without Ca⁺⁺ and Mg⁺⁺
-  Bovine Gelatin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G1393**

Materials for Polydimethylsiloxane (PDMS) clamps fabrication and sterilization

- SYLGARD 184 silicone elastomer (PDMS, Dow, #2646340)
-  Sterilization Pouch **Cardinal Health Catalog #92510**
- Digital lab oven (VWR # 97025-630)
-  Vacuum desiccator **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Z119008**
-  Razor Blades **VWR International (Avantor) Catalog #55411-050**

Materials for TEBV mold chamber and perfusion/viewing chamber

- Six computer numerical control machined (CNC) polycarbonate chamber pieces (Protolabs) and one CNC machined viewing top (Protolabs) for viewing chamber (Fig 1)
- Eight stainless steel hypodermic tubes of 0.8" length (20XTW - New England Small Tube Corporation)
- Four stainless steel hypodermic tubes of 3" length referred to as the TEBV mandrels, will be used in the following sections (23RW - New England Small Tube Corporation)
- 22 mm × 40 mm #1.5 Coverglass (VWR, #48393-172), cut into 22mm × 19mm to fit the viewing chamber
- Buna-N O-Rings
-  Durometer 70A (Medium), 2 mm wide and 28.5 mm inner diameter **McMaster-Carr Catalog #9262K679**
- 1/16" Tube Adapter with 10-32 Threaded Pipe Fitting (McMasterCarr, #2974K123)
- Four pedicle screws with four nuts (6/32 Thread Size 1" long, McMasterCarr, #99607A128 and #91240A007)
- 2 mm  Standard Biopsy Punches **Fisher Scientific Catalog #12-460-399**
-  PTFE tape **Merck MilliporeSigma (Sigma-Aldrich) Catalog #20808-U**
- Epoxy (Henkel, #235033)

Materials for preparing collagen matrix and harvesting hNDFs

-  Collagen I, High Concentration, Rat Tail, 100 mg **Corning Catalog #354249**
- 10X  DMEM - low glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2429**
- 1M  Sodium hydroxide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S8045** , sterile, prepared by dissolving NaOH pellets in deionized water and sterile filtering after cooling
- 2.5 mL  Microcentrifuge Tubes: 2.0 mL **Fisher Scientific Catalog #05-408-138**
-  LSE™ Mini Microcentrifuge **Corning Catalog #6770**

Materials for TEBV plastic compression/dehydration and endothelialization

- Two  50 mL PP Centrifuge Tubes **Corning Catalog #430828**
- One  1 mL Tuberculin Syringe **BD Biosciences Catalog #309659**
- One  5 mL Luer-Lok™ Syringe **BD Biosciences Catalog #309646**
- 5–6 Kimwipes
- 2 cm × 2 cm Kimwipes (cut)
- One 5 mL syringe adapter (McMasterCarr, #51525K281)
- One 1mL syringe adapter (Strategic Applications Inc., #B23-50)
- As syringe tubing  L/S® Precision Pump Tubing **VWR International (Avantor) Catalog #MFLX96410-13**
- One  Straight Tapered Fine Tip Forceps **Fisher Scientific Catalog #16-100-113**
-  Petri Dishes **Fisher Scientific Catalog #FB0875711**
- Customized rotator with 24 rotations per hour (LABQUAKE SHAKER, Barnstead)
- Mechanical Accurate Countdown Timer with 15 Minute Increments (BN-LINK)

Materials for TEBV Perfusion

- Stainless steel tray (SouthPointe Surgical Supply, #RT-1350S)
- Y-shaped connectors, for 1/1600 tube (McMasterCarr, #5117K65)
-  Masterflex® Microbore Transfer Tubing, Tygon **VWR International (Avantor) Catalog #MFLX06419-03**
-  L/S® Precision Pump Tubing **VWR International (Avantor) Catalog #MFLX96410-13**
-  Masterflex® Ismatec® Pump Tubing, 3-Stop, PharMed® BPT, 0.89 mm **VWR International (Avantor) Catalog #MFLX95714-26**
- Two-way straight connectors (McMaster-Carr, #5117K41)
- Masterflex L/S Precision Modular Drive pump (VWR, #MFLX07557-10)
-  Masterflex® Ismatec® Minicartridge Multichannel Pump Head for Masterflex® L/S® Drives, 8-Channel, 8- **VWR International (Avantor) Catalog #MFLX07623-10**
-  Syringe Filters, 0.22 µm **VWR International (Avantor) Catalog #76479-044**
-  25 mL Round Media Storage Bottles, with GL25 Screw Cap **Corning Catalog #1395-25**
-  Glass Glue **Loctite Catalog #233841**

Materials for evaluating monocyte adhesion in TEBVs and immunofluorescence staining

- Human monocytic THP-1 cells, passages 5–20
- 1X  PBS, pH 7.4 **Gibco - Thermo Fisher Scientific Catalog #10010023**
-  16% Paraformaldehyde **Electron Microscopy Sciences Catalog #15710** , diluted 1:4 in PBS without Ca++ and Mg++

-  Goat serum **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G9023**
- 10% BSA: 5 g  Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9647** in 50 mL PBS.
- 0.1% Triton-X solution: 10 uL  Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #X100** in 10 mL PBS
-  Tween 20 100% Nonionic Detergent **Bio-Rad Laboratories Catalog #1706531**
- 30% sucrose solution: 30 grams of  Sucrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S0389** in 100 mL PBS.
-  Olympus 5ml Centrifuge Tubes, Clear Natural, Polypropylene **Genesee Scientific Catalog #24-285**
-  ImmEdge hydrophobic barrier pap pen **Vector Laboratories Catalog #H-4000**
-  CD54/ICAM-1 (VF27-516) Mouse mAb **Cell Signaling Technology Catalog #62133S**
-  Anti-von Willebrand Factor antibody **Merck MilliporeSigma (Sigma-Aldrich) Catalog #F3520**
-  VE-Cadherin **Cell Signaling Technology Catalog #2500S**
-  Hoechst 33342, Trihydrochloride, Trihydrate - 10 mg/mL Solution in Water **Thermo Fisher Scientific Catalog #H3570**
-  CellTracker™ Red CMTPIX **Thermo Fisher Scientific Catalog #C34552**
-  Microscope Slides, Diamond White Glass, 25 × 75mm, Charged, 90° Ground Edges, White Frosted **Globe Scientific Catalog #1358W**
-  Micro cover glasses **VWR International (Avantor) Catalog #48368-040**
-  Extra Fine Bonn Scissors **Fine Science Tools Catalog #14084-08**
- Tissue-Tek O.C.T. Compound (Sakura Finetech, #4583)
-  Truncated Shape Mold, 12X12mm **Electron Microscopy Sciences Catalog # 70181**
- Leica SP5 inverted confocal microscope
- Zeiss Axio Imager fluorescence microscope (5X, 10X objectives) outfitted with an AxioCam MRm digital camera with ZEN Pro software




Fig 1. Design and schematic of the TEHV chambers. The mold, perfusion and viewing chamber are assembled by combining various top and bottom pieces with the core chamber component.

Troubleshooting

Procedure 0: PDMS clamps fabrication and sterilization of TEBV chamber components (Timing: 2-3 days before the TEBV fabrication)

21h 40m

1 PDMS clamp mold fabrication

Print molds for the PDMS clamps using a Stratsys J750 3D printer with the proprietary VeroPure White filament acrylic formulation. Clean molds thoroughly with sodium hydroxide prior to casting with PDMS.

2 PDMS clamps fabrication



Note

CRITICAL STEP Clean the master mold by removing excess PDMS to ensure a flat fit on the dish surface.

Prepare SYLGARD 184 silicone Elastomer (PDMS) mixture by adding 5 parts of the elastomer base and 1 part of the curing agent into a 50 mL falcon tube. Inverting the tube several times to ensure thorough mixing of the two components.

2.1 Centrifuge the tube (1400 rpm, 00:05:00) and degas the mixture for 01:30:00 .

1h 35m



Note

CRITICAL No visible bubbles should be present in the PDMS solution.

2.2 Slowly pour the PDMS into the dish containing the master mold of TEBV clamps. Ensure there is at least 5 mm of PDMS above the top of the molds.

2.3 Place PDMS-coated molds under vacuum for at least 02:00:00 to ensure thorough degasification.

2h

2.4 After degassing, remove molds from vacuum. Using forceps, gently press down so that each is flat against the bottom of the petri dish.

2.5 Place PDMS-coated molds in an oven set at  50 °C and let them incubate

8h

 Overnight



Note

CRITICAL Ensure the oven temperature does not exceed  55 °C to avoid reaching the glass transition temperature of the VeroPure White molds.

2.6 Once cured, retrieve the molds from the oven and collectively cut them out of the PDMS using a scalpel.

2.7 Separate each mold with a straight edge razor, eliminating excess PDMS to leave only 1-2 mm on each side of the molds.

2.8 With fine-tip forceps, delicately detach the PDMS from all four sides of each mold, taking care to avoid splitting or breaking any part of the PDMS.

2.9 After separating PDMS from all sides, slowly lift one end of the mold out of the PDMS encasement using fine-tip forceps.



Note

CRITICAL Proceed with caution to maintain fidelity of the clamp grooves and prevent PDMS breakage or splitting.

2.10 Upon complete removal of the mold, trim away the excess PDMS layer from the top of the clamp.

2.11 Use a straight edge razor to longitudinally cut the clamps on each side, exposing the grooves. Exercise care by positioning the razor appropriately and pressing straight down, ensuring the clamps are approximately 5 mm in width and avoiding excessive cuts on each side.

2.12 Utilizing a straight edge razor, trim the lateral sides of the clamp, maintaining a gap of approximately 1 mm between the edge of each groove and the clamp's edge.

- 2.13 Rotate the clamp 90 degrees, and employing the same technique used for removing excess PDMS from the longitudinal sides, reduce the height of the clamp to 2.75-3 mm.
- 2.14 The clamps are now ready for steam sterilization.

3 Fabrication of PDMS side port caps

Prepare SYLGARD 184 silicone Elastomer (PDMS) mixture by adding 5 parts of the elastomer base and 1 part of the curing agent into a 50 mL falcon tube. Inverting the tube several times to ensure thorough mixing of the two components.

- 3.1 Centrifuge the tube ( 1400 rpm, 00:05:00) and slowly pour the PDMS into a petri dish.

5m



- 3.2 Degas the dish with PDMS mixture for  01:30:00 -  02:00:00 .

2h



Note

CRITICAL No visible bubbles should be present in the PDMS solution.

- 3.3 Preheat the oven to  50 °C and let PDMS mixture incubate  Overnight .

8h



- 3.4 After curing, carefully remove the dish from the oven. Use a blade to cut out a piece of PDMS at least 0.5 cm in depth.
- 3.5 Employ a 2 mm biopsy punch to create a space for covering side ports in the PDMS material.
- 3.6 The side port caps are now ready for steam sterilization.

TEBV mold chamber assembly and materials preparation

21h 40m

4 Wrap Teflon tape around the two side ports and install them into the side holes on each side of the core chamber component.

4.1 Use epoxy to attach the eight 0.8" stainless steel tubes to the core chamber component and allow the epoxy to cure at  Room temperature .



Note

CRITICAL Ensure that the distance between the two ends of the tube on the opposing sides is approximately 1.2 cm, as this determines the length of the TEBVs.

4.2 Fit the silicone O-rings into the rectangular groove on each side of the core chamber component.



Note

CRITICAL STEP Ensure that the O-rings are sealed into the grooves on either side of the core chamber component. Replace O-rings and epoxy adhesive after 3-4 uses or if they appear worn.

4.3 Assemble the top and bottom pieces of the TEBV mold chamber to the assembled core chamber component and then insert four TEBV mandrels into the chamber.



Note

CRITICAL Make sure that the mandrels do not have any sharp areas that may tear the TEBVs.

4.4 Sterilize all the following materials prior to use (Table 1). All the materials are placed into autoclave pouches unless otherwise indicated.

	A	B
	Procedure 1	Assembled TEBV mold chamber (mold chamber top and bottom, core chamber component, four screws and nuts, four TEBV mandrels, and two PDMS side port caps)

A	B
Procedure 2	Viewing/ Perfusion chamber top and bottom; four PDMS clamps; small tubes with varying lengths to replace dysfunctional TEBVs; one pair forceps; stainless steel surgical tray; small square Kimwipes for dehydration; 1 mL syringe with a 23G needle; 5 mL syringe connected to tubing through a luer adapter; Kimwipes for cleanup
Procedure 3	1 mL syringe with a 23G needle; one pair forceps
Procedure 4	Vessel loop tubing; side loop tubing; media reservoir; two pairs of forceps; 5 mL syringe connected to tubing through a luer adapter; Kimwipes for cleanup
Procedure 5	Two pairs of forceps; 5 mL syringe connected to tubing through a luer adapter; Kimwipes for cleanup

Table 1: Materials list for sterilization.

Note

All the following procedures should be performed in Class II Biosafety cabinets unless specified.

Procedure 1: Harvest hNDFs and prepare collagen mixture, Day 0 (Timing 1 h)

1h 10m

- 5 Prepare a collagen solution with a final concentration of 7 mg/ml . To achieve this, calculate the required volumes of collagen, 10× DMEM, 1M NaOH, and hNDF media needed for dilution.



The formula for 1 mL collagen and cell mixture for the four TEBV system is

- Volume of 10× DMEM = $94 \mu\text{L}$
- Volume of collagen solution (V_c) = $1 \text{ mL} \times (7 \text{ mg/ml} / \text{collagen stock concentration})$
- Volume 1 M NaOH (V_{NaOH}) = Volume of collagen solution $\times 0.023$
- Volume of hDNF suspension = $60 \mu\text{L}$
- Volume of hNDF media (V_{hNDF}) = $1 \text{ mL} - 94 \mu\text{L} - V_c - V_{\text{NaOH}} - 60 \mu\text{L}$
(if negative then do not add this extra media)

Note

CRITICAL May need to aliquot extra volume (>1 mL) to compensate for pipetting error and possible bubble formation.

- 6 Aliquot the calculated volumes of collagen, 10X DMEM, 1M NaOH solutions in separate sterile 2.5 mL Eppendorf tubes and keep  On ice . 

Note

CRITICAL STEP Avoid making bubbles while pipetting collagen solution.

- 7 Tighten the nuts on the pre-assembled TEBV mold chamber using pliers. Place on autoclaved surgical tray. 

Note

CRITICAL STEP Adjust the tightness of the screws in a diagonal, alternating fashion and avoid touching the mandrels. Uneven tightening will cause deformation of the O-ring.

- 8 Add the calculated volumes of collagen, 10 ×DMEM, 1M NaOH and hNDF media into a sterile 2.5-mL Eppendorf tube and mix gently.   

Note

CRITICAL STEP Avoid making bubbles while pipetting.

- 9 Harvest hNDFs by trypsinization and centrifugation. Treat hNDFs (~80% confluent) with trypsin-EDTA for  00:03:00 , then add  3 mL fresh DMEM media to neutralize trypsin-EDTA. Centrifuge the cells ( 180 x g, 00:07:00). Resuspend the cell pellet in a density of 1×10^6 hNDFs/60 μ L DMEM media.  

- 10 Quickly add  60 μL of hNDF suspension and mix thoroughly by gently swirling the cells in the gel mixture. Then add extra volume of 1M NaOH to adjust pH until the mixture turns magenta (as demonstrated in the video provided in S3).



Note

CRITICAL STEP Use a benchtop mini centrifuge to briefly centrifuge the collagen mixture before adding into the mold chamber. Act quickly to avoid collagen gelation.

- 11 Use 200 μL pipette to slowly add the collagen mixture into each of the 4 vessels of the mold chamber. Stop injecting when the collagen mixture begins to extrude from the contralateral opening.



Note

CRITICAL Using a P200 pipette, load  200 μL of the collagen solution, even though each vessel only requires approximately 180 μL . Loading slightly more than the vessel mold capacity ensures that there is enough material to fully fill the mold and compensate for any material lost due to potential bubble formation during injection. Ensure there are no significant bubbles in the mold.

- 12 Put the mold in a petri dish and seal with parafilm. Incubate in the cell culture incubator for  00:45:00 -  01:00:00 .

1h

?TROUBLESHOOTING

Procedure 2: Plastic compression and dehydration, Day 0 (Timing: ~45 minutes)

1m

- 13 Remove petri dish containing TEBV mold chamber from incubator to biosafety hood. Remove the chamber from the petri dish. Use pliers to unscrew the nuts and peel off the top and bottom pieces of the mold chamber gently.



Note

CRITICAL Avoid touching mandrels on each side during the manipulation to prevent damage to the lumens and contamination.

- 14 Use forceps to gently apply 2 cm x 2 cm pre-cut sterile Kimwipes directly onto the TEBVs. The Kimwipes should become soaked through as part of the dehydration process. Apply double-layer Kimwipes on each side (~3 times per side) followed by single-layer Kimwipes (~1-2 times per side). Stop the dehydration process once the water absorption on the wipe is isolated to individual TEBVs.



Note

CRITICAL STEP Avoid pressing directly on TEBVs during dehydration. Use forceps to gently press the space between TEBVs to ensure uniform fluid dehydration.

?TROUBLESHOOTING

- 15 Assemble the bottom piece of the TEBV viewing chamber and hand tighten the screws. Add  2 mL -  2.5 mL warm media into the chamber using a P1000 pipette to soak the dehydrated TEBVs for at least  00:01:00 . Carefully remove the media from chamber using a P1000.

1m



Note

CRITICAL STEP Put pipette tip at the corner of the chamber and slightly tilt the chamber to pull-out media. Avoid touching TEBVs!

- 16 Firmly secure TEBVs to the mandrels with PDMS clamps using forceps.
- 17 Hand tighten the top layer of the perfusion chamber first and use plier to gradually tighten the four corners evenly. Tighten the screws in a diagonally alternating pattern.

?TROUBLESHOOTING

- 18 Prime a sterile 5 mL syringe with luer adaptor with warm hNDF media. Invert the syringe to ensure that the media coming out from the syringe is bubble-free. Fill the chamber with media through the barbed side port using the primed syringe. Tilt the chamber while injecting the media to expel as many bubbles as possible. Cap the side ports using PDMS caps.



Note

CRITICAL Use forceps to connect the tube on syringe to the side port to avoid contamination.

- 19 Remove all the inserted 23RW mandrels carefully using forceps.



Note

CRITICAL Pull out the mandrel horizontally, slowing down at both ends to avoid ripping the TEBV lumens during mandrel removal.

- 20 Perform a leakage test for all the vessels.



Note

CRITICAL Use a 1mL syringe with a BD syringe adaptor filled with PBS for easy observation. Be attentive to bubbles traversing the TEBV and examine for PBS emerging from the opposite side. If leakage is observed, replace the TEBV with a short section of tubing (see [↩ go to step #22](#) for details). **?TROUBLESHOOTING**

- 21



Note

OPTIONAL Place the mold into a petri dish and seal it using parafilm. Inspect each vessel under a phase contrast microscope (5X objective). Look for clear demarcation of the lumen and an intact collagen layer. Replace any vessels showing signs of rupture or a collapsed lumen, as outlined in the step 22 .

- 22 **TEBV replacement:** Place the chamber on sterile absorbent wipes. Choose the suitable length of sterilized Tygon tubing prior to opening the chamber. Use pliers to open the perfusion chamber. Gently detach the dysfunctional vessel from the connected mandrels using tweezers, ensuring careful handling to prevent displacement of PDMS clamps. Use tweezers to affix the tubing. Seal the perfusion chamber and replenish the chamber with media through the barbed side port.



Note

CRITICAL Execute the replacement procedure quickly to prevent adjacent vessel damage.

Procedure 3: Endothelial cell seeding and rotation, Day 0-1 (Timing: ~ 45 minutes)

8h 8m

- 23 Prepare a HUVEC suspension of 4.5×10^6 cells in  600 μL warm EC media. Harvest the cells by using trypsin-EDTA for  00:03:00 (~80% confluent), followed by the addition of  3 mL of fresh EC media to neutralize the trypsin-EDTA. Centrifuge the cells at  1000 rpm, 00:05:00 .

8m



Note

CRITICAL STEP Prepare extra volume of HUVEC suspension for each chamber to compensate for the dead volume of the syringe (~  100 μL).

- 24 Fill a sterile 1-mL syringe with HUVEC suspension and tap to avoid any bubbles until a meniscus is formed at the top of the syringe.
- 25 Attach the syringe to one end of the vessel mandrel using forceps. Gently push the HUVEC suspension through, perfusing each vessel with ~  100 μL of the HUVEC suspension (seeding density: 1.5×10^5 cell/cm²). The cloudy HUVEC suspension drops emerging from the opposite side of the mandrel indicate the presence of HUVECs.
- 26 Repeat the step 25 for all the vessels.



Note

CRITICAL STEP Tap the syringe for before each injection to ensure that the HUVEC suspension is well mixed.

- 27 Place the chamber in a dish and seal it with parafilm. Securely affix the chamber to the rotator using lab tape, ensuring that the vessels are parallel to the rotator pole.

28 Set the timer to 12 hours and keep the chamber in the incubator  Overnight .

8h



?TROUBLESHOOTING

Procedure 4: Perfusion set up, Day 1 (Timing: ~ 3.5 hours)

29 There are two perfusion loops for the TEBV chamber as shown in the Fig 2: the vessel loop and the side loop. All the numbers in the parentheses refer to the components in the Fig 2.

Note

CRITICAL All tubing should be sterilized before TEBV perfusion.




Fig. 2 Perfusion System Configuration. The tubing and pump arrangement for TEBV perfusion. Ensure all components undergo sterilization and leak testing prior to perfusion. For the output vessel loop, substituting the Ismatec 3-stop Microbore Tubing with L/S® Precision Pump Tubing. Refer to the materials section for detailed product specifications.



30 Side Loop Tubing Assembly:

Prepare two pieces of L/S Precision pump tubing with desired length to connect to the MasterFlex pump head. Attach 3-stop microbore tubing to each end of the tubing using a two-way straight connector. Connect the prepared pump tubing to the two-way straight connector (Fig 2, part 1-2-3).

31 Vessel Loop Tubing Assembly:

Prepare eight pieces of L/S Precision pump tubing and four pieces of Tygon tubing with desired length. Connect a 3-stop microbore tubing with a two-way straight connector on each end of the tubing (Fig 2, part 11-12). Connect one end of the two-way straight connector with a L/S Precision pump tubing (Fig 2, part 10-11-12). Assemble branched tubing by connecting each branch of the Y-shaped connector with L/S Precision pump tubing (Fig 2, part 10-13). Insert a portion of the small metal mandrel into each Tygon tubing and connect the end branch of the L/S Precision pump tubing, repeating this step four times (Fig 2, part 14-15). Connect branched tubing to the other two-way straight connector. For the outlet tubing of the vessel loop, replace the 3-stop microbore tubing with a L/S Precision pump tubing.

32 Side Perfusion Loop Setup:



32.1 Connect one end of the side tubing to the chamber side port first, then prime side loop tubing with warm EC media.

Note

CRITICAL Clamp the other end of tube securely with forceps to avoid bubbles.

32.2 Use forceps to connect the side loop tubing to the side port.

33 Vessel Perfusion Loop Setup:

- 33.1 First, prime the outlet tubing by filling it with media and then connect it to the perfusion chamber.
- 33.2 Prepare the inlet tubing by filling it with media and secure it by clamping it at the Y-shaped tubing. Any small bubbles trapped in the tubing can be eliminated by either pulsating the syringe plunger or gently squeezing the tubing.
- 33.3 Connect each branch of the inlet tubing to the mandrels one by one using forceps.
- 33.4 Use forceps to insert the other end of the inlet tubing into the media reservoir.
- 33.5 Attach a 0.22 μm syringe filter to the male luer connector on the cap of media reservoir.
- 33.6 Set up the pump and start the perfusion at the flow rate of $70 \mu\text{L} / \text{min}$ per TEBV (the readout pump number is 5 rpm). Monitor the media flow from each vessel and replace the vessel if any abnormalities are observed. Adjust the flow rate hourly until it reaches $500 \mu\text{L} / \text{min}$ per TEBV (the readout pump number is 40 rpm).



Note

CRITICAL STEP Make sure the flow direction is from the inlet tubing to the outlet.
? TROUBLESHOOTING

- 33.7 Calculate the wall shear stress τ_w in the TEBV for a given volumetric flow rate ($Q, \text{m}^3/\text{s}$) using the equation $\tau_w = 4\mu Q/\pi R^3$, where R (m) is average TEBV inner lumen radius and μ ($\text{Pa}\cdot\text{s}$) is Newtonian fluid of viscosity for cell culture media at 37°C .

Procedure 5: Daily media change (Timing: 5-10 minutes)

- 34 Disconnect the tubes from the pump and carefully transfer the system into hood.

- 35 Disconnect one end of the side loop from the barbed side port and place the tube in the waste container.
- 36 Fill one 5 mL syringe with warm EC media and tap the syringe to remove any bubbles.
- 37 Use forceps to connect the syringe to the barbed side port and gradually push fresh EC media through the chamber. Meanwhile, expel any bubbles that may form in the chamber during perfusion.
- 38 Securely clamp the other end of the tube with forceps to prevent the backflow of media. Use forceps to reconnect the tube to the side port.
- 39 Clamp the inlet tube and aspirate the media in the reservoir.
- 40 Add - fresh EC media to the reservoir and then remove the clamp. 
- 41 Reconnect the tubes to the pump and monitor the perfusion system to ensure no bubbles are visible in the perfusion loop. 

Note

CRITICAL Make sure to avoid introducing bubbles into the perfusion loops during media changes as this can result in EC detachment and possible TEBV loss.

Procedure 6: Downstream experiments

4h 35m

42

Note

The following procedures can be selected and modified according to specific experiment needs. They are not included in the TEBV fabrication and daily maintenance procedures.

43 **Activation of endothelium and THP-1 monocytes perfusion (Timing: 4.5 hours)**

4h

Change the perfusion media to fresh EC media containing 200U/mL TNF- α , then perfuse the TEBVs at  500 μ L /min per TEBV for  04:00:00 .

43.1 During the 4-hour endothelium activation step, stain 3×10^6 THP-1 cells for each chamber of TEBVs using 1 μ M cell tracker red-CMTPX in warm RPMI media without serum for  00:30:00 .

30m

43.2 **CRITICAL** Following incubation, centrifuge THP-1 cell suspension ( 1000 rpm,  00:05:00) and wash cells with fresh RPMI media to remove any residual staining reagent.

5m



43.3 Change the media in the TEBV reservoir before adding stained monocytes. Add THP-1 cells into the main vessel loop at a concentration of 1×10^6 cells/mL, ensuring a total volume of  3 mL in the media reservoir.



Note

CRITICAL Avoid adding monocytes in a large volume to prevent cell deposition in the reservoir.

TEBV fixation (Timing: 1.5 hours)

1h 5m

44 Disconnect the chamber from perfusion pump and place the chamber in a chemical fume hood.

45 Load a 5 mL syringe with a luer adaptor attached to a tube with 4% PFA and substitute the medium in the perfusion chamber through the side port. Incubate the vessels for  00:05:00 while exchanging the medium in the reservoir with 4% PFA.

5m



46 After the 5-minutes incubation, reconnect the TEBVs to the pump in the chemical fume hood and perfuse vessel lumen with 4% PFA for an additional 5 minutes. Gently open the perfusion chamber and carefully detach the TEBVs from the mandrels using tweezers.

47 Transfer the TEBVs to a 6-well plate and continue fixation (4% PFA) for  01:00:00 .
Rinse the TEBVs three times with PBS and store in PBS at  4 $^{\circ}$ C .

1h



Procedure 6: Downstream experiments

48 Embedding and sectioning of TEBVs (Timing: 2 days)

Dehydrate each TEBV in a 5 mL centrifuge tube filled with 30% glucose solution until the vessel sinks at  4 °C .

48.1 Remove excess collagen at ends of TEBV then use a 1 mL syringe with a blunt needle to inject a small amount of OCT into the vessel lumen. Stereoscope may be used here to facilitate observation.

48.2 Transfer the vessels into an embedding cup and hold the vessel vertically using forceps. Gradually fill the cup with OCT. 

Note

CRITICAL Do not clamp on the vessels while transferring the vessel into the block which leads to closed lumen. Hold the vessel only at the remaining clamped end using tweezers.

48.3 Carefully position the block flat on dry ice and wait till the OCT turns white. 

Note

CRITICAL Monitor the vessels and use tweezers to hold them in place as vessels may bend while freezing.

48.4 Store blocks at  -80 °C until ready for cryosectioning.

Immunofluorescence staining of TEBVs (Timing: 2 days)

17h 40m

49 For *en face* staining, use a blade to cut a small section of TEBV and delicately insert fine-tip scissors into the vessel lumen to make the incision. Utilize square glass slides during this process to facilitate cutting and opening of the vessels. For TEBV sections, warm the 

slides to  Room temperature and rinse off OCT using PBS. Outline the sections using an ImmEdge™ pen.

Note

CRITICAL Avoid excessive manipulation during cutting the vessels *en face* as the tools used can damage the integrity of endothelium.

50 For *en face* staining, permeabilize with 0.1% Triton-X in PBS for  00:30:00 . For TEBV section staining, permeabilize for  00:10:00 . Then rinse samples three times with PBS.

40m



51 Block the TEBV samples with blocking buffer (10% goat serum + 10% BSA + PBS with 0.1% Tween) for  08:00:00 at  Room temperature on a shaker.

9h

- For section slides, block for  01:00:00 in a humidified chamber at  Room temperature .

52 Add primary antibody at desired dilution ratio in blocking buffer and incubate at  4 °C  Overnight .

8h



53 After primary antibody staining, wash the samples three times with PBS. Subsequently, add secondary antibody at desired dilution ratio in blocking buffer.



54 Rinse the samples 3 times with PBS before being placed on a glass slide and image using fluorescence or confocal microscopy.

