Jun 29, 2023

C TELEvir Field Protocol

DOI

dx.doi.org/10.17504/protocols.io.n2bvj694xlk5/v1

Anna S Fomsgaard^{1,2}

¹Statens Serum Institut; ²University of Copenhagen



Anna S Fomsgaard



DOI: dx.doi.org/10.17504/protocols.io.n2bvj694xlk5/v1

Protocol Citation: Anna S Fomsgaard 2023. TELEvir Field Protocol . protocols.io https://dx.doi.org/10.17504/protocols.io.n2bvj694xlk5/v1

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: April 06, 2022

Last Modified: June 29, 2023

Protocol Integer ID: 60387

Funders Acknowledgement: Horizon 2020 Grant ID: 773830

Disclaimer

This protocol was prepared for a workshop held the 23-24 June 2022 at Statens Serum Institut, Denmark. Here 25 TELEVIR partners from 10 European countries participated to train in the usage of a field-deployable point-of-incidence toolbox to identify emgerging virus threats. The protocol was used to detec SARS-CoV-2 and Human Papillomavirus a RNA and DNA virus detection test-of-concept.

Abstract

The TELE-Vir project key aim is to develop a fast point-of-incidence (poi) toolbox for identification and characterisation of emerging virus threats for humans and/or domestic and wildlife animals.

Presented here is a protocol including sample pretreatment, NA extraction and random amplification for metagenomic virus detection using the MinION device (Oxford Nanopore Technologies)



protocols.io Part of SPRINGER NATURE

Protocol materials

X DNase Set | Zymo Research Catalog #E1010 Step 1

X MagNA Pure LC Total Nucleic Acid Isolation Kit Roche Catalog #03038505001 In 2 steps

- 🔀 REPLI-g Cell WGA & WTA Qiagen Catalog #150052 Step 30
- 🔀 Rapid Sequencing Kit Oxford Nanopore Technologies Catalog #SQK-RAD004 Step 37

Pretreatment

1

X DNase Set I Imerco Catalog #E1010

Equipment

BD Disposable Syringe with Luer-Lok Tip (5ml)	NAME
Syringe	TYPE
Becton Dickinson	BRAND
309649	SKU

Equipment	
Syringe Filter pore size 0.2 µm	NAME
Filter	TYPE
Minisart	BRAND
16532	SKU
https://www.sartorius.com/shop/ww/en/usd/applications-laboratory-filtra sterile-filtration/minisart-syringe-filter%2C-polyethersulfone-%28pes%29% size-0-22-%C2%B5m%2C-ethylene-oxide%2C-female-luer-lock%2C-male-lue lock%2C-pack-size-50/p/16532-	ation- LINK 2C-pore- er-
Polyethersulfone (PES), Pore Size 0.22 μm, Ethylene Oxide, Female Luer Lock, Male Luer Lock, Pack Size 50	SPECIFICATIONS

🕹 50 µL DNase I

 \underline{A} 50 µL DNase Digestion buffer

3 00:15:00 at room temperature

Note

If 400 μL sample is not available, downscale the reagents to fit the lower sample volume.

2 Add [500 µL PBS (or equivalent to 1 ml sample material in total)

15m

15m

protocols.io Part of Springer NATURE

- 3 Poor the diluted sample into the lid of the Syringe Filter (Minsart). The filter can be placed on a sterile surface meanwhile e.g. the Syringe Filter paper lid (inner side).
- 4 Suck in air corresponding to app 1 ml in a 5 ml syringe
- 5 Extract the 1 ml sample material from the petri dish to the 5 ml syringe
- 6 Attach the 0,22 μ M Minisart syringe filter to the 5 ml syringe with Luer-Lok
- 7 Filter the sample and air directly into a 4.5 ml cryotube containing 1 ml MPLB-buffer. Put lid on, mix by turning tube.

X MagNA Pure LC Total Nucleic Acid Isolation Kit Imerco Catalog #03038505001

Inactivation

8 Incubate

O0:20:00 at room temperature for viral inactivation

CITATION

Vinner L, Fomsgaard A (2007). Inactivation of orthopoxvirus for diagnostic PCR analysis.. Journal of virological methods.

CITATION

Rosenstierne MW, Jensen CE, Fomsgaard A (2018). Rapid, Safe, and Simple Manual Bedside Nucleic Acid Extraction for the Detection of Virus in Whole Blood Samples.. Journal of visualized experiments : JoVE. LINK

https://doi.org/10.3791/58001

Extraction

9 Kit used for hand held NA extraction:

X MagNA Pure LC Total Nucleic Acid Isolation Kit Imerco Catalog #03038505001

10m

20m

20m

Equipment	
4.5 ml cryotube	NAME
tube	TYPE
Thermo Fisher Scientific	BRAND
363452	SKU

https://www.thermofisher.com/order/catalog/product/363452LINK

Equipment	
3.6 ml cryotube	NAME
Tube	TYPE
Thermo Fisher Scientific	BRAND
379189	SKU
https://www.thermofisher.com/order/catalog/product/363452 ^{LINK}	

Cylinder magnet

A rubber band to attach the magnet to a tube.

A 50 ml Nunc tube or similar to pour excess solutions into for disposing.

Following the protocol by:

CITATION

Rosenstierne MW, Jensen CE, Fomsgaard A (2018). Rapid, Safe, and Simple Manual Bedside Nucleic Acid Extraction for the Detection of Virus in Whole Blood Samples.. Journal of visualized experiments : JoVE. LINK

https://doi.org/10.3791/58001

10 Per X number of samples, aliquot in cryotubes or Eppendorf tubes

▲ 960 µL Magnetic Glass Particles (MGPs) in an Eppendorf tube

4 mL Wash buffer I in a 4.5 ml cryotube

🕹 1.5 mL Wash buffer II in an Eppendorf tube

- 🗸 3 mL Wash buffer III in a 3.6 ml cryotube
- Δ 100 μL Elution buffer in an Eppendorf tube

Prepare one sample at the time.

10m

Approximate time 🚫 00:10:00 per sample

- 11 Add 960 µl MGPs to 2 ml solution of sample and MPLB-buffer
- 12 Put lid on and turn tube in hand gently until well mixed
- 13 Attach magnet and trap beads.
- 14 Pour solution into trash tube. The MGPs should stay behind trapped on the side of the tube.
- 15 Pour Wash buffer I in the sample tube
- 16 Put a lid on, remove magnet and turn the tube in hand gently until well mixed
- 17 Attach magnet and trap beads.
- 18 Pour solution into trash tube. The MGPs should stay behind trapped on the side of the tube.
- 19 Pour Wash buffer II in the sample tube
- 20 Put a lid on, remove magnet and turn the tube in hand gently until well mixed
- 21 Attach magnet and trap beads.
- 22 Pour solution into trash tube. The MGPs should stay behind trapped on the side of the tube.
- 23 Pour Wash buffer III in the sample tube
- 24 Put a lid on, remove magnet and turn the tube in hand gently until well mixed
- 25 Attach magnet and trap beads.
- 26 Pour solution into trash tube. The MGPs should stay behind trapped on the side of the tube.
- 27 Pour Elution buffer in the sample tube

28 Put a lid on, remove magnet and turn the gently until beads are eluted.

29 Transfer 🚡 10 µL Elution Buffer with MGPs to **2 x** 0.2 ml PCR tube

Whole Transcriptome and Genome Amplification

30 Equipment to enable portability

- A MiniPCR to run isothermal incubations

- A powerbank that can provide 20V. Use one to power the miniPCR and another to power the

MinION for approximately 5 hours. Use brand of choice.

- A salad swing as a hand driven centrifuge. Use brand of choice.

- A coolbox for reagents and cooling samples. Use brand of choice.

Equipment	
All-in1 Laptop Powerbank 24000	NAME
Power Bank	TYPE
Sanberg	BRAND
420-57	SKU
https://sandberg.world/da-dk/product/all-in-1-laptop-powerbank-2400	LINK
Any powerbank that can produce 20V is acceptable. Have 2-3 to run the MiniPCR and the MinION	SPECIFICATIONS

NAME

TYPE

BRAND

SKU

LINK

Equipment	
miniPCR® mini8 thermal cycler	
Thermal cycler	
miniPCR®	

mini8		
https://www.minipcr.com/		



4h

Equipment

A Salad Spinner Centrifuge

Centrifuge

undefined

undefined

https://www.oxo.com/salad-spinner.html? queryID=975936c98686e2b0325b684dd613daff&objectID=2251&indexName=magento2_pro



Kit used for WTA and WGA

X REPLI-g Cell WGA & WTA Imerco Catalog #150052

Quantiscript RT Enzyme Mix, Ligase Mix, and REPLI-g SensiPhi DNA Polymerase should be thawed on ice.

All other components can be thawed at room temperature (15-25°C)

31 The two aliquots of same sample enables uniform whole genome amplification (WGA) and whole transcriptome amplification (WTA) in parallel reactions.

32 Cleanup for WTA

10m

Add $\[\] _$ 2 μ L gDNA Wipeout Buffer to PCR tube with 10 μ l eluted sample material

€ 00:10:00 **\$** 42 °C on MiniPCR

33 Repair (WGA) and Reverse Transcription (WTA)

33.1 For WGA, mix following and add 🚨 10 µL Repair reagent mix to sample

_	Reagents	Volume (µL) per sample
	RT/Polymerase Buffer	4
	gDNA Wipeout buffer	2
_	H2O sc	1
_	Random primer	1
	Random Hex-P primer (20µM)	1
	WGA Ready Enzym	1
	Total	10

This modified protocol uses 5'-phosphorylated random hexamers (P-N6) instead of kit provided oligo-dT primers according to Rosenstierne et al (2014).

CITATION

Rosenstierne MW, McLoughlin KS, Olesen ML, Papa A, Gardner SN, Engler O, Plumet S, Mirazimi A, Weidmann M, Niedrig M, Fomsgaard A, Erlandsson L (2014). The microbial detection array for detection of emerging viruses in clinical samples--a useful panmicrobial diagnostic tool.. PloS one.

LINK

https://doi.org/10.1371/journal.pone.0100813

33.2 For WTA, mix following and add 📕 8 µL RT reagent mix to sample

_		
	Reagents	Volume (µL) per sample
	RT/Polymerase Buffer	4
	H2O sc	1
	Random primer	1
	Random Hex-P primer (20µM)	1
	Quantiscript RT Enzym Mix	1
	Total	8

33.3 Incubate

Image: Online of the second second

34 Ligation

Mix following and add 🗸 10 µL Ligation mix to each tube

Reagents	Volume (µL) per sample
Ligation buffer	8
Ligase Mix	2
Total	10

Incubate

00:30:00

🖇 24 °C This is the Repli-G recommended temperature, however, the miniPCR only goes to 25C. Instead of incubation in the MiniPCR, try table-top incubation.

€ 00:05:00 **§** 95 °C

Hold at 📲 4 °C

35 Amplification

Mix following and add 🛛 30 µL amplification mix to each tube

Reagents	Volume (µL) per sample
Repli-g reaction buffer	29

2h 5m

35m

1h 3m

Reagents		Volume (µL) per sample
Repli-g SensiPhi DNA polymerase		1
Total		30
Incubate		
02:00:00	₿ 30 °C	
00:05:00	₿ 65 °C	
Hold at 📲 4 °C	C	

36 For library preparation use amplified material **without** MGPs.

Library Preparation

37

Multiple rapid library preparation kits are available through **Oxford Nanopore Technology.**

For Workshop we use.

X Rapid Sequencing Kit Imerco Catalog #SQK-RAD004

Equipment	
MinION Mk1C	NAME
Sequencer	TYPE
Oxford Nanopore Technology	BRAND
MIN-101C	SKU
https://store.nanoporetech.com/minion-mk1c.html	LINK



protocols.io Part of SPRINGER NATURE

Citations

Step 33.1

Rosenstierne MW, McLoughlin KS, Olesen ML, Papa A, Gardner SN, Engler O, Plumet S, Mirazimi A, Weidmann M, Niedrig M, Fomsgaard A, Erlandsson L. The microbial detection array for detection of emerging viruses in clinical samples--a useful panmicrobial diagnostic tool. https://doi.org/10.1371/journal.pone.0100813

Step 8

Vinner L, Fomsgaard A. Inactivation of orthopoxvirus for diagnostic PCR analysis.

Step 8

Rosenstierne MW, Jensen CE, Fomsgaard A. Rapid, Safe, and Simple Manual Bedside Nucleic Acid Extraction for the Detection of Virus in Whole Blood Samples. https://doi.org/10.3791/58001

Step 9

Rosenstierne MW, Jensen CE, Fomsgaard A. Rapid, Safe, and Simple Manual Bedside Nucleic Acid Extraction for the Detection of Virus in Whole Blood Samples. https://doi.org/10.3791/58001