

Oct 24, 2023

Ten(10)X-compatible Combinatorial Indexing ATAC sequencing (txci-ATAC-seq)


DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp3o68vzp/v1>

Hao Zhang¹, Ryan Mulqueen², Andrew Adey², Darren Cusanovich¹

¹University of Arizona, Department of Cellular and Molecular Medicine, Tucson, AZ;

²Oregon Health & Science University



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3o68vzp/v1>

Protocol Citation: Hao Zhang, Ryan Mulqueen, Andrew Adey, Darren Cusanovich 2023. Ten(10)X-compatible Combinatorial Indexing ATAC sequencing (txci-ATAC-seq). [protocols.io https://dx.doi.org/10.17504/protocols.io.dm6gp3o68vzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp3o68vzp/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 22, 2023

Last Modified: October 24, 2023

Protocol Integer ID: 89713

Keywords: Single-cell, ATAC-seq, Chromatin accessibility, Combinatorial indexing, Molecular hashing, compatible combinatorial indexing atac, cell atac, nuclei across multiple sample, atac, combinatorial indexing, seq technique, nuclei, microfluidic barcoding, seq, tn5, based microfluidic barcoding

Funders Acknowledgements:

NIH/NIGMS

Grant ID: R35GM137896

NIH/NIGMS

Grant ID: R35GM124704

NIH/NIMH, BICCN

Grant ID: RF1MH128842

Abstract

The txci-ATAC-seq method is a large-scale single-cell ATAC-seq technique that combines the Tn5-based pre-indexing with the 10X Chromium-based microfluidic barcoding. This molecular hashing strategy enables the profiling of up to 200,000 nuclei across multiple samples in a single emulsion reaction.

Materials

Loading Tn5

■ Annealing Buffer:

	Reagent	Final Concentration	Per 10 ml
	1M Tris-HCl, pH8.0	40 mM	400 µl
	5M NaCl	50 mM	100 µl
	H2O		9.5 ml

■ Sequences of Tn5 linker oligos (The 'N' bases shown in the Tn5ME-B sequence represent the Tn5 barcodes):

	Linker Oligo	Sequence 5' → 3'
	Tn5ME-A	TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG
	Tn5ME-B	CGTGTGCTCTTCCGATCTNNNNNNNNNAGATGTGTATAAGAGACA G
	Tn5MErev	[phos]CTGTCTCTTATACACATCT

■ Sequences of Tn5ME-B barcode:

	Row	1	2	3	4	5	6	7	8	9	10	11	12
	A	GAAC CGC G	AGGT TATA	TCAT CCTT	CTGC TTCC	GGTC ACGA	AACT GTAG	GTGA ATAT	ACAG GCG C	CATA GAGT	TGCG AGAC	GACG TCTT	AGTA CTCC
	B	TGGC CGGT	CAAT TAAC	ATAA TGTG	GCG GCAC A	CTAG CGCT	TCGA TATC	CGTC TGCG	TACT CATA	ACGC ACCT	GTAT GTTC	CGCT ATGT	TATC GCAC
	C	TCTG TTGG	CTCA CCAA	TATT AGCT	CGC CGAT C	TCTC TACT	CTCT CGTC	CCAA GTCT	TTGG ACTC	GGCT TAAG	AATC CGGA	TAAT ACAG	CGGC GTGA
	D	ATGT AAGT	GCAC GGAC	GGTA CCTT	AACG TTCC	GCA GAAT T	ATGA GGC C	ACTA AGAT	GTCG GAGC	CCG CGGT T	TTAT AACC	GGAC TTGG	AAGT CCAA
	E	ATCC ACTG	GCTT GTCA	CAAG CTAG	TGGA TCGA	AGTT CAG G	GACC TGAA	TGAC GAAT	CAGT AGGC	AGCC TCAT	GATT CTGC	TCGT AGTG	CTAC GACA
	F	TAAG TGGT	CGG ACAA C	ATAT GGAT	GCG CAAG C	AAGA TACT	GGAG CGTC	ATGG CATG	GCAA TGCA	GTTC CAAT	ACCT TGGC	CTTA TCGG	TCCG CTAA
	G	GCTC ATTG	ATCT GCCA	CTTG GTAT	TCCA ACGC	CCGT GAAG	TTAC AGGA	GGC ATTC T	AATG CCTC	TACC GAGG	CGTT AGAA	CACG AGCG	TGTA GATA

	Row	1	2	3	4	5	6	7	8	9	10	11	12
	H	GATC TATC	AGCT CGCT	CGG AACT G	TAAG GTCA	TTGC CTAG	CCAT TCGA	ACAC TAAG	GTGT CGG A	TTCC TGTT	CCTT CACC	GCCA CAGG	ATTG TGAA

Isolation of nuclei from cell lines

Buffers to make beforehand

- Omni Resuspension Buffer (RSB; filter and store at 4°C):

	Reagent	Final Concentration	Per 50 ml
	1M Tris-HCl, pH 7.5	10 mM	500 µl
	5M NaCl	10 mM	100 µl
	1M MgCl ₂	3 mM	150 µl
	H ₂ O		49.25 ml

- Omni TD Buffer (filter and store at -20°C):

	Reagent	Final Concentration	Per 50 ml
	1M Tris-HCl, pH 7.5	20 mM	1 ml
	1M MgCl ₂	10 mM	0.5 ml
	Dimethyl Formamide	20%	10 ml
	Sterile water		38.5 ml

- Freezing buffer stock solution (FB stock; filter and store at -20°C):

	Reagent	Final Concentration	Per 50 ml
	1M Tris-HCl, pH 8.0	50 mM	2.5 ml
	1M Mg(OAc) ₂	5 mM	0.25 ml
	50% Glycerol	25%	25 ml
	0.5M EDTA	0.1 mM	0.01 ml
	Sterile water		22.24 ml

Buffers to make on the day of the experiment



■ RSB Lysis Buffer (200 µl per sample):

	Reagent	Final Concentration	Per 200 µl
	RSB	~1x	194 µl
	10% Igepal-CA630	0.1%	2 µl
	1% Digitonin	0.01%	2 µl
	10% Tween-20	0.1%	2 µl

■ RSB Washing Buffer (1 ml per sample):

	Reagent	Final Concentration	Per 1 ml
	RSB	~1x	990 µl
	10% Tween-20	0.1%	10 µl

■ Freezing buffer working solution (FBW; 1 ml for every 3 millions of nuclei):

	Reagent	Final Concentration	Per 1 ml
	FB stock	~1x	975 µl
	1M DTT	5 mM	5 µl
	Protease Inhibitors (Sigma P8340)	2% (v/v)	20 µl

txci-ATAC-seq protocol*Buffers to make beforehand*

■ TMG washing buffer (50 ml):

	Reagent	Final Concentration	Per 50 ml
	0.2M Tris-acetate pH 7.8	10 mM	2.5 ml
	1M Magnesium acetate	5 mM	0.25 ml
	50% Glycerol	10%	10 ml
	Sterile water		37.25 ml

■ Loading Buffer was made by mixing buffer 1 (5x) and buffer 2 below (5x):

1. Buffer1 (5x):

	Reagent	Final Concentration	Per 1000 μ l
	0.2M Tris-acetate pH 7.6	50 mM	250 μ l
	1M Magnesium acetate	25 mM	25 μ l
	Dimethyl Formamide	50%	500 μ l
	H2O		225 μ l

2. Buffer2 (5x):

	Reagent	Final Concentration	Per 1000 μ l
	100% Glycerol	50%	500 μ l
	5M NaCl	100 mM	20 μ l
	1M Tris-HCl, pH 7.5	50 mM	50 μ l
	0.5M EDTA	0.1 mM	0.2 μ l
	1M DTT	1 mM	1 μ l
	H2O		428.8 μ l

3. Loading Buffer:

	Reagent	Final concentration	Per 250 μ l
	5x Buffer1	1x	50 μ l
	5x Buffer2	1x	50 μ l
	H2O		150 μ l

■ Omni Resuspension Buffer (RSB; filter and store at 4C):

	Reagent	Final Concentration	Per 50 ml
	1M Tris-HCl, pH 7.5	10 mM	500 μ l
	5M NaCl	10 mM	100 μ l
	1M MgCl ₂	3 mM	150 μ l
	H2O		49.25 ml

■ 10% (100 mg/mL) BSA:

Dissolve 1 g powdered Fraction V or molecular biology grade BSA in 10 mL of distilled H₂O. Then, store in aliquots at -20°C.

Buffers to make on the day of the experiment

- PBSB (containing 0.04% BSA) (2 ml):

	Reagent	Final concentration	Per 2 ml
	PBS	~1x	1960 µl
	20 mg/ml BSA	0.4 mg/ml	40 µl

- RSB washing buffer (4 ml per sample):

	Reagent	Final Concentration	Per 4 ml
	RSB	~1x	3920 µl
	10% Tween-20	0.1%	40 µl
	10% BSA	0.1%	40 µl

- Loading buffer supplemented with SBS oligo (LBS):

	Reagent	Final Concentration	Per 150 µl
	Loading Buffer	~1x	140 µl
	75 µM SBS Oligo	5 µM	10 µl

Note: SBS oligo sequence (5'-3') is CGTGTGCTCTTCCGATCT

- 300 µM DAPI solution:

	Reagent	Final Concentration	Per 200 µl
	10.9 mM DAPI	300 µM	5.5 µl
	H2O		194.5 µl

Note: Add 1 µl of 300 uM DAPI to each 100 µl nuclei to make a final concentration of 3 µM for staining.

Loading Tn5

- 1 Resuspend Tn5ME-A, Tn5ME-B, and Tn5MErev in the annealing buffer to a final concentration of 100 μ M.
- 2 Prepare annealed linker A: Mix one volume of Tn5ME-A with one volume of Tn5MErev in a PCR tube.
e.g. 100 μ l Tn5ME-A + 100 μ l Tn5MErev.
- 3 Prepare annealed linker B: Mix one volume of each barcoded Tn5ME-B with one volume of Tn5MErev on a 96-well plate.
e.g. 10 μ l Tn5ME-B (Index A1) + 10 μ l Tn5MErev.



Table1_Barcoded_Tn5MEB.xlsx 12.8KB

- 4 Mix briefly by pipetting and run the following PCR program in a thermocycler for annealing oligos.

	Temperature	Time
	95 °C	5 min
	Slowly Cool down to 65 °C	-0.1 °C/sec
	65 °C	5 min
	Slowly Cool down to 4 °C	-0.1 °C/sec

The annealed oligos can be kept at -20°C for long-term storage.

- 5 Add 1 μ l of each annealed linker (A and B) to 20 μ l of the Tn5 stock (0.3 mg/ml) on a 96-well plate with a unique annealed linker B in each well.
- 6 Mix briefly by pipetting, and then incubate at 23°C for 30 minutes in a thermomixer at 350 rpm.
- 7 Store at -20°C.

Isolation of nuclei from cell lines

- 8 Remove approximately 10×10^6 cells from culture.



Note

The nuclei isolation protocol was adapted from Corces MR, et al. 2017.

Citation

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY (2017) . An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>

LINK

- 9 Pellet the cells at 500 RCF at 4°C for 5 min in a swinging-bucket centrifuge.
- 10 Aspirate supernatant.
- 11 Resuspend pellet in 200 µl RSB Lysis Buffer.
- 12 Incubate on ice for 3 minutes.
- 13 Add 1 ml RSB Washing Buffer.
- 14 Take 10 µl nuclei and dilute it with 40 µl of Omni TD buffer, then count the nuclei on a hemocytometer by adding 50 µl Trypan blue solution to the diluted nuclei (Note: we found that adding the RSB-resuspended nuclei straight to Trypan blue solution will cause inflation of nuclei, and diluting nuclei in Omni TD buffer before exposure to Trypan blue improves the nuclei integrity).



- 15 Pellet the remaining nuclei in RSB Washing Buffer at 500 RCF for 10 min at 4°C in a fixed-angle centrifuge.
- 16 Resuspend nuclei pellet in FBW at ~3 million nuclei/ml.
- 17 Snap-freeze nuclei in liquid nitrogen, and then transfer the cryovials to a liquid nitrogen dewar (or -80°C) for long-term storage.

Isolation of nuclei from lung tissue

- 18 The following protocol can be used to isolate nuclei from lung tissues.

Citation

Nikita Joshi, Alexander Misharin (2026)
. Single-nucleus isolation from frozen human lung tissue for single-nucleus RNA-seq.
protocols.io.

<https://protocols.io/view/single-nucleus-isolation-from-frozen-human-lung-ti-zu8f6zw>
LINK

txci-ATAC-seq: Preparing nuclei

- 19 Take out flash-frozen nuclei (~3 million in 1 ml) from liquid nitrogen for each sample and thaw in a water bath at 37°C for about 1 min.
- 20 Add 3 ml RSB washing buffer to an empty 15 ml tube for each sample.
- 21 Transfer 1 ml nuclei stored in the freezing buffer to the 15 ml tube containing 3 ml RSB washing buffer.
- 22 Pellet the nuclei at 500 RCF for 10 min at 4°C.



- 23 Resuspend nuclei with 1 ml RSB washing buffer and then transfer to a 1.5 ml LoBind tube through Flowmi (40 micron).
- 24 Pellet the nuclei at 500 RCF for 5 min at 4°C in a fixed-angle centrifuge.
- 25 Resuspend nuclei with 100 µl of 1X PBSB for each sample.
- 26 Count nuclei with DAPI:
 - 26.1 Add 1 µl of 300 µM DAPI to 100 ul nuclei;
 - 26.2 Incubate on ice for 5 mins;
 - 26.3 Add 10 µl stained nuclei to the countess slide to count nuclei.

txci-ATAC-seq: 96 barcoded Tn5 transposition

- 27 Prepare TD mix:

	Reagent	Final Concentration	Per Rxn	X120 rxn (in a 2ml tube)
	2X Nextera TD Buffer*	1X	12.5 µl	1500
	1% Digitonin	0.01%	0.25 µl	30
	10% Tween-20	0.1%	0.25 µl	30

*Omni TD buffer can be used to replace the Illumina Nextera TD buffer.



Citation

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY (2017)
. An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>

LINK

- 28 Thaw a 96-well plate preloaded with 5 μ l of barcoded Tn5 on ice. Mix by brief shaking at 1400 rpm for 30 seconds, spin for a minute at 2000 RCF at 4°C, and carefully unseal the aluminum foil seal.
- 29 Dilute nuclei to 2857 nuclei/ μ l in PBSB and then mix 7 μ l diluted nuclei with 13 μ l TD mix for each well.
- 30 Add 20 μ l nuclei/TD mix mixture to each well of the 96-well plate containing 5 μ l of barcoded Tn5 per well (total 25 μ l).
- 31 Seal the plate using Bio-Rad Microseal B film.
- 32 Mix by shaking at 1000 rpm for one minute.
- 33 If liquid splashes to the seal, briefly spin at 500 RCF for 10 sec.
- 34 Incubate at 37°C for 60 min in a thermocycler block with a heated lid (47°C).
- 35 Thaw TMG washing buffer on ice.
- 36 Remove the plate from the thermocycler.



- 37 Briefly centrifuge at 500 RCF for 10 sec at 4°C.
- 38 Incubate the plate on ice for 5 min.
- 39 Pool nuclei in a LoBind 12-tube strip and then transfer them to a 15 ml conical tube preloaded with 400 μ l of TMG.
- 40 Add 50 μ l/well of TMG to the first row of the plate and pipette them throughout the whole plate to wash out the residual nuclei remaining in the plate.
- 41 After washing the last row of the plate, the TMG was transferred to the same conical tube that was used to collect the barcoded nuclei.
- 42 Centrifuge nuclei at 500 RCF for 10 min at 4°C.
- 43 Remove most of the supernatant.
- 44 Resuspend nuclei with 500 μ l of TMG, then transfer to a 1.5 ml LoBind Tube through Flowmi.
- 45 Centrifuge at 500 RCF for 5 min at 4°C.
- 46 Remove most of the supernatant and resuspend the nuclear pellet with 30 μ l of LBS.
- 47 Count nuclei with a hemocytometer.
- 48 Take the volume of solution containing the desired number of nuclei and dilute it with the LBS to make a total of 15 μ l.
- 49 Use the 15 μ l dilution as input into the 10X Chromium droplet generator – follow Step 2, page 24 of the Chromium Single Cell ATAC kit instructions (**10x Document CG000209 Rev D**) to complete the assay.



txci-ATAC-seq: Modification of 10X Chromium protocol

50 For Step 2.5. GEM Incubation:

50.1 a. Incubate in a thermal cycler with the following protocol (Lid temperature at 105°C).

Temperature	Time
72 °C	00:05:00
98 °C	00:00:30
98 °C	00:00:10
59 °C	00:00:30
72 °C	00:01:00; Go to step 3, repeat 7X (Total 8 cycles)
15 °C	Hold

50.2 b. Store at 15°C for up to 18 h or at -20°C for up to 1 week, or proceed to the next step.

51 For Step 4.1 Sample Index PCR

51.1 c. Add 2.5 µl of customized i7 TruSeq primer (25 µM) containing an 8 bp custom barcode to each 10X library. Record assignment. Pipette mix and centrifuge briefly.

 Table2_TrueSeq_i7_Primer.xlsx 10.2KB

51.2 d. Incubate in a thermal cycler with the following protocol (Lid temperature at 105°C).

Temperature	Time
98 °C	00:00:45
98 °C	00:00:20
67 °C	00:00:30
72 °C	00:00:20; Go to step 2, repeat 4X (Total 5 cycles)



	Temperature	Time
	72 °C	00:01:00
	4 °C	Hold

Sequencing

52

	Sequencing Read	Cycles
	Read 1N	51 bp
	i7 index (I1)	8 bp
	i5 index (I2)	16 bp
	Read 2N	78 bp

Citations

Step 18

Nikita Joshi, Alexander Misharin. Single-nucleus isolation from frozen human lung tissue for single-nucleus RNA-seq

<https://protocols.io/view/single-nucleus-isolation-from-frozen-human-lung-ti-zu8f6zw>

Step 27

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY. An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>

Step 8

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ, Chang HY. An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues.

<https://doi.org/10.1038/nmeth.4396>