



Sep 08, 2025

Tethering and lipid transfer assay

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly4nkplx9/v1>

Elisabeth Holzer¹

¹Laboratory of Sascha Martens, Max Perutz Labs, University of Vienna, Austria



Elisabeth Holzer

Max Perutz Labs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly4nkplx9/v1>

Protocol Citation: Elisabeth Holzer 2025. Tethering and lipid transfer assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly4nkplx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 01, 2025

Last Modified: September 08, 2025

Protocol Integer ID: 223995

Keywords: ASAPCRN, process of liposome tethering, liposome tethering, lipid transfer, tethering

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Abstract

This protocol describes the process of liposome tethering and lipid transfer.

Materials

Reaction Mix

	A	B
	Donor liposomes (of 0.5 mg/mL total lipids)	2.5 µL
	Acceptor liposomes (of 0.5 mg/mL total lipids)	2.5 µL
	Proteins (final concentrations as per figure legends):	
	ATG2A	250 nM
	WIPI4	100 nM
	LC3B	100 nM

Equipment

Glass-bottom 384-well microplates

NAME

SENSOPLATE

BRAND

781892

SKU

<https://shop.gbo.com/en/india/products/bioscience/microplates/sensoplate-glass-bottom-plates/384-well-sensoplate/781892.html>

LIN
K



Plate Setup and Reaction Conditions

- 1 Perform reactions in glass-bottom 384-well microplates (Greiner Bio-One, Cat# 781892).
- 2 Carry out all assays at  Room temperature . 

Reaction Mix

- 3 Each  50 μL reaction (in SEC300 buffer) includes: 

A	B
Donor liposomes (of 0.5 mg/mL total lipids)	2.5 μL
Acceptor liposomes (of 0.5 mg/mL total lipids)	2.5 μL
Proteins (final concentrations as per figure legends):	
ATG2A	250 nM
WIPI4	100 nM
LC3B	100 nM

NBD Fluorescence Measurement (Lipid Transfer)

2h 20m 20s

- 4 Use a Tecan SPARK Multimode Microplate Reader (RRID:SCR_021897).
- 5 Record NBD fluorescence:
 - Excitation: 485 nm
 - Emission: 535 nm
 - Measurement interval: every  00:00:20
 - Duration:  01:30:00
- 6 Add 0.5% (v/v) DDM to each well to fully solubilize lipids. 



-  2.5 μL of 10% (v/v) DDM to  50 μL reactions.

7 Measure maximum NBD fluorescence  00:50:00 after DDM addition.

50m

8 Normalize fluorescence:

- Subtract the initial fluorescence value.
- Divide by the maximal DDM-treated fluorescence.

Liposome Tethering (Optional Parallel Readout)

9 Measure absorbance at 405 nm every  00:00:20 in the same wells.



10 This quantifies vesicle clustering (tethering) under identical conditions using the same plate reader.