

Jun 30, 2025

TH-TdTomato expression measurement by FACS

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvjbk7xgk5/v1>

Elena Coccia¹

¹Icahn School of Medicine at Mount Sinai



Elena Coccia

Icahn School of Medicine at Mount Sinai

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvjbk7xgk5/v1>

Protocol Citation: Elena Coccia 2025. TH-TdTomato expression measurement by FACS. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvjbk7xgk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 18, 2025



Last Modified: June 30, 2025

Protocol Integer ID: 220539

Keywords: ASAPCRN, tdtomato expression measurement by fac, tdtomato expression measurement, analysis of midbrain neuron, midbrain neuron, cytometry, td, neuron

Abstract

This protocol outlines the preparation and flow cytometry-based analysis of midbrain neurons differentiated from Td-Tomato reporter hPSCs.



Cell Culture and Preparation

- 1 Differentiate TH-TdTomato lines into midbrain neurons in **6-well plates** until the desired **differentiation day**.
- 2 Confirm cell morphology and fluorescence (if applicable) before harvesting.

Cell Detachment

- 3 Remove media and wash cells once with PBS
- 4 Gently lift neurons using Accutase supplemented with ROCK inhibitor (Y-27632). 1ml per well of a 6w plate
- 4.1 Incubate as needed (typically ~5 minutes) at RT until cells detach.
- 5 Collect cells into tubes with twice the volume of culture medium containing ROCK inhibitor.
- 6 Centrifuge the cell suspension at $300 \times g$ for 5 minutes at room temperature.
- 7 Carefully remove the supernatant.
- 8 Resuspend the pellet in FACS buffer: PBS + 1x B-27 supplement + ROCK inhibitor (Y-27632) + DAPI (final concentration $1\mu\text{g/ml}$) for live/dead discrimination
- 8.1 Keep cells on ice from this step forward
- 8.2 Make sure to include a sample with no DAPI for proper gating
- 9 For gating: collect TdTomato-negative line, or undifferentiated iPSC as a negative control to establish the TdTomato gate.





- 10 Pass the resuspended cells through a 35 μ m cell strainer into FACS-compatible tubes

Flow cytometry (FACS)

- 11 Analyze TdTomato fluorescence using a flow cytometer with appropriate excitation/emission settings (e.g., 554/581 nm) and DAPI (Ex/Em ~358/461 nm).
 - 11.1 Record at least 10,000 events from live population
- 12 Gating Strategy:
 - 12.1 FSC/SSC: Gate to exclude debris and cell clumps.
 - 12.2 DAPI-negative cells: Gate for live cells excluding DAPI+ dead cells
 - 12.3 TdTomato fluorescence:
Set gate using Td-tomato negative line
Identify TdTomato+ neurons and MFI from live population.
- 13 Analysis with proper software (FlowJo)

