



Th9 Polarization of Mouse CD4+ Cells V.3

DOI

dx.doi.org/10.17504/protocols.io.798hr9w



Sam Li¹

¹BioLegend

BioLegend

Tech. support email: tech@biolegend.com



Sam Li

BioLegend

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.798hr9w>

External link: <https://www.biolegend.com/protocols/th9-polarization-of-mouse-cd4--cells-protocol/4243/>

Protocol Citation: Sam Li . Th9 Polarization of Mouse CD4+ Cells. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.798hr9w>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: October 15, 2019



Last Modified: October 15, 2019

Protocol Integer ID: 28704

Keywords: Th9, polarization, polarization of mouse cd4, mouse cd4, cell, mouse

Materials

MATERIALS

✕ MojoSort™ Mouse CD4 T Cell Isolation Kit **BioLegend Catalog #480005, 480006, 480033**

✕ RBC Lysis Buffer **BioLegend Catalog #420301**

✕ Anti-mouse CD3ε clone 145-2C11 (Ultra-LEAF™ format) **BioLegend Catalog #100339**

✕ Anti-mouse CD28 clone 37.51 (Ultra-LEAF™ format) **BioLegend Catalog #102116**

✕ Anti-mouse IFN-γ clone XMG1.2 (Ultra-LEAF™ format) **BioLegend Catalog #505834**

✕ Recombinant human TGF-β1 (carrier-free) **BioLegend Catalog #580702**

✕ Monensin Solution **BioLegend Catalog #420701**

✕ Recombinant mouse IL-2 (carrier-free) **BioLegend Catalog #575402**

✕ Recombinant mouse IL-4 (carrier-free) **BioLegend Catalog #574302**

- Sterile PBS
- Cell culture medium (IMDM supplemented with 10% FBS)
- Sterile plastic petri dishes
- PMA (Phorbol 12-myristate 13-acetate) (Cat. No. P8139 from Sigma)
- Ionomycin (Cat. No. I0634 from Sigma)



Troubleshooting



Isolation of CD4+ Cells From Lymph Nodes

- 1 Harvest lymph nodes (superficial cervical, mandibular, axillary, inguinal, and mesenteric) from mice.
- 2 Tease lymph nodes through a sterile 70- μ m nylon cell strainer to obtain single-cell suspensions in complete IMDM containing 10% FCS (complete medium).
- 3 Resuspend cells in complete medium and use your favorite method to isolate CD4⁺ cells. Consider using our [**MojoSort™ Mouse CD4 T Cell Isolation Kit**](#).

Th9 Polarization of CD4+ Cells:

- 4 On day 0, coat 60 × 15mm of plastic petri dishes with anti-mouse CD3 ϵ , clone 145-2C11 (5 μ g/ml). Incubate at 37°C for 2 hours or 4°C overnight. Aseptically decant antibody solution from the plate. Wash plate 3 times with sterile PBS. Discard liquid.
 02:00:00
- 5 Plate CD4⁺ cells at 10×10^6 /5 ml/dish. Culture cells for 3 days in the presence of anti-mouse CD28, clone 37.51 (5 μ g/mL), recombinant human TGF- β 1 (10ng/mL), recombinant mouse IL-4 (10ng/mL), recombinant mouse IL-2 (20ng/ml), and anti-mouse IFN- γ , clone XMG1.2 (10 μ g/mL).
- 6 On day 3, wash cells once and then restimulate in complete medium with 500ng/ml PMA and 500ng/mL ionomycin, in the presence of monensin for 6 hours.  06:00:00
- 7 After harvesting, the cells are ready for staining.

Note: Recombinant human TGF- β is effective for stimulating mouse cells.