


Feb 26, 2019

Version 1

Thawing and Seeding Frozen Cells V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.xkmfku6>



Kenneth Schackart¹

¹University of Arizona

Yoon Lab



Kenneth Schackart

University of Arizona

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xkmfku6>

Protocol Citation: Kenneth Schackart 2019. Thawing and Seeding Frozen Cells. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.xkmfku6>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 30, 2019

Last Modified: February 26, 2019

Protocol Integer ID: 19821

Keywords: thaw cells from the liquid nitrogen storage, seeding frozen cell, thaw cell, tissue culture flask, liquid nitrogen storage, cell, tissue

Abstract

How to thaw cells from the liquid nitrogen storage and seed into a tissue culture flask

Guidelines

Wear gloves. Spray anything entering or exiting the biosafety cabinet with 70% Ethanol.

Materials

(1) T-75 or T-25 flask per frozen cell vial (or more if plating at a lower density)

(1) 15 mL centrifuge tube per frozen cell vial

(1) 10 mL serological pipette tip per T-75 flask or (1) 5 mL serological pipette tip per T-25 flask

Warmed cell culture media

1000 µL filter pipette tips

Before start

UV appropriate number of flasks, 15 mL centrifuge tubes, waste beaker, and serological pipette tips.

Warm appropriate cell culture media.



Thaw Cells

- 1 Thaw cells by suspending in  37 °C water bath until completely thawed, but no longer than necessary

Transfer cell suspension

- 2 Within biosafety cabinet, transfer cell suspension to 15 mL centrifuge tube using 1000 µL pipette.

Dilute freezing medium

- 3 Add  1 mL warmed cell culture medium to cell suspension *dropwise*.

Note

Adding the initial cell culture medium slowly helps prevent cell death caused by a rapid change in osmotic pressure.

- 4 Add an additional  3 mL warmed cell culture medium to cell suspension slowly.

Centrifuge cell suspension

- 5 Centrifuge the cell suspension at 1.5 kRPM for  00:03:00 .

Resuspend Cells

- 6 Remove bulk of supernatant with serological pipette, then remove remainder with 1000 µL pipette.

**Note**

For small cell pellets, you are better off leaving a small amount of media than disturbing the cell pellet.

- 7 Add 1 mL warmed cell culture media to cell pellet.

Note

Allowing the cell pellet rest in media for about 2 minutes will help with resuspension.

- 8 Gently pipette mix the cell pellet into the solution.

- 9 Add an additional 7 mL warmed cell culture media if using a T-75 flask.
Add an additional 3 mL warmed cell culture media if using a T-25 flask.

Seed Cells

- 10 If using a T-75 flask, first add  2 mL warmed cell culture media to the flask.
Using a serological pipette, transfer the cell suspension to the tissue culture flask.

Label Flask

- 11 Label the flask with:
- Cell line
 - Passage number
 - Date
 - Your initials

Incubate

- 12 Transfer flask to CO₂ incubator.



Documentation

- 13 Don't forget to remove the vial you used from the frozen storage inventory.