

Sep 25, 2019

Version 1

# The 'Three Peaks' faecal DNA extraction method for long-read sequencing V.1

DOI

[dx.doi.org/10.17504/protocols.io.584g9yw](https://dx.doi.org/10.17504/protocols.io.584g9yw)

Josh Quick<sup>1</sup>

<sup>1</sup>Sam Nicholls [University of Birmingham], Nicholas Loman [University of Birmingham]

High molecular weight D...

Long Read Club



Josh Quick

University of Birmingham

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.584g9yw>

External link: <https://www.slideshare.net/scalene/the-three-peak-challenge-for-longread-ultradeep-stool-metagenomics-on-the-promethion>

**Protocol Citation:** Josh Quick 2019. The 'Three Peaks' faecal DNA extraction method for long-read sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.584g9yw>

**Manuscript citation:**<https://doi.org/10.1093/gigascience/giz043>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** August 07, 2019

**Last Modified:** September 25, 2019

**Protocol Integer ID:** 26620

**Keywords:** Microbiome research, metagenomics, long-read sequencing, nanopore sequencing, genomics, DNA extraction, sample preparation, faecal dna extraction method, dna

## Materials

### MATERIALS

☒ MetaPolyzyme Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG

### STEP MATERIALS

🧴 500 µL PBS pH 7.5

☒ DNA/RNA Shield Zymo Research Catalog #R1100-50

☒ MetaPolyzyme Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG

Resuspend the contents of the bottle in 🧴 500 µL PBS pH 7.5 and aliquot for optimal activity

## Protocol materials

☒ MetaPolyzyme Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG

☒ DNA/RNA Shield Zymo Research Catalog #R1100-50

☒ MetaPolyzyme Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG

☒ OMNIgene•GUT DNA Genotek Catalog #OMR-200

☒ DNA/RNA Shield Zymo Research Catalog #R1100-50

☒ MetaPolyzyme Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG

## Troubleshooting



- 1 Add  100 mg fresh stool or  50  $\mu\text{L}$  OMNIgene GUT to  200  $\mu\text{L}$  DNA/RNA Shield, vortex briefly and place on a tube rotator at  20 rpm  00:10:00 .

 OMNIgene•GUT **DNA Genotek Catalog #OMR-200**

 DNA/RNA Shield **Zymo Research Catalog #R1100-50**

- 2 Centrifuge at  5000 x g  00:05:00 and retain up to  200  $\mu\text{L}$  supernatant depending on size of pellet.

- 3 Add  100  $\mu\text{L}$  PBS and resuspend material by pipetting up and down, centrifuge  5000 x g  00:05:00 and retain up to  100  $\mu\text{L}$  supernatant depending on size of pellet.

- 4 Add  1 mL PBS and resuspend material by pipetting up and down, centrifuge  5000 x g  00:05:00 and discard supernatant.



- 5 Add  100  $\mu\text{L}$  PBS and  5  $\mu\text{L}$  MetaPolyzyme, mix by pipetting and incubate at  35  $^{\circ}\text{C}$   02:00:00 . Gently mix by pipetting up and down to ensure solution is homogeneous before proceeding.

 MetaPolyzyme **Merck MilliporeSigma (Sigma-Aldrich) Catalog #MAC4L-5MG**

- 6 Add  100  $\mu\text{L}$  DNA/RNA Shield,  10  $\mu\text{L}$  [M] 10 % (w/v) SDS and  10  $\mu\text{L}$  [M] 20 mg/mL Proteinase K and mix by pipetting. Incubate at  55  $^{\circ}\text{C}$   00:30:00 with  300 rpm mixing. Gently mix by pipetting up and down to ensure solution is homogeneous before proceeding.

- 7 Centrifuge  5000 x g  00:05:00 and retain up to  200  $\mu\text{L}$  supernatant depending on size of pellet.

- 8 Resuspend pellet in  750  $\mu\text{L}$  Lysis Solution and transfer to a ZR BashingBead Lysis Tube. Bead-beat on a FastPrep instrument for 1 cycle of  00:00:40 at  6 m/s



## Equipment

FastPrep-24 5G

NAME

Bead-beater

TYPE

MP Biomedicals

BRAND

116005500

SKU

<https://uk.mpbio.com/fastprep-24-5g-instrument>

LINK



9 Centrifuge tube at  10000 x g for  00:01:00 and retain  400  $\mu$ L supernatant.

10 Pool the supernatants retained at each of the steps in a 2 ml Eppendorf tube and measure the volume using a P1000 pipette. Add 2 volumes of Genomic Lysis Buffer and  2 mL MagBinding beads.