

Feb 13, 2026

The establishment of an in vitro immunotherapy model based on MTSs

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1kxnygr2/v1>

Guangli Suo¹, Zhe Zhao¹, Mengdie Yang¹, Xuan Xiong¹

¹CAS Key Laboratory of Nano-Bio Interface, Suzhou Institute of Nano-Tech and Nano-Bionics, Chinese Academy of Sciences, Jiangsu, 215123, China.

MTS



Guangli Suo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1kxnygr2/v1>

Protocol Citation: Guangli Suo, Zhe Zhao, Mengdie Yang, Xuan Xiong 2026. The establishment of an in vitro immunotherapy model based on MTSs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1kxnygr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 13, 2026

Last Modified: February 13, 2026

Protocol Integer ID: 243169

Keywords: immunotherapy model, characterization of human pbmc, mtss isolation, human pbmc, mt

Funders Acknowledgements:

National Natural Science Foundation of China

Grant ID: 32470814

Abstract

Isolation and characterization of human PBMCs or CD8+ T lymphocytes, followed by their co-culture with MTS in a MWC-chip.

The isolation of PBMC

- 1 The peripheral blood mononuclear cells (PBMCs) were isolated from a healthy donor's whole blood sample using Ficoll-Paque density gradient medium (GE Healthcare) according to the manufacturer's instructions. Briefly, the whole blood was diluted with DPBS and layered onto Ficoll-Paque medium, followed by centrifugation at 400 g for 30 minutes. After centrifugation, the mononuclear cells located at the interface were transferred to a new tube, diluted with DPBS, and pelleted at 400 g for 10 minutes. The cells were washed once with DPBS before downstream applications. Human cytotoxic T lymphocytes (CTLs) were isolated from PBMCs using the EasySep Human CD8⁺ T Cell Enrichment Kit (STEMCELL Technologies) following the manufacturer's instructions. Briefly, PBMCs derived from BRC patients were incubated with 50 μ L/mL of Enrichment Cocktail for 10 minutes at room temperature. Subsequently, the sample was incubated with magnetic particles at a concentration of 150 μ L/mL for 5 minutes at room temperature. The sample was then diluted with EasySep Buffer (STEMCELL Technologies) and subjected to cell separation using an EasySep Magnet (STEMCELL Technologies). CTLs, or CD8⁺ T cells were recovered by transferring the cell suspension into a new tube and centrifuging at 300 g for 5 minutes, and then cultured in T cell medium (STEMCELL Technologies).

The establishment of the co-culture system with MTSs and PBMC

- 2 To establish the co-culture system with MTSs and PBMC, the collected MTSs were seeded into a new MWC-chip, then PBMC was stained with DID and seeded in the same MWC-chip at the ratio of PBMC / MTSs = 5 / 1. After that, 1.5 mL SG-Medium mixed with Annexin V-FITC (Vazyme) was added in the co-cultural system and maintained at 37°C in 5% CO₂ for 48 hours and imaged by Confocal (OLYMPUS). To establish the co-culture system with MTSs and CTLs (CD8⁺ T Cell), the sorted MTSs and CTLs were mixed with the ratio of CTLs / MTSs = 5 / 1 and seeded in a new MWC-chip. Then, 1.5 mL SG-Medium mixed with Annexin V-FITC was added in the co-cultural system and maintained at 37°C in 5% CO₂ for 48 hours and imaged by Confocal (OLYMPUS). Before assessing the anti-cancer efficacy of Nivolumab (PD-1 antibody, an immune checkpoint inhibitor) on MTSs, PD-L1 expression in MTSs was assessed by staining them with a PD-L1 antibody on a MWC-chip. In the CTLs + MTSs co-culture system for immune cell potency assay, Nivolumab and annexin V-FITC were added to the medium, followed by incubation at 37°C in 5% CO₂ for 48 hours. Confocal imaging (OLYMPUS) was performed to visualize the results. Control groups included MWC-chips seeded with MTSs or MTS + CTLs without Nivolumab treatment.