

Jun 13, 2025

Version 2

Thick 3D Tetraspeck Beads sample preparation V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.n92ldnkkov5b/v2>



Elena Carlson¹

¹Fred Hutchinson Cancer Center



Elena Carlson

Fred Hutchinson Cancer Center

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldnkkov5b/v2>

Protocol Citation: Elena Carlson 2025. Thick 3D Tetraspeck Beads sample preparation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldnkkov5b/v2> Version created by **Elena Carlson**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 13, 2025

Last Modified: June 13, 2025

Protocol Integer ID: 220162

Keywords: tetraspeck, fluoresence, microscopy, point spread function, agarose, beads, fluorescent beads, thick 3d tetraspeck beads sample preparation, fluorescent bead, sample preparation, microscopy troubleshooting, thick sample, thick sample of agar

Funders Acknowledgements:

National Institutes of Health

Grant ID: P30CA015704 (RRID:SCR_022616)

National Institutes of Health

Grant ID: P30CA015704 (RRID:SCR_022609)

National Institutes of Health

Grant ID: P30CA015704 (RRID:SCR_022610)

Abstract

Create a thick sample of agar and fluorescent beads in a well plate, to use for microscopy troubleshooting (e.g. to determine the point spread function of your system).

Materials

- Tetraspeck beads, 0.5 um, ThermoFisher Scientific cat. no: T7281. TetraSpeck™ Microspheres, 0.5 μm, fluorescent blue/green/orange/dark red
- Agarose HOEFER | SKU: SKU:GR140-25
- P1000 and P200 micropipettes + tips
- vortexer (any)
- rice cooker or microwave
- Thermometer
- precision scale
- 15 mL falcon tubes
- glass bottom well plate, e.g. CLS-1811 CLS-1811 - GLASS BOTTOM CULTURE DISHES, TC TREATED, NEST-Chemglass Life Sciences



Prep 1% [w/v] Agarose solution

- 1 Heat DPBS to approximately $100\text{ }^{\circ}\text{C}$, using a rice cooker, pressure cooker, or microwave
use a thermometer to verify temperature
- 2 Measure out 104 mg agar to 10.4 mL hot DPBS and mix vigorously with vortexer. If using a larger well plate, you can adjust these amounts accordingly. Verify concentration made using this website [Percent \(%\) Solutions Calculator - PhysiologyWeb](#) .

Mix in the beads

- 3 Let Agar solution cool to $40\text{-}50\text{ }^{\circ}\text{C}$
Meanwhile, vortex the stock tetraspeck vial

Equipment

Vortex mixer

NAME

Any

BRAND

xx


SKU

- 3.1 When DPBS-Agar mix has cooled to $40\text{-}50\text{ }^{\circ}\text{C}$, use a P200 micropipette to add $60\text{ }\mu\text{L}$ of the stock beads to the 10 mL of hot Agar-DPBS. For a sparser population of beads, add 30 microliters instead.

Vortex the resulting mixture vigorously.

Well plate

- 4 Use a simple transfer pipette to transfer the beads-agarose solution to a glassbottom well plate. Fill the plate to the very top, avoiding the formation of bubbles. Before completely filling the plate, tap the bottom of the plate gently to release any bubbles.
- 5 Carefully place a coverslip on top of the well plate, and let the agarose cool and harden. Shield the sample from light.



Agarose-Beads suspension in a thick well plate, with a glass coverslip on top. This is then stored in a larger well plate dish.



Store the sample in a larger well plate or similar box to keep the glass from being knocked off. Add tape.

Protocol references

Method adapted from [GitHub - amsikking/fluorescent_beads_in_agarose: How to make a 3D sample of beads for PSF measurements with ~water refractive index \(1.33\)](#)

Above was adapted from:

Yiming Li, Yu-Le Wu, Philipp Hoess, Markus Mund, and Jonas Ries, "Depth-dependent PSF calibration and aberration correction for 3D single-molecule localization," Biomed. Opt. Express **10**, 2708-2718 (2019)

Acknowledgements

This research was supported by the Preclinical Imaging, Comparative Medicine, and Cellular Imaging Shared Resources ([RRID:SCR_022616](#) , [RRID:SCR_022609](#), [RRID:SCR_022610](#)) of the Fred Hutch/University of Washington/Seattle Children's Cancer Consortium (P30 CA015704).