



Dec 06, 2024

## Thymidine analog labeling of proliferating cells

DOI

[dx.doi.org/10.17504/protocols.io.261gerdmol47/v1](https://dx.doi.org/10.17504/protocols.io.261gerdmol47/v1)

Ashley Seifert<sup>1</sup>

<sup>1</sup>University of Kentucky



Ashley Seifert

University of Kentucky

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gerdmol47/v1>

**Protocol Citation:** Ashley Seifert 2024. Thymidine analog labeling of proliferating cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.261gerdmol47/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** December 06, 2024

**Last Modified:** December 06, 2024

**Protocol Integer ID:** 114530

**Keywords:** ASAPCRN, cycling cells with thymidine analog, thymidine analog labeling, thymidine analog, assessing cell proliferation, cycling cells before the injury, proliferating cell, cell proliferation, cells this protocol, cell proliferation in response, cycling cell, seifert lab, label, injured tissue, cell, labeling, lab, institutional animal care

## Abstract


This protocol describes the procedure used by the Seifert Lab to label cycling cells with thymidine analogs to track cell proliferation in injured tissue. EdU is used to label cycling cells before the injury, and BrdU after, thereby assessing cell proliferation in response to injury. It has been successfully used in rodents (*Mus musculus* and *Acomys dimidiatus* - spiny mouse), zebrafish (*Danio rerio*) and axlotl (*Ambystoma mexicanum*).


Note that any protocol involving animals should be reviewed and approved by your Institutional Animal Care and Use Committee (IACUC) before use.

## Materials

1. Thymidine analogs diluted in DMSO (EdU and BrdU)
2. Sterile saline solution (0.9% NaCl)
3. 1.5 mL Eppendorf tubes
4. Syringe (2 mL capacity, suitable gauge for injection)
5. Appropriate injection needles (e.g., 25-30 gauge for intraperitoneal injection)
6. Animals, e.g. *Acomys dimidiatus* (spiny mice), *Mus musculus*, *Ambystoma mexicanum* (axolotls), *Danio rerio* (zebrafish)








## Protocol materials

 BrdU (5-Bromo-2'-Deoxyuridine) **Thermo Fisher Catalog #B23151**

 EdU (5-ethynyl-2'-deoxyuridine) **Thermo Fisher Catalog #A10044**

## Troubleshooting



- 1 *For Acomys, Mus, or Ambystoma:*  
Dilute  EdU (5-ethynyl-2'-deoxyuridine) **Thermo Fisher Catalog #A10044**  
from stock in sterile saline to a final concentration of  10 µg/g body weight
- 2 *For Acomys, Mus, or Ambystoma:*  
Dilute  BrdU (5-Bromo-2'-Deoxyuridine) **Thermo Fisher Catalog #B23151**  
from stock in sterile saline to a final concentration of  1 µg/g body weight
- 3 *For Acomys, Mus, or Ambystoma:*  
Transfer the thymidine analog solution to a 2 mL syringe appropriate for intraperitoneal injection (25-30 gauge)
- 4 *For Acomys, Mus, or Ambystoma:*  
Inject EdU intraperitoneally three times at 3-hour intervals (i.e., 0h, 3h, and 6h)
- 5 *For Acomys, Mus, or Ambystoma:*  
Collect tissue at 9h (for assessment of baseline cell proliferation)
- 6 *For Acomys, Mus, or Ambystoma:*  
At desired time following injury, inject BrdU intraperitoneally three times at 3-hour intervals (i.e., 0h, 3h, and 6h)
- 7 *For Acomys, Mus, or Ambystoma:*  
Collect tissue at 9h (for assessment of injury-induced cell proliferation)
- 8 *For Danio:*  
Prepare a final solution of  0.4 millimolar (mM) EdU or BrdU in zebrafish water 
- 9 *For Danio:*  
Place zebrafish in the prepared water for a  09:00:00 exposure period to allow analog uptake before collecting tissue 