

Feb 25, 2019

Tissue lysis and digestion for MS analysis

DOI

dx.doi.org/10.17504/protocols.io.yguftww

David Kotol¹

¹Department of Protein Science, Royal Institute of Technology, Stockholm, Sweden

Human Cell Atlas Metho...

Human Protein Atlas



David Kotol

Science fo Life Laboratory / KTH - Royal Institute of Techno...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yguftww

Protocol Citation: David Kotol 2019. Tissue lysis and digestion for MS analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yguftww>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 23, 2019

Last Modified: February 25, 2019

Protocol Integer ID: 20724

Abstract

Here, a workflow for sample preparation and quantification of brain and pancreatic tissue proteins with use of heavy labelled protein standards (QPrESTs) is described. The standards are added before the digestion ensuring high reproducibility and elimination of nuances appearing during proteolytical cleavage of proteins and further steps such as solid phase extraction. QPrESTs are heavy labelled recombinant protein fragments covering regions of proteins with low homology and span over 50-150 aminoacids. They were produced as a part of Human Protein Atlas project in auxotrophic E.coli and individually purified, serving as an excellent resource for MS-based protein quantification.



Materials

MATERIALS

☒ Methanol (MeOH) **Fisher Scientific Catalog #10499560**

☒ Acetonitrile (ACN) **Fischer Scientific Catalog #10660131**

☒ Sodium deoxycholate (SDC) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #30970**

☒ Urea **Merck MilliporeSigma (Sigma-Aldrich) Catalog #U5378**

☒ Thiourea **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8656**

☒ DL-Dithiothreitol (DTT) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #43815**

☒ 2-Chloroacetamide (CAA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #22790**

☒ EDTA-free Protease Inhibitor Cocktail **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11836170001**

☒ Acetone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34850**

☒ Pierce Trypsin Protease **Thermo Fisher Scientific Catalog #90057**

☒ Formic acid (FA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #15657520**

☒ Trifluoroacetic acid (TFA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #74564**

☒ Triethylammonium bicarbonate (TEAB) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T7408**

☒ Bio-Rad Protein Assay **Bio-Rad Laboratories Catalog #5000001**

☒ Ammonium hydroxide **Alfa Aesar Catalog #L13168**

STEP MATERIALS

☒ Bio-Rad Protein Assay **Bio-Rad Laboratories Catalog #5000001**

Buffers were prepared using MilliQ water if not stated differently

Tissue lysis and digestion

Lysis buffer, 1mL (7M Urea, 2M Thiourea, 2% SDC)

m (Urea) = 0.42 g

m (Thiourea) = 0.15 g

m (SDC) = 0.02 g

1/5 tablet of protease inhibitors

300mM DTT, 500 µL

m (DTT) = 23 mg

500mM CAA, 1 mL

m (CAA) = 46.75 mg

Dissolved in 100mM TEAB

100mM TEAB, 1 mL

100 µL 1M TEAB

Solvent A, 100 mL (3% ACN, 0.1% FA)

3 mL 100% ACN

100 µL 100% FA

Solid Phase Extraction using Strong Cation Exchange

Wash Buffer, 1 mL (30% MeOH, 0.1% FA)

300 µL 100% MeOH

1 µL 100% FA

Elution Buffer, 1 mL (30% MeOH, 1.25% NH₄OH)

300 µL 100% MeOH

45 µL 28% NH₄OH

Protocol materials

⊗ 2-Chloroacetamide (CAA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #22790

⊗ Bio-Rad Protein Assay Bio-Rad Laboratories Catalog #5000001

⊗ Urea Merck MilliporeSigma (Sigma-Aldrich) Catalog #U5378

⊗ Bio-Rad Protein Assay Bio-Rad Laboratories Catalog #5000001

⊗ Acetone Merck MilliporeSigma (Sigma-Aldrich) Catalog #34850

⊗ Triethylammonium bicarbonate (TEAB) Merck MilliporeSigma (Sigma-Aldrich) Catalog #T7408

⊗ Ammonium hydroxide Alfa Aesar Catalog #L13168

⊗ Thiourea Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8656

⊗ DL-Dithiothreitol (DTT) Merck MilliporeSigma (Sigma-Aldrich) Catalog #43815

⊗ EDTA-free Protease Inhibitor Cocktail Merck MilliporeSigma (Sigma-Aldrich) Catalog #11836170001

⊗ Acetonitrile (ACN) Fischer Scientific Catalog #10660131

⊗ Formic acid (FA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #15657520

⊗ Sodium deoxycholate (SDC) Merck MilliporeSigma (Sigma-Aldrich) Catalog #30970

⊗ Pierce Trypsin Protease Thermo Fisher Scientific Catalog #90057

⊗ Methanol (MeOH) Fisher Scientific Catalog #10499560


⊗ Trifluoroacetic acid (TFA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #74564

⊗ Bio-Rad Protein Assay Bio-Rad Laboratories Catalog #5000001

Safety warnings

- ❗ Chloroacetamide (2-chloroacetamide) is a chlorinated organic compound used for alkylating reduced cysteine residues and is suspected of damaging fertility.

Protein extraction

- 1 Incubate adapter of the Tissue Lyser on dry ice for  00:30:00

Equipment

Tissue Lyser LT

NAME

QIAGEN

BRAND

85600

SKU

<https://www.qiagen.com/us/shop/automated-solutions/sample-disruption/tissuelyser-lt/#orderinginformation>

LINK

Bead mill

SPECIFICATIONS

- 2 Add one 3mm bead to the tissue sample and incubate on dry ice for  00:30:00

Equipment

Tungsten Carbide Beads, 3 mm

NAME

QIAGEN

BRAND

69997


SKU

<https://www.qiagen.com/se/shop/automated-solutions/accessories/beads/#orderinginformation>

LINK

Beads for use with Tissue Lyser systems

SPECIFICATIONS

3 Disrupt the tissue using the Tissue Lyser for  00:02:00

Equipment

Tissue Lyser LT

NAME

QIAGEN

BRAND

85600


SKU


<https://www.qiagen.com/us/shop/automated-solutions/sample-disruption/tissuelyser-lt/#orderinginformation>

LINK

Bead mill

SPECIFICATIONS

4 Add  250 μ L of Lysis Buffer (7M Urea, 2M Thiourea, 2% SDC)

5 Homogenize the tissue using the Tissue Lyser for  00:02:00



Equipment

Tissue Lyser LT

NAME

QIAGEN

BRAND

85600










SKU

<https://www.qiagen.com/us/shop/automated-solutions/sample-disruption/tissue-lyser-lt/#orderinginformation>



LINK

Bead mill

SPECIFICATIONS









- 6 Centrifuge for  00:30:00 at 20,000 × g at  4 °C
- 7 Precipitate the proteins with 4 volumes of ice-cold acetone and incubate in  -20 °C over night
- 8 Centrifuge at 20,000 × g at  4 °C for  00:30:00 and discard the supernatant
- 9 Wash the pellet two times with  400 µL of ice-cold acetone followed by centrifugation at 20,000 × g at  4 °C for  00:10:00  00:10:00
- 10 Remove the acetone excess and dry the pellet



- 11 Dissolve the pellet in  200 μL of Pellet Solubilizer (7M Urea, 2M Thiourea)
- 12 Dilute 8 times with  100 millimolar (mM) TEAB and measure protein concentration using Bio-Rad Protein assay according to the manufacturers instructions

 Bio-Rad Protein Assay **Bio-Rad Laboratories Catalog #5000001**

Protein Digestion

- 13 Take volume that corresponds to  15 μg of proteins and add QPrESTs mixture
- 14 Add DTT to a final concentration of  10 millimolar (mM) and incubate  01:00:00 at  30 $^{\circ}\text{C}$
- 15 Add CAA to a final concentration of  50 millimolar (mM) and incubate  00:30:00 in the dark at room temperature
- 16 Add  300 ng of Pierce Trypsin protease and incubate over night at  37 $^{\circ}\text{C}$
- 17 Quench the digestion with TFA to a final concentration of 0.5% (v/v)

Solid phase extraction using Strong Cation Exchange

- 18 Prepare one StageTip per sample by inserting 3 layers of Empore Cation Extraction Discs in a 200 μL tip



Equipment

Empore Discs

NAME

47mm cation extraction disks

TYPE















3M

BRAND

2251-1

SKU

<https://www.sterlitech.com/empore-disks-cation-exchange-sr-2251-1.html>^{LINK}

- 19 Activate membrane with  50 μL 100% MeOH, spin at 1000 x g for  00:01:00
- 20 Equilibrate StageTip with  50 μL Wash Buffer (30% MeOH, 0.1% FA), spin at 1000 x g for  00:01:00
- 21 Apply sample to the StageTip, spin at 1000 x g for  00:01:00
- 22 Wash two times with  30 μL Wash Buffer (30% MeOH, 0.1% FA), spin at 1000 x g for  00:01:00  00:01:00
- 23 Elute two times with  20 μL Elution Buffer (30% MeOH, 1.65% NH_4OH), spin at 1000 x g for  00:01:00  00:01:00
- 24 Vacuum dry for  00:30:00 at  42 °C and store at  -20 °C until MS-analysis



25 For MS analysis, dissolve samples in Solvent A and inject 1.75 ug