

Jan 30, 2025

Tissue Mito-IP: Immunopurification of Mitochondria from MitoTag Mice

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvo98rdv4o/v1>

Christos Themistokleous¹, Miratul M K Muqit¹

¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Str, Dundee DD1 5EH, UK



Christos Themistokleous

University of Dundee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvo98rdv4o/v1>

Protocol Citation: Christos Themistokleous, Miratul M K Muqit 2025. Tissue Mito-IP: Immunopurification of Mitochondria from MitoTag Mice. [protocols.io](https://dx.doi.org/10.17504/protocols.io.5qpvo98rdv4o/v1) <https://dx.doi.org/10.17504/protocols.io.5qpvo98rdv4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 06, 2025

Last Modified: July 01, 2026

Protocol Integer ID: 118216

Keywords: Tissue Mito-IP, Immunopurification of Mitochondria, MitoTag Mice, immunopurification of mitochondria, mitochondrial isolation, purifying mitochondria, mitochondria from tissue, mitochondrial proteome, changes in the mitochondrial proteome, intact mitochondria, mitochondria, outer mitochondrial membrane protein omp25, immunopurification of other organelle, mitochondrial viability, mitotag mice cell organelle, isolation of other organelle, compromising organelle purity, components of conventional organellar isolation buffer, mitotag mice, organelle purity, tissue mito, conventional organellar isolation buffer, cell profiling, immunopurification, mass spectrometry, other organelle, rapid immunopurification, intact organelle, release of intact organelle, protocol details the immunoprecipitation, organelle, chimeric protein, success of the mito, including lysosome

Abstract

Cell organelles represent a minor fraction of the total cellular content, making whole-cell profiling inadequate for monitoring changes in the mitochondrial proteome, metabolome, and lipidome. Traditional techniques for purifying mitochondria have inherent limitations, often compromising organelle purity, isolation time, or viability. Additionally, the components of conventional organellar isolation buffers, such as sucrose, can interfere with mass spectrometry (MS) profiling. To overcome these challenges, a novel method called 'Mito-IP,' was developed which facilitates the rapid immunopurification of pure and intact mitochondria. This method enables mitochondrial isolation within 10 minutes and supports various downstream applications, including immunoblotting, proteomic, metabolomic, and other -omic analyses. The method employs a chimeric protein comprising of three HA epitope tags fused to the outer mitochondrial membrane protein OMP25 (rat, aa109-145). An LC/MS-compatible buffer (termed 'KPBS') was also developed containing only KCl and KH_2PO_4 , significantly improving performance and mitochondrial viability. Following the success of the Mito-IP method, the same epitope-tagged concept is being extended to the isolation of other organelles, including lysosomes (Lyso-IP), Golgi (Golgi-IP), and peroxisomes (Peroxo-IP). The following optimised protocol details the immunoprecipitation of mitochondria from tissues of MitoTag mice. The same steps apply to the immunopurification of other organelles when the HA-epitope tag is present on the organelle of interest.




Fig. 1. Schematic of the Mito-IP method.

Tissues are collected and homogenized using a dounce homogenizer in KPBS buffer. After the release of intact organelles from the cells, the homogenate is incubated with anti-HA magnetic beads for 5 minutes, allowing the beads to bind to the mitochondria. Using a magnetic separator, the mitochondria-bound beads are drawn to the side of the tube. Following three washes, an appropriate lysis buffer is added for downstream analysis.

Materials

-  DPBS no calcium no magnesium **Gibco - Thermo Fisher Scientific Catalog #14190169**

- KPBS Buffer:

	A	B
	KCl	136mM
	KH ₂ PO ₄	10mM

Adjust to pH 7.25 with KOH and filter (0.22µm conring). On the day of use, add

-  cComplete™ EDTA-free Protease Inhibitor Cocktail **Roche Catalog #11873580001** and Roche PhosSTOP tablet (#04906837001).

- Lysis Buffer
- MS grade water (Fisher, 11947199)
- Pierce™ BCA Protein Assay Kit (#23227, lot# VA294738).

-  Pierce BCA protein assay **Thermo Scientific Catalog #23227**

- NuPAGE 4xLDS sample buffer (Invitrogen, #1941674).
-  1% β-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**

Equipment:

Equipment

VWR GLASS TISSUE GRINDER 2ML

NAME

Homogenizer

TYPE

neobits

BRAND

89026-386-EACH

SKU

https://www.neobits.com/vwr_89026_386_each_vwr_glass_tissue_grinder_2ml_p3461292.html?srltid=AfmBOopMUFg2KjilCr83nFVhOC03xxabSjciG3yYHqPT5OE6wEAVZpBQ

LINK

Equipment

VWR PL TISSUE GRINDER 2ML (Each)

NAME

plain plunger

TYPE

neobits

BRAND

89026-398-EACH

SKU

https://www.neobits.com/vwr_89026_398_each_vwr_pl_tissue_grinder_2ml_each__p3461298.html^{LINK}

•

Equipment

The Belly Dancer Shaker (Orbiter)

NAME

The Belly Dancer®

BRAND

BDRAA115S

SKU

<https://www.ibisci.com/products/belly-dancer-shaker>^{LINK}

- DynaMag™ Magnet. (Thermofisher scientific, #12320D or #12321D).

Equipment

DynaMag™- Spin Magnet

NAME

Invitrogen™

BRAND

12320D

SKU

<https://www.thermofisher.com/order/catalog/product/12320D>

LINK

Equipment

DynaMag-2

NAME

Magnet

TYPE

Invitrogen

BRAND

12321D

SKU

<https://www.thermofisher.com/order/catalog/product/12321D#/12321D>

LINK

•

Equipment

Microcentrifuges, Micro Star 17R (VWR #521-1647)	NAME
Microcentrifuges	TYPE
Micro Star 17R	BRAND
521-1647	SKU
https://in.vwr.com/store/product/8306728/microcentrifuges-ventilated-refrigerated-micro-star-17-17r ^{LINK}	

- Bioruptor NGS Sonication System (UCD-600).
- Dissecting tools.

Consumables:

-  Pierce & Warriner; Anti-HA Magnetic Beads **Thermo Fisher Catalog #88837**
-  SafeSeal reaction tube 1.5 ml PP PCR Performance Tested Low protein-binding **Sarstedt Catalog #72.706.600**
-  15 ml centrifuge tubes **greiner bio-one Catalog #188271**
-  50 ml centrifuge tubes **greiner bio-one Catalog #227261**
- TipOne bevelled 1000µl, 200µl and 20µl pipette tips (Starlab).



(A) Pre-clearing of anti-HA magnetic beads

5m

1 Resuspend beads by shaking the bottle until there is a homogeneous suspension.

2 Pipette  200 μL of anti-HA bead slurry into a 1.5ml tube (100 μL /sample).



Note

The volume of beads used per sample can be adjusted.

3 Place the tube onto a magnet for  00:00:30 and remove the overlying solution.

30s

4 Remove the tube from the magnet, add  1 mL of cold KPBS and gently resuspend 3 times to disperse any clumps.



5 Repeat the wash step (3+4) 2 more times.



6 After the last wash, resuspend the beads in  200 μL KPBS, split them into two tubes of 100 μL (one for the Mito-IP and one for the Control-IP) and keep them  On ice .



(B) Tissue collection and homogenisation

10m

7 Euthanize the mouse by cervical dislocation.

8 Dissect and collect the tissue of interest in a tube containing cold PBS.

Note

Make sure to remove any fat tissue connected to it if applicable.

9



**Note****Continue in the cold room until the detergent lysis step**

Place the tissue in the prechilled homogeniser and add  1 mL of KPBS.

- 10 Homogenise the tissue with 25 strokes without rotation. 

Note

Wash the douncer 10 times with milliQ water before using it again and always start with the Control-IP to avoid any contamination.

- 11 Move  25 μ L into a new tube (**whole-cell input**) and keep it  On ice . 

- 12 Spin down at  2000 x g, 4°C, 00:02:00 . This step is to remove any non-lysed cells, nuclei and debris such as cell membranes and leave the intact organelles in the supernatant.  2m

Note

Critical step: Optimise the homogenisation process by following the steps of section D [**Optimization of tissue homogenisation**].

- 13 Move  25 μ L to a new tube (**IP input**) and keep it  On ice . 

(C) Immunoprecipitation

- 14 Move the homogenate into the tube containing the HA-beads. Ensure that the bead clumps are dispersed by gently pipetting up and down 2-3 times. 

**Note****Critical step:**

Be very careful not to take any of the pellet as it will affect the purity of the IP. For some tissue homogenates, the separation between the pellet and the supernatant is not very clear. In this case, is preferred to spin down the homogenate at 2000xg instead and take less supernatant. See section D [**Optimization of tissue homogenisation**].

- 15 Incubate the homogenate with the beads for  00:05:00 on an orbiter.

5m

**Note**

Incubation time can be adjusted.

- 16 Separate the beads by putting the tube on the magnet for  00:00:30 .

Note

Residual material left on the cap can be brought down with a pulse spin ( 00:00:01) before placing the tubes on the magnet.

- 17 Collect  25 μL of flow-through into a new tube (**FT input**) and keep it  On ice .



- 18 Wash beads with  1000 μL of cold KPBS 3 times.

**Note**

Move the beads to a new tube after every wash to minimise contamination from cell extracts sticking to the tube.

- 19 The isolated mitochondria on the beads can be either stored at  -80 °C for later use or eluted off the beads using an appropriate lysis buffer for:





- Immunoblot analysis.
- Proteomic analysis.
- Lipidomic analysis.
- Metabolomic analysis.
- Other.

(D) Elution and lysis

40m 30s

- 20 Incubate the beads with  50 μ L lysis buffer for  00:10:00 . 10m 
- 21 Place the tube on the magnet for  00:00:30 and collect the supernatant into a new tube. Repeat this step one more time to ensure total removal of the beads. 30s
- 22 For the input samples, add  100 μ L of lysis buffer, resuspend and incubate for  00:20:00 . 20m  
- 23 Leave samples  On ice and continue with the next round of IP. 
- 24 Spin down input samples at  13000 x g, 4°C,  00:10:00 and move supernatant into a new tube. 10m  
- 25 For samples intended for immunoblot or proteomic analysis: Sonicate the IP and input samples for maximal protein extraction (15 cycles – 30sec on, 30sec off). Determine the protein concentration of the IP and inputs using a micro BCA Protein Assay Kit. Use 2 μ l of undiluted IP samples and 2 μ l of diluted (1:10) input samples. Quantification should be done in triplicate or duplicate.
- 25.1 Samples can be stored at  -80 °C for future application. 
- 25.2 For immunoblot analysis: samples are diluted into 4xLDS loading buffer supplemented with fresh 5% (by vol) 2-mercaptoethanol prior to analysis on SDS-polyacrylamide gel electrophoresis and immunoblot analysis. Loading  1 μ g -  2 μ g of sample is enough to detect the mitochondrial markers (such as HSP60, CS, VDAC, CISD1, OPA1, MFN1/2).
- 26 **For samples intended for proteomic analysis, use the protocol**



'Sample preparation for proteomic analysis of isolated mitochondria and whole-cell extracts'

[dx.doi.org/10.17504/protocols.io.3byl4wrj8vo5/v1](https://doi.org/10.17504/protocols.io.3byl4wrj8vo5/v1)

(E) Optimization of tissue homogenisation

8m

27 Determining the minimum amount of homogenization necessary is critical. Poor homogenization leads to inadequate amounts of mitochondria being available for immunocapture. However, excessive homogenization can damage the organelles released from cells. This section provides troubleshooting guidance for various mouse tissues to ensure the successful isolation of pure mitochondria.

28 Brain

28.1 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS. 

28.2 Homogenise the tissue using 40 strokes while avoiding introducing any bubbles.

28.3 Pour the homogenate into a 1.5ml tube slowly, as it will be viscous.

28.4 Spin down at  2000 x g, 4°C, 00:02:00 .

2m



28.5 Move  750 μ L of supernatant into the tube with the beads. There will be a big pellet. 

28.6 Continue with the steps of section C [Immunoprecipitation].

Note

Homogenising just one hemisphere or both separately, improves the homogenisation efficiency and Mito-IP purity.

29 Lung



29.1 Place the tissue at the bottom of the homogeniser and add  200 μL KPBS. 

29.2 Homogenise the tissue with an initial 5-7 strokes until there are no big tissue pieces.

29.3 Add the remaining  800 μL KPBS and perform another 25 strokes. With each stroke, move the probe all the way out of the buffer and then all the way back to the bottom of the homogeniser. 

29.4 Pour the homogenate into a 1.5ml tube.

29.5 Spin down at  2000 x g, 4°C, 00:02:00 .

2m



29.6 Move  900 μL of supernatant into the tube with the beads. 

29.7 Continue with the steps of section C **[Immunoprecipitation]**.

30 Heart

30.1 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS. 

30.2 Homogenise the tissue using 25 strokes. Make sure you push the probe all the way to the bottom of the homogeniser with the first stroke forcing the tissue to the sides of the probe.

30.3 Pour the homogenate into a 1.5ml tube slowly.

30.4 Spin down at  2000 x g, 4°C, 00:02:00 .



30.5 Move  900 μL of supernatant into the tube with the beads. 



30.6 Continue with the steps of section C [**Immunoprecipitation**].

31 Kidney

31.1 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS. 

31.2 Homogenise the tissue using 40 strokes while avoiding introducing any bubbles.

31.3 Pour the homogenate into a 1.5ml tube slowly.

31.4 Spin down at  2000 x g, 4°C, 00:02:00 .  

31.5 Move  900 µL of supernatant into the tube with the beads. 

31.6 Continue with the steps of section C [**Immunoprecipitation**].

Note

Homogenising only one kidney or both separately, improves the homogenisation efficiency and Mito-IP purity.

32 Spleen

32.1 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS. 

32.2 Homogenise the tissue using 25 strokes. Make sure you push the probe all the way to the bottom of the homogeniser with the first stroke forcing the tissue to the sides of the probe.

32.3 Pour the homogenate into a 1.5ml tube slowly.



32.4 Spin down at  2000 x g, 4°C, 00:02:00 .



32.5 Move  900 μ L of supernatant into the tube with the beads.



32.6 Continue with the steps of section C [**Immunoprecipitation**].

33 Liver

33.1 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS.



33.2 Homogenise the tissue using 25 strokes. Make sure you push the probe all the way to the bottom of the homogeniser with the first stroke forcing the tissue to the sides of the probe.

33.3 Pour the homogenate into a 1.5ml tube slowly.

33.4 Spin down at  2000 x g, 4°C, 00:02:00 .

2m



33.5 Move  900 μ L of supernatant into the tube with the beads.



33.6 Continue with the steps of section C [**Immunoprecipitation**].

Note

Homogenising only one liver lobe or more separately, improves the Mito-IP purity.

34 Skeletal muscle (gastrocnemius and soleus)



34.1 Place the tissue at the bottom of the homogeniser and add  200 μL KPBS. 

34.2 Homogenise the tissue with an initial 5-7 strokes until there are no big tissue pieces.

34.3 Add the remaining  800 μL KPBS and perform another 40 strokes. With each stroke, move the probe all the way out of the buffer and then all the way back to the bottom of the homogeniser. 

34.4 Pour the homogenate into a 1.5ml tube.

34.5 Spin down at  2000 x g, 4°C, 00:02:00 .

2m



34.6 Move  900 μL of supernatant into the tube with the beads. 

34.7 Continue with the steps of section C [**Immunoprecipitation**].

Note

Homogenising the skeletal muscles of both the hindlimbs together can improve the Mito-IP efficiency. Skeletal muscles are more difficult to homogenise and will require to apply more force, especially during the initial strokes.

35 Intestine (small or large)

35.1 Dissect a 5-6 cm segment and remove its contents using an angled probe. Intestines are surrounded by membranes. If not removed, they will not be able to be pelleted down and will affect the Mito-IP purity.

35.2 Place the tissue at the bottom of the homogeniser and add  1 mL KPBS. 

35.3 Homogenise the tissue using 40 strokes. Make sure you push the probe all the way to the bottom of the homogeniser with the first stroke forcing the tissue to the sides of the probe.



35.4 Pour the homogenate into a 1.5ml tube slowly.

35.5 Spin down at  5000 x g, 4°C, 00:02:00 .



35.6 Move  800 μ L of supernatant into the tube with the beads.



35.7 Continue with the steps of section C [**Immunoprecipitation**].