

Jun 30, 2025

Tissue processing for immunofluorescence

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl41m42lo5/v1>

Benjamin David Webster Belfort¹, Anthony Insalaco¹

¹Baylor College of Medicine



Benjamin David Webster Belfort

Baylor College of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl41m42lo5/v1>

Protocol Citation: Benjamin David Webster Belfort, Anthony Insalaco 2025. Tissue processing for immunofluorescence. protocols.io <https://dx.doi.org/10.17504/protocols.io.3byl41m42lo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We have published data resulting from this protocol.

Created: June 26, 2025



Last Modified: June 30, 2025

Protocol Integer ID: 221065

Keywords: ASAPCRN, murine olfactory bulb tissue processing, tissue processing for immunofluorescence protocol, immunofluorescence protocol, tissue processing, confocal microscopy application, cryopreservation, tissue fixation, omp labeling, including tissue fixation, tissue

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 10000400

Abstract

Protocol for murine olfactory bulb tissue processing, including tissue fixation, cryopreservation, and OMP labeling. This protocol is intended for confocal microscopy applications.

Olfactory bulb harvesting, fixation, and cryopreservation

- 1 Begin by deeply anesthetizing the mouse using any preferred method (following institutionally approved protocols). Confirm surgical plane of anesthesia by firm toe pinch and proceed with transcardiac perfusion of at least 10 mL of ice-cold 1x PBS followed by 10 mL of 4% paraformaldehyde (16% PFA diluted in PBS; Electron Microscopy Sciences Cat# 15710).
- 2 Harvest olfactory bulbs (OBs) and fix in 4% PFA at 4°C overnight.
- 3 Move OBs to 30% sucrose at 4°C for 36h for cryoprotection.
- 4 Freeze OBs in OCT (Fisher HealthCare Cat# 23-730-571).

Tissue sectioning and OMP labeling

- 5 Section OBs at 40 μ m on a cryostat (Leica Cat# CM1860), collecting every third section for analysis.
- 6 Place tissue sections in blocking buffer (5% Donkey Serum, 1X PBS, 0.1% Triton) for 1 h at room temperature.
- 7 Wash sections once in PBST (0.1% Triton).
- 8 Incubate sections with anti-OMP antibody (1:20,000; FUJIFILM Wako Chemicals U.S.A. Corporation Cat# 544-10001-WAKO) in blocking buffer at 4°C overnight.
- 9 Wash tissue sections three times with PBST (0.1% Triton).
- 10 Incubate sections with secondary antibody (1:1,000 anti-goat, Alexa Fluor 488 Invitrogen Cat# A-11055) in blocking buffer for 1 h at room temperature.
- 11 Wash sections three times with 0.1% PBST.



- 12 Incubate sections with Hoechst (1:1,000, Thermo Fisher Cat# 62249) in 0.1% PBST for 15 min at room temperature.
- 13 Wash sections once more with 0.1% PBST.
- 14 Mount sections on slide with Fluoromount-G (Southern Biotech Cat# 0100-01).

Guidelines and Warnings

- 15 This protocol requires prior approval by the users' Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee.