

Aug 09, 2023

TMT labelling

DOI

dx.doi.org/10.17504/protocols.io.n92ldmp69l5b/v1

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis

WEHI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmp69l5b/v1>

Protocol Citation: Louise Uoselis 2023. TMT labelling. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldmp69l5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 86211

Keywords: ASAPCRN, labelling tmt, tmt

Abstract

TMT labelling

Troubleshooting



- 1 Reconstitute lyophilised peptides in 25 μL of 100 mM TEAB by vortexing each sample for ~ 00:00:10 , and then leaving to sit at Room temperature for 00:15:00 . Each sample was then sonicated for 00:05:00 in a waterbath sonicator in an ice slurry. 20m 10s
- 2 Determine the peptide concentration of each sample spectroscopically
- 3 Aliquot 10 μg of peptides for each sample into a LoBind microfuge tube and add 100 millimolar (mM) triethylammonium bicarbonate to a final volume of 100 μL for each sample.
- 4 For the pooled batch control, divide 100 μg by the total number of samples, and aliquot that amount of peptides from each sample into one tube. Make up the pooled control volume to 100 μL with 100 millimolar (mM) triethylammonium bicarbonate (can lyophilise the samples to reduce the volume if the total pooled volume exceeds 100 μL).
- 5 Bring the TMT labels to Room temperature .
- 6 Reconstitute each label in acetonitrile as per the manufacturer's instructions.
- 7 Add the required volume of the designated label to each sample, and vortex each sample for ~ 00:00:05 to mix. 5s
- 8 Leave samples to incubate at Room temperature for 01:00:00 (static). 1h
- 9 Add hydroxylamine to a final concentration of 0.26%/sample, vortex each sample for ~ 00:00:05 , and incubate at Room temperature for 00:15:00 to quench the TMT labelling reaction. 15m 5s



- 10 Combine the samples into their designated batches, and lyophilise the pooled samples. Seal each tube with parafilm and store at 🌡️ -80 °C for downstream processing.