

Jan 12, 2026

Tn5 loading scalebarcodes

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk1w1wg5r/v1>

Koen Theunis^{1,2,3,4}

¹Laboratory of Computational Biology, VIB Center for AI & Computational Biology, Leuven, Belgium;

²VIB-KU Leuven Center for Brain & Disease Research, Leuven, Belgium;

³Department of Human Genetics, KU Leuven, Leuven, Belgium;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, United States

Aertslab



Koen Theunis

Laboratory of Computational Biology, VIB Center for AI & Com...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk1w1wg5r/v1>

Protocol Citation: Koen Theunis 2026. Tn5 loading scalebarcodes. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlk1w1wg5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 09, 2026

Last Modified: January 12, 2026

Protocol Integer ID: 238301

Keywords: ASAPCRN, Tn5, scaleATAC, scATAC-seq, tn5 transposase, scalebarcodes this protocol, barcoded adapter complex, loading scalebarcode, tn5, adapter complex

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000430

Aligning Science Across Parkinson's

Grant ID: ASAP-025179

Abstract

This protocol describes how to prepare barcoded adapter complexes that can be loaded onto Tn5 transposase.

Attachments



[Tn5_loading_scalebar...](#)

13KB

Guidelines

All mixes contain glycerol so make sure to pipet slowly and sufficient for proper mixing.

Materials

	A	B	C	D	E	F
	Oligo anneal buffer B	Stock conc		Final conc		V (μl)
	Tris-HCl (pH 8)	1000	mM	40	mM	40,0
	NaCl	5000	mM	50	mM	10,0
	UltraPure Water					950,0
	Total Volume					1000,0

	A	B	C	D
	Name	Truseq_R2_SBS12_partial	Tn5_96plx_idx	ME
	R2_18_ME_scl_4	CGTGTGCTCTTCCGATCT	CTGCTTCC	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_5	CGTGTGCTCTTCCGATCT	GGTCACGA	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_6	CGTGTGCTCTTCCGATCT	AACTGTAG	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_16	CGTGTGCTCTTCCGATCT	GCGGCACA	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_17	CGTGTGCTCTTCCGATCT	CTAGCGCT	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_18	CGTGTGCTCTTCCGATCT	TCGATATC	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_28	CGTGTGCTCTTCCGATCT	CGCCGATC	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_29	CGTGTGCTCTTCCGATCT	TCTCTACT	AGATGTGTATAAGAGACAG
	R2_18_ME_scl_30	CGTGTGCTCTTCCGATCT	CTCTCGTC	AGATGTGTATAAGAGACAG



	A	B	C	D
	R2_18_ME_scl_40	CGTGTGCTCTTCCGATCT	AACGTTCC	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_41	CGTGTGCTCTTCCGATCT	GCAGAATT	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_42	CGTGTGCTCTTCCGATCT	ATGAGGCC	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_52	CGTGTGCTCTTCCGATCT	TGGATCGA	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_53	CGTGTGCTCTTCCGATCT	AGTTCAGG	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_54	CGTGTGCTCTTCCGATCT	GACCTGAA	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_64	CGTGTGCTCTTCCGATCT	GCGCAAGC	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_65	CGTGTGCTCTTCCGATCT	AAGATACT	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_66	CGTGTGCTCTTCCGATCT	GGAGCGTC	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_76	CGTGTGCTCTTCCGATCT	TCCAACGC	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_77	CGTGTGCTCTTCCGATCT	CCGTGAAG	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_78	CGTGTGCTCTTCCGATCT	TTACAGGA	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_88	CGTGTGCTCTTCCGATCT	TAAGGTCA	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_89	CGTGTGCTCTTCCGATCT	TTGCCTAG	AGATGTGTATAAGAGACA G
	R2_18_ME_scl_90	CGTGTGCTCTTCCGATCT	CCATTCGA	AGATGTGTATAAGAGACA G
	Name	Nextera_R1_A14		ME
	R1-A-oligo	TCGTCGGCAGCGTC		AGATGTGTATAAGAGACA G
	Name			ME'

	A	B	C	D
	Nextera_Mosaic_End_REV COMP			/5Phos/CTGTCTCTTATAC ACATCT

Order HPLC purified oligo's

	A	B	C	D	E	F
	Tn5 dilution buffer	Stock conc		Final conc		V (μl)
	Tris-HCl (pH 7.4)	1000	mM	50	mM	25 0,0
	NaCl	5000	mM	100	mM	10 0,0
	EDTA	50	mM	0,1	mM	10, 0
	DTT	1000	mM	1	mM	5,0
	Glycerol	100	%	50	%	25 00, 0
	Water					213 5,0
	Total Volume					50 00, 0

Add DTT fresh



Prepare oligo's

- 1 Dissolve oligo's in Oligo annealing buffer (40mM Tris-HCl (pH 8), 50mM NaCl) to a stock of [M] 100 micromolar (μM) .
Work in 96-well plates for easy pipetting.

Anneal oligo's

- 2 Mix oligo's equimolar for each barcode:

	A	B	C	D
		μL	Stock	
	A=R1	3	100	μM
	B=R2-barcodedoligo1-24	3	100	μM
	ME'	6	100	μM

- 3 5min 95°C, ramp (-1°C/s) to 65°C, 5min 65°C, ramp (-1°C/s) to 4°C.

A	B
Temperature (°C)	Time
95 °C	5min
ramp (-1°C/s)	
65 °C	5min
ramp (-1°C/s)	
4 °C	Hold

Adapters are now [M] 50 micromolar (μM) .

Dilute oligo's and Tn5

- 4



A	B	C	D	E	F
	Stock conc		Final conc		V (μl)
oligo mix	50	μM	24	mM	12,0
Tn5 dil buffer					13,0
Total Volume					25,0

5

A	B	C	D	E
	Storage	Vol (μL)	Conc (ng/μL)	Molarity (μM)
EMBL Tn5	-80	20	1000	18

A	B	C	D	E	F
Tn5 (diluted in fresh dilution buffer)	Stock conc		Final conc		V (μl)
Tn5	18	μM	6	μM	20,0
Tn5 dil buffer					40,0
Total Volume					60,0

Tn5 is now 333,33 ng/μL; 6 μM

Load Tn5

6



A	B	C	D	E	F
Tn5 loading on ice	Stock conc		Final conc		V (μl)
Tn5 (diluted in fresh dilution buffer)	6	μM	5	μM	5,0
adapter1-24	24	μM	4	μM	1,0
Total Volume					6,0

Repeat this for every barcoded adapter mix.
Adapter loaded Tn5 is now max at 4 μM = 222,22ng/μL functional Tn5.

7 Incubate for  00:45:00 @  23 °C on a heater/shaker (300 RPM).

Dilute Tn5 to stock concentration and store

- 8 Dilute to 62,5 ng/μL for 50 μL reactions
Dilute to 25 ng/μL for 15 μL tagmentation reactions (like 10x Genomics ATAC)
- 9 Store at -20 °C.

Protocol references

Our barcoded oligo's are modified versions of the oligonucleotide list that was published in: Mulqueen, R.M., Pokholok, D., O'Connell, B.L. *et al.* High-content single-cell combinatorial indexing. *Nat Biotechnol* **39**, 1574–1580 (2021). <https://doi.org/10.1038/s41587-021-00962-z>

DIY Tn5 protocols would not have been possible without the work of: Simone Picelli, Åsa K. Björklund, Björn Reinius, Sven Sagasser, Gösta Winberg and Rickard Sandberg. Tn5 transposase and tagmentation procedures for massively scaled sequencing projects. <http://www.genome.org/cgi/doi/10.1101/gr.177881.114>.

The source of our Tn5 is EMBL and they describe it in the following paper:

Hennig BP, Velten L, Racke I, Tu CS, Thoms M, Rybin V, Besir H, Remans K, Steinmetz LM. Large-Scale Low-Cost NGS Library Preparation Using a Robust Tn5 Purification and Tagmentation Protocol. *G3 (Bethesda)*. 2018 Jan 4;8(1):79-89. doi: 10.1534/g3.117.300257. PMID: 29118030; PMCID: PMC5765368.