


Jan 05, 2026

Tomogram reconstruction and pre-processing using standardized workflow in Scipion-EM

DOI

<https://dx.doi.org/10.17504/protocols.io.36wggpo45vk5/v1>



Benoit Gallet¹, Guy Schoehn¹

¹Univ. Grenoble Alpes, CEA, CNRS, IBS, F-38000 Grenoble



Benoit Gallet

IBS Grenoble

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wggpo45vk5/v1>

Protocol Citation: Benoit Gallet, Guy Schoehn 2026. Tomogram reconstruction and pre-processing using standardized workflow in Scipion-EM. protocols.io <https://dx.doi.org/10.17504/protocols.io.36wggpo45vk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 11, 2025

Last Modified: January 05, 2026

Protocol Integer ID: 227018

Keywords: Cryo-electron tomography, scipion-em, image processing, segmentation, Cryo-tomogram, GPU processing, traceable workflows for electron microscopy, modular workflow for tomogram reconstruction, electron microscopy, em electron tomography, tomogram reconstruction, microscopy, architecture of biological specimen, standardized workflow in scipion, nanometer resolution, source image processing platform, biological specimen, cell biology, subcellular structure, traceable workflow, specimen, implementing standardized pipeline

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Electron tomography (ET) is a powerful technique for determining the three-dimensional (3D) architecture of biological specimens at nanometer resolution. It is widely used in structural biology and cell biology to analyze subcellular structures in near-native conditions. However, tomogram reconstruction is a complex, multi-step process involving data pre-processing, alignment, reconstruction, and post-processing. The diversity of tools and the flexibility of workflows can lead to inconsistent results, making reproducibility and traceability a significant challenge. To address this, we present a standardized and modular workflow for tomogram reconstruction using the Scipion-EM framework. Scipion-EM is an open-source image processing platform that integrates multiple software packages within a unified environment, enabling reproducible and traceable workflows for electron microscopy (EM) data processing. The platform supports automation, provenance tracking, and user-friendly graphical interfaces, making it particularly suitable for implementing standardized pipelines.

Guidelines

This protocol outlines a standardized tomogram reconstruction workflow within Scipion-EM, with an emphasis on:

- Reproducibility and traceability of results;
- Interoperability between software packages (e.g., IMOD, Dynamo, RELION, Warp/M, etc.);
- Flexibility for customization based on project-specific requirements;
- Best practices for data management and metadata handling.

The goal is to enable both novice and experienced users to execute robust tomographic reconstructions, reduce methodological variability, and facilitate the sharing and publication of processing pipelines.

This protocol is intended for:

- Researchers and technicians engaged in electron tomography data processing;
- Cryo-EM core facility personnel developing training resources;
- Developers interested in expanding or benchmarking workflows within Scipion-EM;
- Collaborators seeking to replicate or audit published results using standardized methods.

While some knowledge with tomographic concepts and EM software is prerequisites, this guide is designed to be accessible to users new to Scipion-EM.

Materials

Software/tools referenced on these pages: Scipion-EM GUI, AreTomo (ProtAreTomoAlignRecon), IMOD/3dmod (tomogram viewing, dose filtering and apply transformation), MotionCorr/ProtTsMotionCorr, CistemProtTsCtffind (CTF estimation), CTFEstimationTomoViewer, ProtImodDoseFilter (IMOD - Dose filter), ProtImodApplyTransformationMatrix (IMOD - Apply transformation), ProtImodTomoReconstruction (IMOD - Tomo reconstruction), CryoCARE (ProtCryoCARETraining and ProtCryoCAREPrediction/CryoCARE Prediction), Warp (ProtWarpDeconvTomo - warp deconvolve tomograms), tomo3D (ProtTomo3dProtDenoiseTomogram - denoise tomogram, Edge Enhancing Diffusion (EED) option), MemBrain (ProtMemBrainSeg - membrain - tomogram membrane segmentation), ChimeraX (3D rendering).

Hardware Requirements

While Scipion-EM is platform-independent, the following hardware is recommended for efficient tomographic reconstruction:

- CPU: Multi-core processor (8 cores or more recommended)
- GPU: CUDA-compatible GPU (for GPU-accelerated tools like MotionCorr2, AreTomo, Membrain-seg, CryoCARE or RELION)
- RAM: Minimum 64 GB or more preferred for large datasets
- Storage: SSD or high-speed disk with sufficient capacity for intermediate files and reconstructed tomograms (8Tb minimum) and SATA discs for storage (100Tb is a good start)

Safety warnings

- ! Ensure the following are installed and properly configured:
- Scipion-EM (v3.0 or higher) with tomography-related plugins:
- scipion-em-xmipp
 - scipion-em-xmipptomo
 - scipion-em-imod
 - scipion-em-motioncor
 - scipion-em-aretomo
 - scipion-em-relion (for CTF estimation or subtomogram averaging)
 - scipion-em-dynamo (for subtomogram averaging workflows)
 - scipion-em-tomo
 - scipion-em-warp
 - scipion-em-tomo3d
 - scipion-em-membrain
 - scipion-em-cryocare
- Python: Included with Scipion; no separate installation needed
- CUDA Toolkit (if using GPU-accelerated tools)
- Environment variables set appropriately (e.g., for IMOD's IMOD_DIR, CUDA paths)

Before start

To ensure consistency and reproducibility:

- Organize your data: Maintain a consistent folder structure (e.g., one folder per tilt-series with raw data and metadata)
- Backup raw data: Before starting any processing, create a read-only copy of your original files
- Version control: Record software versions and parameters, especially if you are working in a collaborative environment
- Check plugin updates: Scipion plugins evolve rapidly; always check for the latest stable versions with bug fixes or new features
- GPU vs CPU paths: Be aware of which steps benefit from GPU acceleration; plan your workflow accordingly if using a shared cluster

Organize your data in the following:

i. Tilt-Series Data

- Format: MRC, TIFF, or other supported image stacks
- Acquisition Metadata: Must include tilt angles, pixel size, defocus (if applicable), and microscope parameters
- Fiducials: Gold fiducials (10–20 nm) are recommended for robust alignment unless using fiducial-less methods

ii. Microscope and Acquisition Parameters (in .mdoc file)

- Voltage (e.g., 300 kV)
- Spherical aberration (Cs, 2.7 for Krios)
- Detector type (e.g., K3, Falcon 4)
- Exposure time and dose per tilt (Exposure dose)
- Acquisition scheme (dose-symmetric, bidirectional, etc.)
- Number of frames per movie
- Tilt angle (~85°)
- Pixel size (Å/pix2)

Accurate metadata is critical for CTF estimation, tilt alignment, and dose compensation steps, all these informations are available in the .mdoc file.


Start Scipion

- 1 Go to Document and create a new folder with the name of your experiment, copy your raw data inside (.mrc or .tiff/.eer AND .mdoc files). In the same folder, in which the movies are, you need the .gain file of the day of data collection.


Start scipion3 in a terminal.

Create new project

- 2 Create new project and give the same name that your data folder




On the left panel select the tomography tree for availability of adapted protocols



Workflow import and setup

- 3 Go to import workflow
Go to ScipionUserData/Workflows_template
Select the Workflow adapted to your project and data and import it



NB: you can always start a blank workflow and add each step as a new protocol from the list.



Import data

- 4 What to do before: open a .mdoc file and identify the mandatory information: kV, Cs, pixel size, movies frames, electron dose.

4.1 In the case of .mrc stack files

Protocol to use: tomo - import tilt-series

You have to specify:

- Folder with your data
- Pattern of the .mdoc file containing information of the data (*.mdoc or Position*.mdoc or TS*.doc)



Protocol Run: ProtImportTs

tomo - import tilt-series finished Cite Help

Run

Run name: **tomo - import tilt-series** edit Comment: edit

Run mode: ☒ Continue ☐ Restart ? Use a queue engine?: ☐ Yes ☒ No edit ?

Wait for: ?

Expert Level: ☐ Normal ☒ Advanced

Import

Files directory: ?

Pattern: ?

Tomo5 mdoc?: ☐ Yes ☒ No ?

Exclusion words: ?

Acquisition values provided below will override the corresponding mdoc values

Have images been CTF corrected?: ☐ Yes ☒ No ?


Have images been aligned?: ☐ Yes ☒ No ?

Import action on files: ☐ Copy files ☐ Absolute symlink ☒ Relative symlink ?

Acquisition info - override mdoc values if provided

Microscope voltage (kV)	300.0	<a>?
Spherical aberration (mm)	2.7	<a>?
Amplitude Contrast	0.1	<a>?
Magnification rate	33000	<a>?
Pixel size (sampling rate) $\text{\AA}/\text{px}$	3.769	<a>?
Tilt axis angle (deg.)	83.7	<a>?
Dose (electrons/sq, $\text{\AA}$)	0.0	<a>?
Initial dose	0.0	<a>?
Dose per tilt image	2.23311	<a>?

Close Save Execute



4.2 In the case of individual tilt movies

Protocol to use: tomo - import tilt-series movies

Protocol Run: ProtImportTsMovies

tomo - import tilt-series movies finished [Cite](#) [Help](#)

Run

Run name [Edit](#) Comment [Edit](#)

Run mode ☒ Continue ☐ Restart [?](#) Use a queue engine? ☐ Yes ☒ No [Edit](#) [?](#)

Wait for: [?](#)

Expert Level ☒ Normal ☐ Advanced

Import

Files directory [Browse](#) [?](#)

Pattern [?](#)

Tomo5 mdoc? ☐ Yes ☒ No [?](#)

Acquisition values provided below will override the corresponding mdoc values

Acquisition info - override mdoc values if provided

Microscope voltage (kV) [?](#)

Spherical aberration (nm) [?](#)

Amplitude Contrast [?](#)

Pixel size (sampling rate) Å/px [?](#)

Tilt axis angle (deg.) [?](#)

Dose (electrons/sq.Å) Initial dose Dose per tilt image [?](#)

Gain image [Browse](#) [?](#)

[Close](#) [Save](#) [Execute](#)

Tilt series viewer

[Save](#) [Help](#)

Filter:


Tilt series	Order	Tilt angle	Excluded	Dose	Path
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 315, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 315, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 315, 3,75 Å/px)
▶ lamella1					TiltSeriesM (51 items, 4096 x 4096 x 531, 3,75 Å/px)

[Close](#)

Correct motion

- 5 In the case of .eer or .tiff movies of individual tilt, motion needs to be corrected using motion correction.

Protocol to use: motioncorr - align tilt-series movies



NB: In the case of using CryoCARE afterwards, it is mandatory to split odd/even frames!

Protocol Run: ProtTsMotionCorr

motioncorr - align tilt-series movies finished [Cite](#) [Help](#)

Run

Run name: [Edit](#) Comment: [Edit](#)

Run mode: ☒ Continue ☐ Restart [?](#) Use a queue engine? ☐ Yes ☒ No [Edit](#) [?](#)

Parallel Scipion threads: [?](#) Wait for: [?](#)

GPU IDs: [?](#)

Expert Level: ☒ Normal ☐ Advanced

Input Motioncorr params Gain and defects EER Mag. correction

Motioncorr params

Number of patches X: Y: [?](#)

Patches overlap (%): [?](#)

Group N frames global align: local align: [?](#)

[Close](#) [Save](#) [Execute](#)

Protocol Run: ProtTsMotionCorr

motioncorr - align tilt-series movies finished [Cite](#) [Help](#)

Run

Run name: [Edit](#) Comment: [Edit](#)

Run mode: ☒ Continue ☐ Restart [?](#) Use a queue engine? ☐ Yes ☒ No [Edit](#) [?](#)

Parallel Scipion threads: [?](#) Wait for: [?](#)

GPU IDs: [?](#)

Expert Level: ☒ Normal ☐ Advanced

Input Motioncorr params Gain and defects EER Mag. correction

Gain and defects

Rotate gain reference: [?](#)

Flip gain reference: [?](#)

Camera defects file: [?](#)

Camera defects map: [?](#)

[Close](#) [Save](#) [Execute](#)

Protocol Run: ProtTsMotionCorr

motioncorr - align tilt-series movies

finished Cite Help

Run

Run name motioncorr - align tilt-s

Comment

Run mode ☒ Continue ☐ Restart

Use a queue engine? ☐ Yes ☒ No

Parallel Scipion threads 1

Wait for

GPU IDs 0

Expert Level ☒ Normal ☐ Advanced

Input

Motioncorr params

Gain and defects

EER

Mag. correction

EER

These options are ignored for non-EER movies.

EER fractionation 10

EER upsampling ☒ 1x ☐ 2x ☐ 4x

Close

Save

Execute

Protocol Run: ProtTsMotionCorr

motioncorr - align tilt-series movies

finished Cite Help

Run

Run name motioncorr - align tilt-s

Comment

Run mode ☒ Continue ☐ Restart

Use a queue engine? ☐ Yes ☒ No

Parallel Scipion threads 1

Wait for

GPU IDs 0

Expert Level ☒ Normal ☐ Advanced

Input

Motioncorr params

Gain and defects

EER

Mag. correction

Mag. correction


Correct anisotropic magnification? ☐ Yes ☒ No

Close

Save

Execute

Open the results:



Align tilt series and evaluate the quality of the data

6 AreTomo, IMOD patch, IMOD fiducials

Protocol Run: ProtAreTomoAlignRecon

aretomo - tilt-series align and reconstruct finished Cite Help

Run

Run name aretomo - tilt-series ali Comment

Run mode ☒ Continue ☐ Restart Use a queue engine? ☐ Yes ☒ No

Scipion threads 2 Wait for

GPU IDs 0

Expert Level ☐ Normal ☒ Advanced

Input Extra options

Input

Input set of Tilt-Series tomo - import tilt-series,outputTiltSeries

Skip alignment? ☐ Yes ☒ No

Reconstruct the tomograms? ☒ Yes ☐ No

Reconstruct the odd/even tomograms? ☐ Yes ☒ No

Binning 8

Volume height for alignment (voxels) 300

File with volume height for alignment per tilt-series

Tomogram thickness unbinned (voxels) 1200

Refine tilt angles? Measure only

Refine tilt axis angle? Refine and use the refined value for the

Generate extra IMOD output? No

CTF

Estimate the CTF? ☒ Yes ☐ No

Do phase shift estimation? ☐ Yes ☒ No

Close Save Execute

Select the imported data or the motion corrected data

- Skip alignment: No
- Reconstruct the tomogram: Yes
- Reconstruct Odd/Even: No

NB: In the case of the use of CryoCARE later on, you need to always activate the Odd/Even option, in all the steps.

- Binning: 8
- Volume height for alignment: try between 200 and 400
- Tomogram thickness unbinned : 1200
- Refine tilt angle: No
- Refine tilt axis angle: Refine and use the refine value for the entire "stack"
- Generate extra IMOD output: No
- Estimate the CTF: Yes
- Do phase shift estimation: No

Select "Extra Options" Tab:

Protocol Run: ProtAreTomoAlignRecon

aretomo - tilt-series align and reconstruct finished Cite Help

Run

Run name aretomo - tilt-series ali Comment

Run mode ☒ Continue ☐ Restart Use a queue engine? ☒ Yes ☐ No

Parallel Scipion threads 10 Wait for

GPU IDs 0

Expert Level ☒ Normal ☐ Advanced

Input Extra options

Extra options

Do dose-weighting? ☒ Yes ☐ No

Reconstruction method ☒ SART ☐ WBP

SART options iterations 10 projections per subset 5

Flip intensity? ☒ Yes ☐ No

Flip volume? ☒ Yes ☐ No

Local motion correction

Sample type Well distributed

Patches X 5 Y 5

Dark tolerance 0.7


Additional parameters TiltCor: 1

Close Save Execute


Select the SIRT method to have more contrast

- You can add the Patch tracking by selecting the option :
- Local motion correction - Well distributed 5 5
- Additional parameter: TiltCor 1

Check the result: (open output: Tomograms)




Open each, one by one, by double clicking on it: (you may need to transfer the files locally).




Check the tomogram and if the alignment is good, keep those parameters, otherwise restart this step trying to optimize it: increase/decrease the volume for alignment, use Patch (Local motion correction),...

Estimate CTF

7 Plug-in to use: cistem - tilt-series ctffind




Open result with CTFEstimationTomoViewer:



You can remove the bad tilt-series, or at least the bad images

Dose weighting

8 Protocol to use: **Imod - Dose filtering**



If the dose was correctly set in the import step you can select the Scipion import as Input dose source. Otherwise you can set the dose as fixed and enter the dose for each image (tilt) in electron / angstrom².


In the case of reconstructing the even and odd tomograms (for the use of CryoCARE denoising) you can apply this dose filtering to the halves, or not.

Tomogram reconstruction

9


9.1 Alignment

Binning 2, 4, 8 and purpose




The recommendations for the binning are the following:


- Start with Binning 8 because it's faster, give more contrast, is useful for quality check of the data and the alignments/reconstructions parameters
- Then do Binning 4 and Binning 2 (and Bin1 for highest resolution reconstruction)
- Bin 4 for denoising and segmentation of membranes, for morphological studies
- Bin 2 for high resolution interpretation of the data, particle picking and STA (Sub-Tomogram Averaging)




Bin8



Bin2



Bin4 (WBP)




Bin4 (SART)

Binning effects on reconstruction contrast and resolution

Binning effects on reconstruction contrast and resolution

9.2 Reconstruction




Plung-in: mod - Tomo reconstruction

Denoising

- 10 We propose three different ways for denoising: CryoCARE, Warp and tomo3D.
- 1) CryoCARE is a training-based approach, it is very long for the training but the result is of very high quality.
 - 2) Warp is a CTF based filtering approach that is very fast and gives very good results
 - 3) tomo3D is more a filtering method that can use either Anisotropic Nonlinear Diffusion (EED method) or TomoBFlow, a nonlinear filtering method, it is also a fast method for denoising.
- 10.1 To be able to use CryoCARE you need to have even and odd reconstructions of your tomograms. This is done at the "motion correction" step level and is also in option in all the subsequent steps by always selecting "Apply to odd/even".

1- Training:

Protocol to use: cryocare - CryoCARE training



Protocol Run: ProtCryoCARETraining

cryocare - CryoCARE Training finished [Cite](#) [Help](#)

Run

Run name: cryocare - CryoCARE Train [Edit](#) Comment: [Edit](#)

Run mode: ☒ Continue ☐ Restart [?](#) Use a queue engine? ☐ Yes ☒ No [Edit](#) [?](#)

EGPU IDs: 0,1 [?](#) Wait for: [?](#)

Expert Level: ☐ Normal ☒ Advanced

Input [Config Parameters](#) [Training Parameters](#) [U-Net Parameters](#)


Are odd-even associated to the Tomograms? ☒ Yes ☐ No

Tomograms: imod - Tomo reconstruction (copy 2).Tomograms [Search](#) [Delete](#) [Refresh](#)

[Close](#) [Save](#) [Execute](#)

2- Prediction:

In case you want to denoise a full set of tomograms from a single data collection session, you can train the CryoCARE NN with one representative tomogram and then use the trained NN to predict all the other tomograms.




3- Results:

Result:

Before

After denoising



CryoCARE denoising


10.2 Warp



Protocol to use : warp - deconvolve tomograms


Input: the reconstructed tomogram in Bin4 (or Bin8)

CTF tomo series: the CTF estimation from AreTomo or gctf



Plug-in: warp - deconvolve tomogram

Results:



Warp denoising

10.3 Protocol to use: tomo3d - denoise tomogram

Protocol Run: ProtTomo3dProtDenoiseTomogram

tomo3d - denoise tomogram finished [Cite](#) [Help](#)

Run

Run name: tomo3d - denoise tomogram [Edit](#) Comment: [Edit](#)

Run mode: ☒ Continue ☐ Restart [?](#) Use a queue engine? ☐ Yes ☒ No [Edit](#) [?](#)

Parallel Threads: 8 [?](#) Wait for: [?](#)

Expert Level: ☒ Normal ☐ Advanced

Input

Set Of Tomograms: imod - Tomo reconstruction.Tomograms [Search](#) [Trash](#) [Eye](#) [?](#)

Denoising method: Edge Enhancing Diffusion (EED) [?](#)


Sigma Gaussian Filter: 0.5 [?](#)

Number of Iterations: 10 [?](#)

[Close](#) [Save](#) [Execute](#)

If you want to explore the best settings you can play either with the Denoising method or with the parameters.

Results:



tomo3d denoising

Segmentation


11 Membrain-seg

Protocol to use : membrain - tomogram membrane segmentation

Membrain is an automated method to segment membranes in denoised cryo-tomograms. It uses a neural network that has been trained on *Chlamydomonas* data mainly. The result will be a .mrc file with the segmented membranes, it can be grey scaled according to connected components

"Input" tab:


Select the data: denoised tomograms (provide denoised tomograms as input in the Input tab of the ProtMemBrainSeg protocol).




The screenshot shows the 'Protocol Run: ProtMemBrainSeg' window. The title bar includes 'MemBrain' and 'membrain - tomogram membrane segmentation'. The status is 'finished' with 'Cite' and 'Help' links. The 'Run' section has fields for 'Run name' (membrain - tomogram membra), 'Comment', 'Run mode' (Continue, Restart), 'Use a queue engine?' (Yes, No), 'Threads' (1), 'Wait for' (679), and 'GPU IDs' (0). The 'Expert Level' is set to 'Advanced'. The 'Input' tab is selected, showing 'Input tomograms' (aretomo - tilt-series align and reconstruct (copy).Tomograms), 'Threshold for segmentation' (0.0), 'Sliding window size' (160), and 'Additional options'. At the bottom are 'Close', 'Save', and 'Execute' buttons.

“Connected components analysis” tab:

You can either choose to label connected component, or not, depending on the software you will use in a second step.




The screenshot shows the 'Protocol Run: ProtMemBrainSeg' window with the 'Connected components analysis' tab selected. The 'Run' section is the same as the previous screenshot. The 'Expert Level' is 'Advanced'. The 'Connected components analysis' section has 'Label connected components?' (Yes, No) and 'Threshold for connected components' (-1). At the bottom are 'Close', 'Save', and 'Execute' buttons.



"Test-time augmentation" tab: always "Yes"

Open the result to check, with 3dmod or with chimera (as volume)



Results of the segmentation

Downstream processes

- 12 Particle picking
- 13 Sub-Tomogram Averaging

Protocol references

<https://scipion-em.github.io/docs/release-3.0.0/docs/user/tutorials/tomo/tomography-intro.html#tomography-intro>

<https://www.sciencedirect.com/science/article/pii/S1047847722000429?via%3Dihub>

Acknowledgements

This work used the platforms of the Grenoble Instruct-ERIC center (ISBG ; UAR 3518 CNRS-CEA-UGA-EMBL) within the Grenoble Partnership for Structural Biology (PSB), supported by FRISBI (ANR-10-INBS-0005-02) and GRAL, financed within the University Grenoble Alpes graduate school (Ecoles Universitaires de Recherche) CBH-EUR-GS (ANR-17-EURE-0003).