

Sep 05, 2019

## Top2 Chromatin Accessibility by Etoposide Cross-linking

DOI

[dx.doi.org/10.17504/protocols.io.2rngd5e](https://dx.doi.org/10.17504/protocols.io.2rngd5e)

Jacob G Kirkland<sup>1</sup>

<sup>1</sup>Stanford University



Jacob G Kirkland

Oklahoma Medical Research Foundation, Oklahoma University He...

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.2rngd5e](https://dx.doi.org/10.17504/protocols.io.2rngd5e)

**Protocol Citation:** Jacob G Kirkland 2019. Top2 Chromatin Accessibility by Etoposide Cross-linking. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.2rngd5e>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** May 09, 2019

**Last Modified:** September 05, 2019

**Protocol Integer ID:** 23054

**Keywords:** TOP2, DNA accessibility, chromatin

### Abstract

**The best way to determine TOP2 specific accessibility to chromatin is via etoposide treatment, which leads to a covalent cross-link between TOP2 proteins and DNA, that can only occur where TOP2 proteins have access to DNA. TOP2 cross-linking to DNA has previously been shown to require DNA accessibility and chromatin regulators that alter DNA accessibility (Dykhuizen and Hargreaves et al. Nature 2013, Miller et al. NSMB 2017).**



## Guidelines

### Buffer A

- 25mM HEPES pH 7.0
- 25mM KCl
- 0.05mM EDTA
- 5mM MgCl<sub>2</sub>
- 10% Glycerol
- 0.1% NP-40

### High Salt RIPA

- 50mM TRIS pH 7.8
- 500mM NaCl
- 1% NP-40
- 0.1% DOC
- 0.1% SDS

### 1% SDS RIPA

- 50mM TRIS pH 7.8
- 150mM NaCl
- 1% NP-40
- 0.1% DOC
- 1% SDS

### Protease Inhibitors (1000x stock in DMSO all from Calbiochem)

- Chymostatin #230790 (10mg/ml)
- Leupeptin, Hemisulfate #108975 (10mg/ml)
- Pepstatin A, Synthetic #516481 (10mg/ml)



- 1 Cell Culture: 15m  
Trypsinize and count cells  
Transfer 0.75e6-1e6 cells to (4) different eppendorfs  
Save (1) eppendorf for your input sample for step 4  
Bring volumes to 1ml with media (If volume is already more than 1ml you can transfer to a 15ml conical and bring volume up to known amount where all cell lines are at the same volume)
- 2 Etoposide Cross-linking 1h  
Add Etoposide (1:1000; final concentration of 100uM) for 0, 15 and 60 minutes  
Add to 60' sample first, and then 15' sample 45' later  
Rotate at RT for 60'
- 3 Washes 10m  
Spin at 500g for 4 minutes at 4°C  
Wash with PBS  
Spin at 500g for 4 minutes at 4°C  
Remove supernatant
- 4 Buffer A -- Making Nuclei 15m  
For all samples + Input  
Resuspend in 1 ml Buffer A + Protease Inhibitors (1:1000) + DTT (1:1000; final concentration:)  
Incubate on ice for 10'  
Spin at 300g for 4 minutes  
Aspirate out supernatant  
Spin briefly and remove remaining supernatant
- 5 High salt RIPA extraction 25m  
Resuspend nuclei in 50-100ul High salt RIPA with 500mM NaCl  
Incubate on ice for 20'  
Spin 12,000g for 5' at 4°C  
Remove supernatant (can save if you want)  
Spin briefly and remove remaining supernatant
- 6 Solubalize chromatin pellet 33m  
Resuspend in 50ul 1% SDS RIPA + PIs + DTT + Benzonase (1:200)  
Be sure to fully resuspend pellet by pipetting up and down ~20 times  
Incubate at room temp for 10'  
Incubate on ice for 20'  
Spin 12,000g for 3 minutes and transfer supernatant to a new tube (there should not be much if any of a visible chromatin pellet after spinning. If there is then try resuspending in a larger volume)
- 7 Proceed to Western Blot adding equal volumes of each sample

