




Feb 01, 2022

Version 2

Total Lipid Extraction from Baker's yeast (*Saccharomyces cerevisiae*) V.2

DOI

dx.doi.org/10.17504/protocols.io.b4h8qt9w



Israel Olayide¹, Monica Rieth²

¹St. Louis University Dept. of Chemistry; ²Southern Illinois University-Edwardsville

Labyrieth



Monica Rieth

Southern Illinois University-Edwardsville

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b4h8qt9w>

Protocol Citation: Israel Olayide, Monica Rieth 2022. Total Lipid Extraction from Baker's yeast (*Saccharomyces cerevisiae*). [protocols.io https://dx.doi.org/10.17504/protocols.io.b4h8qt9w](https://dx.doi.org/10.17504/protocols.io.b4h8qt9w) Version created by **[Monica Rieth](#)**



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 01, 2022

Last Modified: February 01, 2022

Protocol Integer ID: 57632

Keywords: lipids, extraction, extracting total lipid, total lipid extraction, total lipids from baker, saccharomyce, yeast,

Abstract

This protocol outlines a method for extracting total lipids from Baker's yeast, *Saccharomyces cerevisiae*. It has been adapted from Roy et al., *J. Lipid Res.* **2018**. doi: 10.1194/jlr.M088559.

Attachments



[Extraction of lipids...](#)

35KB

Guidelines

All culture ODs should be normalized for comparative analyses prior to the extraction. Weighed Eppendorf tubes to be used before getting started.

Extracted lipids can be stored at -20 degC

Troubleshooting

Safety warnings

! chloroform and methanol should be used in a chemical hood



Day 1-3

3d

- 1 Sterilize the inoculating loop by dipping it in the ethanol and heating it for 30 second in the flame
- 1.1 Allow loop to cool the loop at least for 10 seconds
- 2 Streak plate from prepared glycerol stock (-80 °C) and spread on YPD agar plate
- 3 Allow plate to incubate for 2-3 days at 25-30 °C until colonies are 1-2mm in diameter

Day 4

1d

- 4 Inoculate 5 mL YPD with yeast colony from plate
- 5 Grow overnight at 30°C with shaking at 250 rpm

Day 5

1d

- 6 Check OD₆₀₀ at UV-Vis spectrophotometer
- 7 Harvest the whole cells when OD₆₀₀ is 1.5
Centrifuge at 1000 x g for 2 minutes or until clear to pellet the yeast
- 8 Pour off the supernatant without disturbing cell pellet
- 9 Wash the cells three times by resuspending pellets in 1ml PBS, centrifuge for at 500 x g for 3-5 mins.
- 10 After washing three times with PBS, resuspend in 1ml PBS and transfer this mixture to a 15 ml conical centrifuge tube.



- 11 Add 3.75 mL of chloroform/methanol (1:2) to the 1 mL of resuspended cells
Add 1.25 mL chloroform and 1.25 mL sterile water subsequently.
- 12 Vortex vigorously for 5 mins, and centrifuge at 150 x g for 5 mins at room temperature.
- 13 After centrifugation, a two-phase system is obtained: aqueous top phase and organic bottom phase which contains the lipids.
- 14 Carefully remove the top aqueous layer and the middle insoluble layers (precipitated proteins).
- 15 The organic bottom layer is dried using a speed vacuum or dried under a steady stream of nitrogen.
- 16 Record the weight of the dried sample by pre-weighing an Eppendorf tube or equivalent before or after drying the lipid.