

Apr 05, 2023

Version 1

Total nucleic acid extraction - NucleoMag® DNA/RNA Water Kit (MACHEREY-NAGEL Inc.) V.1

 Forked from [Total nucleic acid extraction - Maxwell\(R\) HT Environmental TNA Kit, custom \(Promega\)](#).

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7yjq9gwz/v1>

Adélaïde guet¹, Melissa Schussman²

¹UWM; ²University of Wisconsin - Milwaukee

mclellan lab



Melissa Schussman

University of Wisconsin - Milwaukee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7yjq9gwz/v1>

Protocol Citation: Adélaïde guet, Melissa Schussman 2023. Total nucleic acid extraction - NucleoMag® DNA/RNA Water Kit (MACHEREY-NAGEL Inc.). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7yjq9gwz/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 05, 2023

Last Modified: April 05, 2023

Protocol Integer ID: 80055

Keywords: wastewater, SARS-CoV-2, total nucleic extraction, RNA, DNA, extraction, purification, Promega, KingFisher, total nucleic acid extraction from wastewater, total nucleic acid extraction, rna water kit, dna, rna, ht environmental tna kit, nagel

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Total nucleic acid extraction from wastewater using Maxwell(R) HT Environmental TNA Kit, custom (Promega)

Guidelines

When work is completed, remove equipment and supplies from the cabinet. Wipe the work area with 10% bleach, let stand for 10 min, rinse with water, then with 70% ethanol, and finally with RNAase AWAY.

Materials

MATERIALS

- Ethanol USP/ACS or molecular biology grade (100%)
- Molecular biology grade water
- Isopropanol** molecular biology grade (100%)
- 6x** KingFisher **96-well plates** (Cat. no.: 95040460)
- 1x** KingFisher **96-tip comb well plate**(Cat. no.: 97002534)
- Screw cap** microcentrifuge tubes

Equipment

Mini-Beadbeater-16

NAME

high-energy cell disrupter

TYPE

BioSpec

BRAND

607

SKU

<https://biospec.com/product/mini-beadbeater-16>

LINK

1 speed

SPECIFICATIONS

Equipment

Kingfisher Flex

NAME

Automated Extraction System

TYPE

ThermoFisher

BRAND

5400630

SKU



Promega_Maxwell_HT_RNA_Waste...



Before start

1. Clean the working area and all equipment: wipe down with 10% bleach and let dry. Wipe down with 70% ethanol and let dry. Then, wipe down using RNase AWAY and let dry.
2. Prepare the 50% ethanol solution (it must be fresh!)
3. Prepare the 6 purification plates:
 - **Wash 1 plate:** Add 100 μ l of 50% ethanol and 900 μ l of wash buffer (WBA) to each well required for purification.
 - **Wash 2 plate** (same as plate Wash 1): Add 100 μ l of 50% ethanol and 900 μ l of wash buffer (WBA) to each well required for purification.
 - **Ethanol Wash plate:** Add 450 μ l of 50% ethanol to each well required for purification.
 - **Elution plate:** Add 100 μ l of 25 mM Tris-HCl (pH 8.0) to each well required for purification.
 - **Tip plate:** Place KingFisher 96-tip comb into an empty KingFisher 96-well plate. *While opening the 96-tip comb plate, pay attention to not touch the tips.*
 - **Lysis and Bind plate:** Add 35 μ L of Resin to each well required for purification (*vortex the bottle at max speed before use*). Add 50 μ l of Alkaline Protease Solution custom (APA) to each well required for purification. Add 250 μ l of cell lysis solution (CLD) to each well required for purification. Add 400 μ l of Isopropanol (100%) to each well required for purification.



Total nucleic acid extraction

2h

- 1 For **HA filter** extraction, let the sample thaw on ice and go to **step 2**.

5m

For **BCoV/BRSV** extraction (in duplicate), add 5 μ L of BCoV/BRSV solution to the 2-mL tube containing 475 μ L MWA1. Vortex for 15 seconds (speed 7 out of 10) and flash freeze the tube. Go to step 4.

For **Direct extraction**, add 150 μ L of wastewater to the 2-mL tube containing 100 μ L CTAB. Vortex for 15 seconds (speed 7 out of 10) and flash freeze the tube. Go to step 4.

2

For **HA filter** extraction, place the 2-mL tubes in the bead beater.

Equipment

Mini-Beadbeater-16

NAME

high-energy cell disrupter

TYPE

BioSpec

BRAND

607

SKU

<https://biospec.com/product/mini-beadbeater-16>

LINK

1 speed

SPECIFICATIONS

- 2.1 Bead beat for  00:02:30

2m 30s

Safety information

Start the bead beating when the beads start to be loose in the tubes.

- 2.2 Cooldown the samples on ice for  00:05:00 .

5m



2.3 Repeat Steps 9.1 and 9.2 once  .

7m 30s

3 Centrifuge at maximum speed for 1 min at room temperature.

1m

 150000 rpm, Room temperature, 00:01:00

4 For **HA filter** extraction, transfer 400 µL of supernatant to the **Lysis/Bind plate**.
For **BCoV/BRSV/Direct extraction**, transfer all supernatant to the **Lysis/Bind plate**.

10m

5 Start the protocol MN_96_Flex.bdz on the KingFisher Flex  01:14:00

1h 14m

Equipment

Kingfisher Flex

NAME

Automated Extraction System

TYPE

ThermoFisher

BRAND

5400630

SKU

6 Transfer the purified sample from the **Elution plate** to the **microcentrifuge tubes**.

10m

Note

The DNA/RNA is now ready for downstream applications. RNA extract may be stored in RNase-free water at -80°C for 1 year.

RT-ddPCR

7 Quantification by Droplet Digital PCR (ddPCR)



dx.doi.org/10.17504/protocols.io.bpg6mjze