




Feb 20, 2023

Version 11

Total RNA and DNA from Microalgae (12 samples per microplate) V.11

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvro85bvmk/v11>



Ying-Yu Hu¹, Zoe V. Finkel¹

¹Dalhousie University

Marine Microbial Macroecology Lab
Tech. support email: ruby.hu@dal.ca



Ying-Yu Hu

Dalhousie University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvro85bvmk/v11>

Protocol Citation: Ying-Yu Hu, Zoe V. Finkel 2023. Total RNA and DNA from Microalgae (12 samples per microplate) . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.6qpvro85bvmk/v11>Version created by **Ying-Yu Hu**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 19, 2022

Last Modified: February 20, 2023

Protocol Integer ID: 71512

Keywords: RNA, DNA, SYBR Green II, DNase, RNase, microalgae, fluorescence, dna from microalgae, microalgae sample, quantifying bulk rna, microalgae, total rna, bulk rna, ng rna, samples per microplate, rna, ng dna, dna, fluorochrome sybr green ii, microplate

Funders Acknowledgements:

Simons Collaborative on Ocean Processes and Ecology



Grant ID: 723789

Abstract

Here we describe a protocol for extracting and quantifying bulk RNA and DNA from microalgae, which is adapted from Berdalet E. et al. (2005).

RNA and DNA are extracted from microalgae samples and then quantified by fluorochrome SYBR Green II.

The level of sensitivity of this method is set at ca. 20 ~300 ng RNA and 10 ~ 100 ng DNA (ml assay)⁻¹.



Citation

Berdalet E, Roldán C, Olivar MP, Lysnes K (2026)
. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part A. Optimisation of the assay. Scientia Marina.

<https://doi.org/10.3989/scimar.2005.69n11>

LINK

Citation

Berdalet E, Roldán C, Olivar MP (2026)
. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part B. Quantification in natural samples. Scientia Marina.

<https://doi.org/10.3989/scimar.2005.69n117>

LINK

Guidelines

Estimation of RNA/DNA in the collected microalgae samples:

Under replete condition, RNA and DNA is about 5.7% and 1% in total dry mass, while Chl-a is about 1.1% in total dry mass. Therefore, $\text{RNA}_{\text{ug/L}} = \text{Chl-a}_{\text{ug/L}} \times (5.7/1.1)$, $\text{DNA}_{\text{ug/L}} = \text{Chl-a}_{\text{ug/L}} \times (1/1.1)$.

Common dilution from sample collected on the filter to assay is 1/40.

Materials

STEP MATERIALS

- ☒ Ribonuclease A from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R6513-50MG**
- ☒ DEOXYRIBONUCLEASE1 RNase and Protease Free **Bioshop Catalog # DRB002.10**
- ☒ Magnesium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #63069-100ML**
- ☒ Calcium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #21115-100ML**
- ☒ SYBR[®]; Green II RNA Gel Stain, 10,000X concentrate in DMSO **Thermo Fisher Catalog #S7564**
- ☒ Tris(hydroxymethyl)aminomethane hydrochloride 1M pH 8.0 RNase free **Fisher Scientific Catalog #AAJ60080AK**
- ☒ Deoxyribonucleic acid from calf thymus **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4522-1MG**
- ☒ N-Lauroylsarcosine sodium salt solution (20% RNase/DNase free) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L744-50mL**
- ☒ EDTA buffer (0.5M DNase/RNase free) **Bioshop Catalog #EDT333.100**
- ☒ UltraPure[™] DNase/RNase-Free Distilled Water **ThermoFisher Catalog #10977023**
- ☒ E. coli Total RNA **Thermo Fisher Scientific Catalog #AM7940**

Protocol materials

- ☒ Calcium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #21115-100ML**
- ☒ SYBR[®]; Green II RNA Gel Stain, 10,000X concentrate in DMSO **Thermo Fisher Catalog #S7564**
- ☒ Deoxyribonucleic acid from calf thymus **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4522-1MG**
- ☒ N-Lauroylsarcosine sodium salt solution (20% RNase/DNase free) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #L744-50mL**
- ☒ EDTA buffer (0.5M DNase/RNase free) **Bioshop Catalog #EDT333.100**
- ☒ UltraPure[™] DNase/RNase-Free Distilled Water **ThermoFisher Catalog #10977023**
- ☒ Ribonuclease A from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R6513-50MG**
- ☒ DEOXYRIBONUCLEASE1 RNase and Protease Free **Bioshop Catalog # DRB002.10**
- ☒ E. coli Total RNA **Thermo Fisher Scientific Catalog #AM7940**
- ☒ Magnesium chloride solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #63069-100ML**
- ☒ Tris(hydroxymethyl)aminomethane hydrochloride 1M pH 8.0 RNase free **Fisher Scientific Catalog #AAJ60080AK**

Safety warnings

- ! No data is available addressing the mutagenicity or toxicity of SYBR[®] Green II Nucleic Acid Gel Stain. Because this reagent binds to nucleic acids, it should be treated as a potential mutagen and used with appropriate care. The DMSO stock solution should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. We strongly recommend using double gloves when handling the DMSO stock solution. As with all nucleic acid stains, solutions of SYBR Green II Nucleic Acid Gel Stain should be poured through activated charcoal before disposal or collected in waste container to be treated later. The charcoal must then be incinerated to destroy the dye.

Day 1: Freeze-dry samples

- 1 Freeze dry samples and blank filters. Freeze at  -80 °C until processed.

Note

1. Freeze-drying should be as short as possible to reduce sample degradation.
2. The exact duration of freeze-drying depends on size of filter, quantity of sample and the size of container.

Equipment

FreeZone® 2.5 L Benchtop Freeze Dryers	NAME
Labconco®	BRAND
700202000	SKU

Day 1: Prepare primary solutions

- 2 Turn on UV light in biosafety cabinet for  00:15:00

- 3 Clean working surface with decontamination solution.

- 4 Prepare Tris buffer  5 mM  8.0

- 4.1 Pour  1 M  8.0 Tris into an RNase free 15 mL Falcon tube

 Tris(hydroxymethyl)aminomethane hydrochloride 1M pH 8.0 RNase free **Fisher Scientific Catalog #AAJ60080AK**

Equipment

Falcon® Centrifuge Tubes	NAME
Polypropylene, Sterile, 15 mL	TYPE
Corning®	BRAND
352096	SKU

- 4.2 Directly add  2.5 mL  1 M  8.0 Tris into 500 mL RNase free water in its original package.



⊗ UltraPure™ DNase/RNase-Free Distilled Water **Thermofisher Catalog #10977023**

Equipment

BT Barrier Pipet Tips	NAME
Pre-Sterile	TYPE
Neptune®	BRAND
BT1250, BT100, BT10	SKU

5 RNA primary standard solution (1M 200 ug/ml)

5.1 In the original package, the **frozen** E. Coli Total RNA is of 1 mg/mL, in which total RNA is 200 ug.

⊗ E. coli Total RNA **Thermo Fisher Scientific Catalog #AM7940**

Note

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/sp_7940.pdf

5.2 Uncap the original package of E. Coli Total RNA and directly add $800\text{ }\mu\text{L}$ Tris buffer (1M 5 mM , pH 8.0).

Cap the package and vortex for a thorough mix.

5.3 Aliquot 30 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at $-80\text{ }^{\circ}\text{C}$

Equipment

Finnpipette Stepper Pipette	NAME
Thermo Scientific™	BRAND
4540000	SKU



Equipment

Finntip stepper pipette tips

NAME

500 ul (sterile)

TYPE

Thermo Scientific

BRAND

Thermo Scientific™ 9404173

SKU

6 DNA primary standard solution (≈ 1M 500 ug/ml)

6.1 Uncap the original package of Deoxyribonucleic acid from calf thymus and add

 2 mL Tris buffer (1M 5 mM , pH 8.0). Deoxyribonucleic acid from calf thymus **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4522-1MG**

Note

https://www.sigmaaldrich.com/content/dam/sigma-aldrich/docs/Sigma/Product_Information_Sheet/d4522pis.pdf

6.2 Cap the package. Do not vortex or sonicate.

6.3 Keep the solution at $0\text{ }^{\circ}\text{C}$ ~ $4\text{ }^{\circ}\text{C}$ overnight to completely solubilize the DNA.
Gentle reversion is recommended.6.4 Aliquot 10 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at $-80\text{ }^{\circ}\text{C}$

Equipment

Finntip stepper pipette tips

NAME

500 ul (sterile)

TYPE

Thermo Scientific

BRAND

Thermo Scientific™ 9404173

SKU

7 RNase primary stock solution (1M 10 mg/ml)

7.1 Uncap the original package of Ribonuclease A from bovin pancreas and add  5 mL

Tris buffer ( 5 mM ,  8.0).

Cap the package and vortex for a thorough mix.

 Ribonuclease A from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R6513-50MG**

7.2 Aliquot 30 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at

 -20 °C

Equipment

Finntip stepper pipette tips

NAME

500 ul (sterile)

TYPE

Thermo Scientific

BRAND

Thermo Scientific™ 9404173

SKU

Equipment

Finntip™ Stepper Pipette Tips

NAME

500 ul (Sterile)

TYPE

Thermo Scientific

BRAND

21-377-149

SKU

8 DNase primary stock solution ( 5 mg/ml = 10,000 U/mL)

8.1 Uncap the original package of Deoxyribonuclease1 and add  1 mL Tris buffer (

 5 mM ,  8.0).

Cap the package and vortex for a thorough mix.

 DEOXYRIBONUCLEASE1 RNase and Protease Free **Bioshop Catalog # DRB002.10**

8.2 Aliquot 70 uL to 600 uL RNase free microtubes (every assay required 60 uL). Keep

frozen at  -20 °C .

Day 2: Exact RNA and DNA

9 Turn on UV light in biosafety cabinet for  00:15:00

10 Clean working surface with decontamination solution.

11 Prepare falcon tubes and tube rack in biosafety cabinet

	Volume of the tube (mL)	Content in the tube
	5	0.5 M EDTA
	5	20% sarcosine
	50	5 mM Tris
	15 or 50	1% STEB

Equipment

Falcon® Centrifuge Tubes

NAME

Polypropylene, Sterile, 15 mL

TYPE

Corning®

BRAND

352096

SKU

Equipment

Falcon® Centrifuge Tubes

NAME

Polypropylene, Sterile, 50 mL

TYPE

Corning®

BRAND

352070

SKU

12 Prepare STEB ([M] 1 %)

Note

Use the following formula to determine the total volume of 1% STEB required:
 (# samples + # blank filters) X (500 ul) + (500 ul) = total volume of 1% STEB required

12.1 Pour sarcosine ([M] 20 %) into an RNase free 5 mL falcon tube.

⊗ N-Lauroylsarcosine sodium salt solution (20% RNase/DNase free) **Merck**
MilliporeSigma (Sigma-Aldrich) Catalog #L744-50mL

12.2 Pour EDTA ([M] 0.5 M) into an RNase free 5 mL falcon tube.

⊗ EDTA buffer (0.5M DNase/RNase free) **Bioshop** Catalog #EDT333.100

12.3 Pour Tris buffer ([M] 5 mM , pH 8.0) into an RNase free 50 mL falcon tube.

12.4 Mix $\text{500 } \mu\text{L}$ sarcosine ([M] 20 %) , $\text{10 } \mu\text{L}$ EDTA ([M] 0.5 M) and 9 mL + $\text{490 } \mu\text{L}$ Tris buffer ([M] 5 mM , pH 8.0) to obtain STEB ([M] 1 %).

13 Prepare ice bath

14 Remove freeze-dried samples from -80°C freezer and place them On ice .

15 Add $\text{500 } \mu\text{L}$ Tris buffer ([M] 5 mM , pH 8.0) and $\text{500 } \mu\text{L}$ STEB ([M] 1 %) to the bead tube. Place tubes On ice .

Equipment

LYSING TUBES

MATRIX D 2 mL/15 mL

MP BIOMEDICALS

116913500/116933050

NAME

TYPE

BRAND

SKU

16 Rinse forceps by [M] 70 % volume ethanol and air dry.

Equipment

Filter forceps

blunt end, stainless steel

Millipore

XX6200006P

NAME

TYPE

BRAND

SKU

17 Transfer sample/blank filter into the bead tube by using clean forceps.

18 Invert immediately then put back On ice .

20s

19 Disrupt samples on the bead mill at 6.5 m/s.

30s

Equipment

Fastprep-24 5G™ Sample Preparation Instrument^{NAME}

MP Biomedicals

BRAND

116005500

SKU

20 Keep tubes  On ice . Check the label on each tube, restore the label if it fades. 1m 30s

21 Disrupt samples on the bead mill at 6.5 m/s. 30s

22 Keep tubes  On ice . Check the label on each tube, restore the label if it fades. 1m 30s

23 Disrupt samples on the bead mill at 6.5 m/s 30s

24 Keep tubes  On ice . Check the label on each tube, restore the label if it fades. 1m 30s

25 Disrupt samples on the bead mill at 6.5 m/s. 30s

26 Continuously shake homogenate in a multi-head vortex at the highest speed for  01:00:00  Room temperature 1h

Note

Vortex mixer should be able to remain stable on the bench under this vortex speed.

27 In the biosafety cabinet, transfer all homogenate into RNase free 2 mL micro-tube.

28 Freeze at  -80 °C until analyzed.

Day 3: Run the assay

29 Prepare ice bath.

30 Turn on UV light in biosafety cabinet for  00:15:00

31 Clean working surface with decontamination solution.

32 Prepare falcon tubes, microtubes and tube racks in biosafety cabinet

	Number of tubes	Type of tubes	Contents
	5	5 mL falcon tubes	1 M MgCl₂
			1 M CaCl ₂
			Working solution A (WS-A)
			Working solution B (WS-B)
			Working solution C (WS-C)
	1	50 mL falcon tube	5 mM Tris buffer
1		15 mL falcon tubes	0.05% STEB
	6	2 mL RNase free tubes	RNase working solution
			RNA secondary standard solution
			DNA tertiary standard solution
			900 mM MgCl ₂
			900 mM CaCl ₂
			Sybr green working solution (SG-II WS)
1		600 uL RNase free tube	DNA secondary standard
24		2 mL RNase free tubes	RNA standard solutions for RNA standard curves
			DNA standard solutions for DNA standard curves
N= total number of samples and blanks		2 mL RNase free tubes	Samples and blanks
3XN		2 mL RNase free tubes	Diluted samples and blanks
4		Microtube racks	Tubes of 2 mL in Set 1
			Tubes of 2 mL in Set A
			Tubes of 2 mL in Set B
			Tubes of 2 mL in Set C
1		Tube racks	Falcon tubes



Equipment

Screw-Cap Centrifuge Tube

NAME

5 mL

TYPE

VWR

BRAND

10002-738

SKU

Day 3: Run the assay (Caution: It is a long procedure!)

33 Organize and label the tubes as shown below

Set 1:

In microtube rack, label 2 mL tubes for samples and blanks to be further diluted.

Blk	1	2	3	4	5	6	7	8	9	10	11
-----	---	---	---	---	---	---	---	---	---	----	----

Set A, B and C:

In microtube rack, label 2 mL tubes for RNA (marked in pink), DNA (marked in blue) standard solutions and samples (marked in yellow)

Set A is for working solution A (WS-A) treatment, i.e. treated with DNase

Set B is for working solution B (WS-B) treatment, i.e. treated with RNase

Set C is for working solution A (WS-A) and C (WS-C) treatment, i.e. treated with DNase and RNase

Set A	R1A	R2A	R3A	R4A	R5A	D1A	D2A	D3A				
	BlkA	1A	2A	3A	4A	5A	6A	7A	8A	9A	10A	11A

Set B	R1B	R2B	R3B	R4B	R5B	D1B	D2B	D3B				
	BlkB	1B	2B	3B	4B	5B	6B	7B	8B	9B	10B	11B

Set C	R1C	R2C	R3C	R4C	R5C	D1C	D2C	D3C				
	BlkC	1C	2C	3C	4C	5C	6C	7C	8C	9C	10C	11C

34 Label tubes for reagents as following.

Follow the sheet, add Tris buffer (1M 5 mM , pH 8.0) to the reagent tubes:

Content	5 mM Tris (uL)
SG-II WS	1000+250
WS-A	2X1000+820
WS-B	2X1000+820
WS-C	2X1000+940
RNase	380
900 mM MgCl ₂	40
900 mM CaCl ₂	40
RNA secondary	990+495
DNA secondary	95
DNA tertiary	960
0.05% STEB	9X1000 + 500

35 Thaw Sybr green II at room temperature

 SYBR[®] Green II RNA Gel Stain, 10,000X concentrate in DMSO **Thermo**
Fisher Catalog #S7564

36 Add  900 μ L Tris buffer ( 5 mM ,  8.0) to each tube in **Set 1**

Note

Depending on the dilution of extracted sample

37 Follow the sheet, add Tris buffer ( 5 mM ,  8.0) to each tube in **Set A, B and C**.
The unit of volume is uL.

Set A	650	640	600	550	500	640	610	550				
	650	650	650	650	650	650	650	650	650	650	650	650

Set B	650	640	600	550	500	640	610	550				
	650	650	650	650	650	650	650	650	650	650	650	650

Set C	600	590	550	500	450	590	560	500				
	600	600	600	600	600	600	600	600	600	600	600	600

38 Prepare STEB ([M] 0.05 %)
Add  500 µL STEB ([M] 1 %) to 0.05% STEB tube, and vortex.

39 Add  250 µL STEB ([M] 0.05 %) to RNA and DNA standards in **Set A, B and C** by reverse pipetting.

Set A, B, C	250	250	250	250	250	250	250	250				

40 Place RNase and DNase primary stock solutions, RNA and DNA primary standard solutions and samples  On ice .

41 Turn on refrigerated centrifuge and set the temperature to  4 °C .

Equipment

CENTRIFUGE 5430 R

NAME

Eppendorf

BRAND

MP2231000510

SKU

42 Turn on shaker/incubator and set temperature to  37 °C .

Equipment

SHAKING INCUBATOR

NAME

71L

TYPE

Corning® LSE™

BRAND

6753

SKU

43 Prepare [M] 900 mM MgCl₂

43.1 Pour [M] 1 M MgCl₂ solution into 5 mL RNase free Falcon tube

 Magnesium chloride solution Merck MilliporeSigma (Sigma-Aldrich) Catalog #63069-100ML

43.2 Transfer  360 µL [M] 1 M MgCl₂ solution into 900 mM MgCl₂ tube

44 Add  60 μL  900 mM MgCl_2 to WS-A

45 Add  60 μL  900 mM MgCl_2 to WS-B

46 Prepare  900 mM CaCl_2

46.1 Pour  1 M CaCl_2 solution into 5 mL RNase free Falcon tube

 Calcium chloride solution Merck MilliporeSigma (Sigma-Aldrich) Catalog #21115-100ML

46.2 Transfer  360 μL  1 M CaCl_2 solution into 900 mM CaCl_2 tube

47 Add  60 μL  900 mM CaCl_2 to WS-A

48 Add  60 μL  900 mM CaCl_2 to WS-B

49 Prepare SG-II WS

49.1 Centrifuge one tube of SG-II concentrate at  Room temperature

 13000 rpm, 00:05:00 to deposit DMSO.

49.2 Wrap SG-II WS tube with foil, transfer **supernatant** of SYBR Green II 10,000X concentrate to SG-II WS tube in biosafety cabinet ( 8.75 μL per 1.25 mL Tris)

Note

Any step involving SYBR Green II should be operated in dark room or at least dim light. Prepare Sybr green II WS at this step to allow enough time for stabilization.

50

Note

Lunch break!

51 Prepare RNase working solution  0.5 mg/ml

Add  20 μL RNase primary stock solution ( 10 mg/ml) to RNase tube

52 Add  60 μL  0.5 mg/ml RNase to WS-B.

Keep WS-B  On ice .

53 Add  60 μL  0.5 mg/ml RNase to WS-C.

Keep WS-C  On ice .



54 Add  60 μL DNase primary stock solution ( 5 mg/ml) to WS-A.

Keep WS-A  On ice .

55 Centrifuge extracted samples  10000 x g, 4°C, 00:04:00

56 Prepare RNA secondary standard solution  2 $\mu\text{g}/\text{ml}$

Add  15 μL RNA primary standard solution to RNA standard tube and mix.

Keep  On ice .

57 Prepare DNA secondary standard solution ( 25 $\mu\text{g}/\text{ml}$)

Add  5 μL primary DNA standard solution ( 500 $\mu\text{g}/\text{ml}$) to DNA secondary tube and mix.

Note

DNA secondary standard will also be used to verify the actual concentration of DNA by using μdrop plate.

Note

Avoid vortexing DNA standard. Mix with pipette tip by aspirating up and down several times.

58 Prepare DNA tertiary standard solution  1 $\mu\text{g}/\text{ml}$

Add  40 μL DNA secondary solution ( 25 $\mu\text{g}/\text{ml}$) to DNA tertiary standard tube and mix.

Keep  On ice .

59 Load  50 μL WS-A to tubes in **Set A**.

Note

From 59 to 62: Reverse pipetting
Decontaminate pipet between different WS.

60 Load  50 μL WS-A to tubes **Set C**.

61 Load  50 μL WS-B to tubes in **Set B**.

62 Load  50 μL WS-C to tubes in **Set C**.

63 Add  100 μL centrifuged samples to its corresponding tubes in **Set 1**.

Gently invert the tube to mix sample.

- 64 From Set 1, transfer  250 μL of diluted samples to each corresponding tubes (marked in yellow) in **Set A, B and C**.
Add RNA secondary standard to tubes (marked in pink) in **Set A, B and C**.
Add DNA secondary standard to tubes (marked in blue) in **Set A, B and C**.
The unit of volume is μL .

Set A, B, C	0	10	50	100	150	10	40	100				
	250	250	250	250	250	250	250	250	250	250	250	250

Note

In order to avoid cross contamination from RNase or DNase, use one tip for each dispensing.
Pipette solution in the tube up and down for mixing.

- 65 Invert each tube to mix well and place all tubes into the shaker/incubator at  37 $^{\circ}\text{C}$,
continuously shaking at 200 RPM for  00:20:00 .

20m

Note

Incubation time is critical. Temperature might be disturbed by door open/close. Don't start the timer until temperature returns to 37 $^{\circ}\text{C}$.

- 66 After incubation, invert each tube for mixing and then place into the fridge to stop the reaction.

Day 3: Verify DNA concentration and SG-II absorbance

- 67 Measure DNA secondary concentration by using μdrop plate (sample volume: 4 μl)
Use Tris buffer ( 5 mM ,  8.0) as blank.

Equipment

$\mu\text{Drop}^{\text{TM}}$ Plates	NAME
Thermo Scientific	BRAND
N12391	SKU



Equipment

Varioskan LUX Multimode Microplate Reader NAMEThermo Fisher BRANDVL0L00D0 SKU

- 68 DNA_primary concentration ($\mu\text{g/ml}$) = $(\text{Abs}_{260} - \text{Abs}_{260}(\text{blank})) \times 50 \mu\text{g/ml} \times (10\text{mm}/0.5\text{mm}) \times \text{DF}$

Where, $\text{DF}=20$

Note

The DNA concentration is around 500 $\mu\text{g/mL}$ but can be much lower (since the small volume of DNA primary aliquot is hard to be mixed), use the measured DNA value to calculate the DNA primary concentration.

- 69 Check absorbance of SG-II WS:

- 69.1 In a transparent microplate, load
(1) 200 μL Tris buffer as blank
(2) 10 μL SG-II WS and 190 μL Tris buffer

- 69.2 Read absorbance at 480 nm, the value after subtracted by blank shall be no higher than 0.21

Day 3: Read fluorescence

- 70 Remove samples out of the fridge and allow to reach  Room temperature for  00:02:00 before loading the microplate.

Note

Since fluorescence decreases with increasing temperature, with percentage changes depending on the fluorophore (Bashford, 1987), the SG-II WS must be kept dark at RT (22°C) and the samples must be equilibrated at RT (c. 2 min).

- 71 Adhere black film on the top of a microplate lid.

Equipment

Black Vinyl Films for Fluorescence and Photoprotection	NAME
VWR	BRAND
89087-692	SKU

Equipment

Microplate Lids	NAME
Polystyrene	TYPE
Greiner Bio-One	BRAND
07000288	SKU

- 72 Load  10 µL SG-II WS to each well in the microplate with 0.5 mL tip of stepper, and cover the plate with the black-film lid.

Equipment

Finttip™ Stepper Pipette Tips	NAME
500 uL	TYPE
Thermo Scientific™	BRAND
9404170	SKU

Equipment

96-Well Black Microplates	NAME
Polystyrene	TYPE
Greiner Bio-One	BRAND
655076	SKU



Note

Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

- 73 Organize tubes in 96-well microtube rack in the same order as how microplates are loaded.

R1A	R1B	R1C	BlkA	BlkB	BlkC	4A	4B	4C	8A	8B	8C
R2A	R2B	R2C									
R3A	R3B	R3C	1A	1B	1C	5A	5B	5C	9A	9B	9C
R4A	R4B	R4C									
R5A	R5B	R5C	2A	2B	2C	6A	6B	6C	10A	10B	10C
D1A	D1B	D1C									
D2A	D2B	D2C	3A	3B	3C	7A	7B	7C	11A	11B	11C
D3A	D3B	D3C									

- 74 Load  190 μ L working sample to the microplate by reverse pipetting.

45m

	1	2	3	4	5	6	7	8	9	10	11	12
A	R1A	R1B	R1C	BlkA	BlkB	BlkC	4A	4B	4C	8A	8B	8C
B	R2A	R2B	R2C	BlkA	BlkB	BlkC	4A	4B	4C	8A	8B	8C
C	R3A	R3B	R3C	1A	1B	1C	5A	5B	5C	9A	9B	9C
D	R4A	R4B	R4C	1A	1B	1C	5A	5B	5C	9A	9B	9C
E	R5A	R5B	R5C	2A	2B	2C	6A	6B	6C	10A	10B	10C
F	D1A	D1B	D1C	2A	2B	2C	6A	6B	6C	10A	10B	10C
G	D2A	D2B	D2C	3A	3B	3C	7A	7B	7C	11A	11B	11C
H	D3A	D3B	D3C	3A	3B	3C	7A	7B	7C	11A	11B	11C

Note


Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

- 75 Shake black film covered microplate at  Room temperature for  01:30:00

1h 30m

Note

Read fluorescence right after 1h30 incubation at room temperature.



- 76 Setup microplate reader:
 Plate: Greiner F bottom chimney well PP 96 well;
 Shake: Continuous 5s at 600 rpm
 Bandwidth: 5 nm
 Endpoint reading: Ex 490 nm/Em 520 nm;

Equipment

Varioskan LUX Multimode Microplate Reader NAME
 Thermo Fisher BRAND
 VL0L00D0 SKU

- 77 Read fluorescence and export data to excel sheet.
- 78 In the fume hood, dispose any waste with SG-II into fluorescence stain waste container (some stain waste has DMSO solvent).

Calculate

- 79 RNA standard curve

79.1 Concentrations of RNA standards in the microplate

	RNA standard	Secondary 2 ug/mL (uL)	Tris (uL)	STEB (uL)	WS (uL)	Sample in microplate (uL)	SG-II (uL)	Conc in microplate (ng/mL)
	R1	0.00	650.00	250.00	50.00	190.00	10.00	0.00
	R2	10.00	640.00	250.00	50.00	190.00	10.00	20.00



	RNA standard	Secondary 2 ug/mL (uL)	Tris (uL)	STEB (uL)	WS (uL)	Sample in microplate (uL)	SG-II (uL)	Conc in microplate (ng/mL)
	R3	50.00	600.00	250.00	50.00	190.00	10.00	100.00
	R4	100.00	550.00	250.00	50.00	190.00	10.00	200.00
	R5	150.00	500.00	250.00	50.00	190.00	10.00	300.00

- 79.2 Slope of fluorescence in Set A vs concentration of RNA standard gives $m_{\text{RNA}+\text{DNase}}$ (≈ 0.03)
 Slope of fluorescence in Set B vs concentration of RNA standard gives $m_{\text{RNA}+\text{RNase}}$

- 79.3 Calculate ρ

$$\rho = \frac{m_{\text{RNA}+\text{RNase}}}{m_{\text{RNA}+\text{DNase}}}$$

- 80 Total RNA of the samples

$$\text{ug Total RNA (ml assay)}^{-1} = 0.001 \times \frac{(RFU_A - RFU_{A\text{Blank}}) - (RFU_C - RFU_{C\text{Blank}})}{(1 - \rho)} / m_{\text{RNA}+\text{DNase}}$$

Where,

RFU_A and RFU_C are the fluorescence in Tube A and Tube C of the same sample.

$RFU_{A\text{Blank}}$ and $RFU_{C\text{Blank}}$ are the fluorescence in Tube A and Tube C of the blank.

- 81 DNA standard curve

- 81.1 Concentrations of DNA standards in the microplate: Use measured DNA primary concentration instead of 500 ug/mL:

	DNA primary Conc (ug/mL)	DNA primary (uL)	Tris (uL)	Conc. DNA secondary (ug/mL)
		5	95	

	DNA secondary Conc. (ug/mL)	DNA secondary (uL)	Tris (uL)	Conc. DNA tertiary (ug/mL)
		40	960	

	DNA standard	DNA tertiary (uL)	Tris (uL)	STEB (uL)	WS (uL)	Sample in microplate (uL)	SG-II (uL)	Conc. in microplate (ng/mL)
	R1	0	650	250	50	190	10	0
	D1	10	640	250	50	190	10	~10
	D2	40	610	250	50	190	10	~40
	D3	100	550	250	50	190	10	~100



- 81.2 Slope of fluorescence in Set A vs concentration of DNA standard gives $m_{\text{DNA+DNase}}$
Slope of fluorescence in Set B vs concentration of DNA standard gives $m_{\text{DNA+RNase}}$
(≈ 0.12)

- 81.3 Calculate δ

$$\delta = \frac{m_{\text{DNA+DNase}}}{m_{\text{DNA+RNase}}}$$

- 82 Total DNA of the samples

$$\mu\text{g Total DNA (ml assay)}^{-1} = 0.001 \times \frac{(RFU_B - RFU_{B\text{Blank}}) - (RFU_C - RFU_{C\text{Blank}})}{(1 - \delta)} / m_{\text{DNA+RNase}}$$

Where,

RFU_B and RFU_C are the fluorescence in Tube B and Tube C of the same sample

$RFU_{B\text{Blank}}$ and $RFU_{C\text{Blank}}$ are the fluorescence in Tube B and Tube C of the blank.

- 83 Dilution factor=40

If,

- Sample is extracted by 1 mL extraction reagent
- In Set 1, sample is diluted to 1/10
- In Set 3, diluted by Tris and all working solutions to 250/950
- In microplate, diluted by SG-II WS to 190/200

Citations

Berdalet E, Roldán C, Olivar MP, Lysnes K. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part A. Optimisation of the assay

<https://doi.org/10.3989/scimar.2005.69n11>

Berdalet E, Roldán C, Olivar MP. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part B. Quantification in natural samples

<https://doi.org/10.3989/scimar.2005.69n117>