

Jan 28, 2022

Total RNA extraction from frozen placenta tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.b4b9qsr6>

Scott Lindsay-Hewett¹

¹University of California, San Diego

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Scott Lindsay-Hewett

University of California San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.b4b9qsr6>

Protocol Citation: Scott Lindsay-Hewett 2022. Total RNA extraction from frozen placenta tissue. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4b9qsr6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 26, 2022

Last Modified: January 28, 2022

Protocol Integer ID: 57441

Keywords: total rna extraction from frozen placenta tissue, quality total rna from frozen placenta tissue, mirvana mirna isolation kit from ambion, total rna extraction, mirvana mirna isolation kit, rna extraction, total rna, frozen placenta tissue, manual for the mirvana mirna isolation kit, rna, placenta, extraction, ambion

Abstract

This protocol describes the isolation of high-quality total RNA from frozen placenta tissue. Tissue is disrupted using a bead beater, and total RNA is isolated using the *mirVana* miRNA Isolation Kit from Ambion.

Written steps are adapted from Ambion's manual for the *mirVana* miRNA Isolation Kit and BioSpec's instructions for the Mini-BeadBeater-16.

Materials

Equipment	
Centrifuge	NAME
Benchtop Centrifuge	TYPE
Eppendorf	BRAND
5405000441	SKU
https://online-shop.eppendorf.us/US-en/Centrifugation-44533/Centrifuges-44534/Centrifuge-5425-PF-243560.html	LINK
Any benchtop centrifuge will suffice	SPECIFICATIONS

Equipment	
Block heater	NAME
--	BRAND
--	SKU

Equipment	
Vortex mixer	NAME
Any	BRAND
xx	SKU

Equipment	
Mini-Beadbeater-16	NAME
high-energy cell disrupter	TYPE
BioSpec	BRAND
607	SKU
https://biospec.com/product/mini-beadbeater-16	LINK
1 speed	SPECIFICATIONS

Equipment	
Fume hood	NAME
Fume hood	TYPE
Generic	BRAND
Unknown	SKU

Equipment

Set of micropipettes with rack: 100-1000 µl, 20-200 µl, 2-20 µl, and 0.5-10 µl ^{NAME}

Pipettor set ^{TYPE}

Pipetman ^{BRAND}

QP-1001-07 ^{SKU}

<https://www.minipcr.com/product/set-of-four-adjustable-volume-micropipettes-rack/> ^{LINK}

Can use equivalent Pipettors ^{SPECIFICATIONS}



Equipment

Bioanalyzer ^{NAME}

Bioanalyzer ^{TYPE}

Agilent ^{BRAND}

G2991AA ^{SKU}

<https://www.agilent.com/en/product/bioanalyzer-automated-electrophoresis/bioanalyzer-instrument/2100-bioanalyzer-instrument-228250> ^{LINK}

Any bioanalyzer will suffice. ^{SPECIFICATIONS}





✕ mirVanaTM miRNA isolation kit

✕ Ethanol Pure 200 proof for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023-500mL**

✕ Acid-Phenol:Chloroform, pH 4.5 (with IAA, 125:24:1) **Thermo Fisher Catalog #AM9720**

✕ Dry Ice

✕ DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021**

✕ Nuclease-free water, not DEPC-treated **Life Technologies Catalog #AM9932**

✕ Filter Tips

✕ Bioanalyzer RNA 6000 Nano Kit **Agilent Technologies Catalog #5067**

Preparation

15m

- 1 Clean workspace, pipettes, and gloves with RNaseZAP.
- 2 Prepare bucket of ice.
- 3 Heat nuclease-free water to 95°C .
- 4 Pre-cool microcentrifuge to 4°C .
- 5 Add 200 proof pure ethanol to *mirVana* wash buffers as instructed.
- 6 Pull frozen placenta samples (~150-200 mg each) from -80°C and place on dry ice until ready to process.

Note

Integrity of isolated total RNA will be greater if tissue was preserved in RNA*later* as soon as possible after harvesting.

Typically, tissue is submerged in RNA*later* overnight. The following day, the RNA*later* is removed and the sample moved to -80°C storage.

Cell lysis and tissue disruption

20m

- 7 Add approximately 1 mL of 1.0mm zirconia/silica beads to each frozen placenta sample, and $700\text{ }\mu\text{L}$ *mirVana* Lysis/Binding Buffer.

**Note**

After addition of beads and lysis buffer, tubes should be almost full. Exclude as much air as possible to reduce foaming.

Safety information

Samples should be in screw-cap microcentrifuge tubes with integral O-rings in the caps. Snap-cap tubes should not be used, unless secured with an adapter.

- 8 Load samples immediately into Mini-BeadBeater-16 vial holder ring. Up to 16 samples can be accommodated.

Safety information

IMPORTANT: Rotate the vial holder ring to a position where the small hole in the vial holder ring engages the anti-rotation pin sticking out of the wiggle mechanism. Slide the vial holder down the pin and seat it flat on the wiggle mechanism. Slide the large, black plastic hold-down cap over the stainless steel center bolt, aligning it so that it too slides down the anti-rotation pin. The hold-down cap must make contact with the top of the aluminum wiggle mechanism - not just the tops of the microcentrifuge tubes. Finally, screw on and hand-tighten the black knob firmly. Repeat: Tighten firmly.

- 9 Switch on the Mini-BeadBeater-16 and run for  00:02:00 .

2m

- 10 Place samples immediately into pre-chilled microcentrifuge and spin for  00:05:00 at maximum speed .

5m

- 11 Remove  500 μ L lysate, being careful not to draw up particulates, and dispense into a fresh labeled microcentrifuge tube.

Organic extraction

40m



12 Add  50 μL *mirVana* miRNA Homogenate Additive (1:10 volume of original lysate) and vortex to mix. Incubate  00:10:00  On ice .

10m

13 Add  500 μL Acid-Phenol:Chloroform (1:1 volume of original lysate) and vortex  00:01:00 to homogenize sample.

1m

Note

Acid-Phenol:Chloroform may appear as a clear, homogeneous phenol phase (lower), overlayed by a small aqueous phase (upper). Pipette from the lower, not the upper, phase.

Safety information

Phenol is very corrosive and will severely burn the skin. Safety precautions such as gloves, protective eyewear, a lab coat, and working in a fume hood are critical. Discard contaminated pipette tips in appropriate waste container.

14 Spin  00:05:00 maximum speed in pre-chilled microcentrifuge.

5m

15 Carefully remove  350 μL of the top aqueous phase and transfer to a fresh labeled microcentrifuge tube.

Note

It is possible to remove a greater volume to maximize RNA yield, but make sure not to disturb the interphase or organic phase.

Safety information

Discard phenol liquid waste and contaminated tubes in appropriate waste containers.

Total RNA isolation

30m



- 16 Add  437.5 μL (1.25 volumes) 200 proof pure ethanol and mix thoroughly.
- 17 Transfer up to  700 μL lysate/ethanol mixture to *mirVana* Filter Cartridge (placed into *mirVana* Collection Tube) and spin  10000 x g, 4°C, 00:00:15 in pre-chilled microcentrifuge. Discard flow-through. Repeat with remaining volume of lysate/ethanol mixture.

15s

Note

mirVana Filter Cartridges can accommodate up to  700 μL volume - do not overfill. To avoid filter damage, do not spin at speeds greater than  10000 x g.

- 18 Add  700 μL *mirVana* Wash Solution 1 and spin  10000 x g, 4°C, 00:00:15 in pre-chilled microcentrifuge. Discard flow-through.

15s

Note

Tip: flow-through can be aspirated into an appropriately-labeled waste container using a vacuum.

- 19 Add  500 μL *mirVana* Wash Solution 2/3 and spin  10000 x g, 4°C, 00:00:15 in pre-chilled microcentrifuge. Discard flow-through. Repeat wash step.

15s

- 20 After discarding flow-through from the last step, spin  10000 x g, 4°C, 00:02:00 in pre-chilled microcentrifuge to remove residual Wash Solution.

2m

- 21 Transfer *mirVana* Filter Cartridge into a fresh *mirVana* Collection Tube. Add  40 μL pre-heated nuclease-free water to the center of the filter and close the cap. Incubate  00:01:00 and then spin  00:00:30 maximum speed in pre-chilled centrifuge to elute RNA.

1m 30s



- 22 Transfer eluate from the *mirVana* Collection Tube to a low bind microcentrifuge tube and store at  -80 °C .

Note

Minimize freeze/thaw cycles.

Tip: Store a small aliquot in a separate tube for QC.

Quality control

1h 30m

- 23 Quantitate RNA using NanoDrop. Pay attention to A260/A280 ratio. For highly pure RNA, a ratio of 1.8-2.1 is expected. If necessary, repurify by adding 1/10th nuclease-free 5M NaCl and 1.38 volumes 200 proof pure ethanol before repassing the sample over a fresh *mirVana* Filter Cartridge. Continue the total RNA isolation from Step 17.

Note

RNA can also be quantitated using Qubit RNA Broad Range Assay.

- 24 Assess RNA quality by running the RNA 6000 Nano Assay on an Agilent 2100 Bioanalyzer. A RIN score of >7.0 is normal for total RNA isolated from placenta.