



Sep 17, 2025

Total RNA isolation from cells using RNAqueous-Micro Total RNA Isolation Kit

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg318zzl25/v1>

Peng Xu^{1,2}, Pietro De Camilli^{1,2}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Peng Xu

Yale University School of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg318zzl25/v1>

Protocol Citation: Peng Xu, Pietro De Camilli 2025. Total RNA isolation from cells using RNAqueous-Micro Total RNA Isolation Kit . protocols.io <https://dx.doi.org/10.17504/protocols.io.kqdg318zzl25/v1>

**Manuscript citation:**

Defect in hematopoiesis and embryonic lethality at midgestation of Vps13a/Vps13c double knockout mice
Peng Xu, Rubia Isler Mancuso, Marianna Leonzino, Caroline J. Zeiss, Diane S. Krause, Pietro De Camilli
bioRxiv 2025.05.09.653147; doi: <https://doi.org/10.1101/2025.05.09.653147>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 17, 2025

Last Modified: September 17, 2025

Protocol Integer ID: 227444

Keywords: micro total rna isolation kit, isolation of total rna, total rna isolation, total rna, small number of cell, isolation

Funders Acknowledgements:

National Institutes of Health

Grant ID: NS36251

National Institutes of Health

Grant ID: DA018343

Aligning Science Across Parkinson

Grant ID: ASAP-000580

Yale Cooperative Center of Excellence in Hematology

Grant ID: U54DK106857

Abstract

For the isolation of total RNA from a small number of cells

Sample disruption

- 1 Resuspend cell pellet (up to ~500,000 cells) by vortexing vigorously in at least 100 μ L Lysis Solution.

RNA isolation

- 2 For a standard prep of 100 μ L of lysate, add 50 μ L of 100% ethanol, and vortex briefly but thoroughly.
- 3 Load the lysate/ethanol mixture (up to 150 μ L) onto a Micro Filter Cartridge Assembly and close the cap.
- 4 Centrifuge for ~10 sec at maximum speed or until all of the mixture has passed through the filter. Longer centrifugation times may be needed to filter the lysate from tissue samples. For lysate/ethanol mixtures >150 μ L, load and filter the first 150 μ L, then repeat with additional aliquots until the entire sample has passed through the filter. The Collection Tube has a capacity of ~700 μ L when assembled with a Micro Filter Cartridge; if more than 150 μ L of lysate/ethanol mixture is filtered, empty the Collection Tube before proceeding. The RNA is now bound to the filter in the Micro Filter Cartridge.
- 5 Open the Micro Filter Cartridge, add 180 μ L of Wash Solution 1 (working solution mixed with ethanol) to the filter and close the cap. Centrifuge for ~10 sec to pass the solution through the filter.
- 6 Open the Micro Filter Cartridge, add 180 μ L of Wash Solution 2/3 (working solution mixed with ethanol) to the filter and close the cap. Centrifuge for ~10 sec to pass the solution through the filter. Repeat with a second 180 μ L aliquot of Wash Solution 2/3.
- 7 Open the Micro Filter Cartridge assembly, remove the filter cartridge from the Collection Tube, and pour out the flow-through. Replace the Micro Filter Cartridge into the same Collection Tube, close the cap, and centrifuge at maximum speed for 1 min to remove residual fluid and dry the filter.
- 8 Label a Micro Elution Tube (1.5 mL tubes provided with the kit) and transfer the Micro Filter Cartridge into it.
- 9 Apply 5–10 μ L of Elution Solution, preheated to 75°C, to the center of the filter. Close the cap and store the assembly for 1 min at room temperature. Centrifuge the assembly for ~30 sec to elute the RNA from the filter.
- 10 Repeat with a second 5–10 μ L aliquot of preheated Elution Solution, collecting the eluate in the same Micro Elution Tube.



Protocol references

https://documents.thermofisher.com/TFS-Assets/LSG/manuals/cms_057352.pdf