

Feb 19, 2019

## Tracking bleach synchronized worms

DOI

[dx.doi.org/10.17504/protocols.io.yb4fsqw](https://dx.doi.org/10.17504/protocols.io.yb4fsqw)

Priota Islam<sup>1</sup>

<sup>1</sup>Imperial College London

Behavioural Genomics



Priota Islam

Imperial College London

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yb4fsqw>

**Protocol Citation:** Priota Islam 2019. Tracking bleach synchronized worms. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.yb4fsqw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 19, 2019

**Last Modified:** February 19, 2019

**Protocol Integer ID:** 20572

**Keywords:** tracking bleach, worm

## Troubleshooting



## Preparing Imaging plates

- 1
  - Prepare and pour low peptone NGM onto 35mm imaging plates
  - Put them in the cold room for at least 2 days before use
  - Seed imaging plates with 50ul of freshly made 1:10 solution (OP50:M9 solution) the day before recording
  - Leave on bench to dry with lid on overnight

## Preparing worms

- 2
  - Chunk the worms on 3 maintenance plates
  - Bleach the worms after 2 days following the protocol for ***Bleach synchronization of c elegans***
  - Refeed the starved L1s on 2 different NGM plates that are OP50 seeded, 72hr prior to the experiment day (Repeat for every strain to be tracked)

## Day of recording

- 3
  - Using a hair pick – pick 5 young adults onto the seeded imaging plates
  - Put the plates on the tracker – agar side up with the lids off
  - Wait 30min for the worms to acclimatise
  - Record for 15min (See protocol ***Imaging on the multiworm tracker*** for detailed tracking instructions)