



Sep 12, 2025

Tragopogon (Asteraceae) transformation and regeneration

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwbo5zvmk/v1>

Shengchen Shan^{1,2}

¹Florida Museum of Natural History; ²University of Florida



Shengchen Shan

University of Florida

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwbo5zvmk/v1>

Protocol Citation: Shengchen Shan 2025. *Tragopogon* (Asteraceae) transformation and regeneration . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwbo5zvmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 02, 2025

Last Modified: September 12, 2025

Protocol Integer ID: 117586

Keywords: tetraploid tragopogon, genome duplication, polyploidy, asteraceae, immediate consequences of polyploidy, tragopogon, regeneration, regeneration protocol, whole genome, evolutionary model system, significant evolutionary force

Funders Acknowledgements:

National Science Foundation

Grant ID: IOS-1923234

National Science Foundation

Grant ID: DEB-2043478

Abstract

Polyploidy or whole-genome duplication is a significant evolutionary force. *Tragopogon* (Asteraceae) includes an evolutionary model system for studying the immediate consequences of polyploidy. Here, we provide an improved transformation and regeneration protocol for diploid and tetraploid *Tragopogon*.

Materials

	afds	hello	C	D

Tragopogon transformation and regeneration protocol

1 Prepare *Tragopogon* material

- 1.1 Remove the pericarp (fruit wall) of *Tragopogon* achenes to expose the seed; sterilize the seeds in 20% (v/v) bleach (containing a few drops of Tween-20) for 20 min (shake slowly). Rinse the seeds 4 times with autoclaved deionized water to remove excess bleach on the seed surface. Blot the seeds on the autoclaved filter paper.
- 1.2 Grow the seeds on **seed germination medium** (also known as **1/2 MS medium**) (2.22 g/L MS basal medium with vitamins, 30 g/L sucrose, 5 g/L agar, pH 5.8) in a sterile ice cream box in an incubator (25 °C, 14 h light/10 h dark, 52.5 $\mu\text{mol m}^{-2} \text{s}^{-1}$; model CU36L4C8 from Percival 26 Scientific, Perry, IA, USA).
- 1.3 Cotyledons and true leaves are harvested approximately 2 weeks after seed germination for *Agrobacterium*-mediated transformation (Step 3.1).

2 *Agrobacterium* culture

- 2.1 Streak the *Agrobacterium tumefaciens* strain EHA 105 onto a **selective LB plate** (40 g/L LB agar, 50 mg/L kanamycin, 50 mg/L rifampicin). Incubate the plate at 28 °C for 2 days.
- 2.2 Pick *Agrobacterium* colonies and inoculate into 20 mL **selective liquid TY medium** (5 g/L tryptone, 3 g/L yeast extract, 50 mg/L kanamycin, 50 mg/L rifampicin, 0.1 mM acetosyringone, pH 5.5). Shake at 100 rpm, 28 °C for 24 h.

3 *Tragopogon* transformation and regeneration

- 3.1 Prepare leaf explants. Chop leaf tissues from Step 1.3 into 1-cm segments. To promote transformation efficiency, the surfaces of the leaf segments are slightly scored to create wounds using a sterile blade. The *Tragopogon* leaf explants are placed on the **callus induction medium** [4.44 g/L MS basal medium with vitamins, 30 g/L sucrose, 0.5 mg/L naphthaleneacetic acid (NAA), 1 mg/L 6-benzylaminopurine (BAP), 6 g/L gelrite, pH 5.8] for 2 days (in the dark at room temperature, at 23 °C) before transformation (a process referred to as "pre-culture").
- 3.2 Prepare *Agrobacterium* suspension. Centrifuge TY medium from Step 2.2 at 2000 rpm for 10 min; remove the TY medium and re-suspend the *Agrobacterium* pellet with **AAM infection medium** (2.95 g/L KCl, 0.25 g/L $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$, 0.15 g/L $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$, 0.15 g/L

NaH₂PO₄·H₂O, 10 mg/L MnSO₄·4H₂O, 3 mg/L H₃BO₃, 2 mg/L ZnSO₄·7H₂O, 0.75 mg/L KI, 0.25 mg/L Na₂MoO₄·2H₂O, 0.025 mg/L CuSO₄·5H₂O, 0.025 mg/L CoCl₂·6H₂O, 0.876 g/L L-glutamine, 0.266 g/L L-aspartic acid, 0.174 g/L L-arginine, 7.52 mg/L glycine, 0.1 g/L myo-inositol, 0.5 mg/L nicotinic acid, 0.5 mg/L pyridoxine hydrochloride, 0.1 mg/L thiamine hydrochloride, 27.8 mg/L FeSO₄·7H₂O, 37.26 mg/L Na₂EDTA, 0.5 g/L casamino acids, 68.5 g/L sucrose, 36 g/L glucose, 0.1 mM acetosyringone, pH 5.2). For each transformation experiment, prepare 20 mL AAM infection medium with OD₆₀₀ = 0.1-0.2.

- 3.3 Co-culture of leaf explants and *Agrobacterium*. Immerse explants in 20 mL AAM infection medium for 30 min; explants are thoroughly blotted dry on autoclaved filter paper and transferred to **co-cultivation medium** (4.44 g/L MS basal medium with vitamins, 30 g/L sucrose, 0.5 mg/L NAA, 1 mg/L BAP, 6 g/L gelrite, 0.1 mM acetosyringone, pH 5.2). Place the plates in the dark for 3 days at room temperature (23 °C).
- 3.4 Transfer the explants to **selective callus induction medium** (4.44 g/L MS basal medium with vitamins, 30 g/L sucrose, 0.5 mg/L NAA, 1 mg/L BAP, 6 g/L gelrite, 300 mg/L timentin, 15 mg/L hygromycin, pH 5.8), and place the plates (covered by foil) inside the incubator (25 °C) for approximately 24 days. Callus will emerge from the cut edge of the explants. Successful transformants will survive on the selection medium and emit green fluorescence under blue light.
- 3.5 Shoot initiation. Calli of *Tragopogon porrifolius* (2x) and *T. mirus* (4x) are transferred to the **selective COCO-1** [2.41 g/L Woody plant basal medium with vitamins, 30 g/L sucrose, 2.5 mg/L BAP, 8.5 g/L agargellan, 10% (v/v) coconut water, 300 mg/L timentin, 15 mg/L hygromycin, pH 5.7] and **COCO-3** (identical to COCO-1, but with 15 g/L agargellan instead of 8.5 g/L) shooting media, respectively. Place the plates inside the same incubator for approximately 5 weeks. The medium is renewed every 2-3 weeks.
- 3.6 Shoot elongation. Calli with tiny shoots are transferred to **selective WPM-3 medium** (2.41 g/L Woody plant basal medium with vitamins, 30 g/L sucrose, 8.5 g/L agargellan, 300 mg/L timentin, 15 mg/L hygromycin, pH 5.7) for approximately 2 weeks.
- 3.7 Rooting. The elongated shoots are gradually transferred (up to approximately 4 weeks on selective WPM-3 medium) to **selective COCO-R medium** [2.41 g/L Woody plant basal medium with vitamins, 30 g/L sucrose, 1.5 mg/L indole-3-butyric acid (IBA), 8.5 g/L agargellan, 300 mg/L timentin, 15 mg/L hygromycin, pH 5.7] for rooting.
- 3.8 After approximately 6 weeks, rooted shoots are gradually transferred to soil and continue their growth inside a growth chamber (14 h light/10 h dark, 25 °C; model AR-1110 from Percival Scientific). After approximately three months of growth in soil, the regenerated plants are subjected to a three-month cold treatment (8 h light/16 h dark, 7 °C). The plants are then moved back to the growth chamber (14 h light/10 h dark, 25 °C), and after approximately 3-4 weeks, most plants begin to bolt and flower.