

Oct 19, 2023

Transfection by electroporation of GFP-LRRK2 and Immunofluorescent imaging of MEFs VPS35 (D620N) mutants stably expressing LysoTag

DOI

dx.doi.org/10.17504/protocols.io.n2bvj3dwxlk5/v1


Rotimi Y. Fasimoye¹, Dario R. Alessi¹

¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK



Francesca Tonelli

MRC-PPU at The University of Dundee



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj3dwxlk5/v1>

Protocol Citation: Rotimi Y. Fasimoye, Dario R. Alessi 2023. Transfection by electroporation of GFP-LRRK2 and Immunofluorescent imaging of MEFs VPS35 (D620N) mutants stably expressing LysoTag. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvj3dwxlk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 26, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 88587

Keywords: ASAPCRN, mouse embryonic fibroblast, lysotag transfection of foreign dna, expressing lysotag transfection, subcellular localisation of protein, embryonic fibroblast, immunofluorescent imaging of mefs vps35, lrrk2 into mef, colocalization of gfp, transfection by electroporation, subcellular localisation, immunofluorescent imaging, protein, microscopy, immunofluorescent, molecular biology, tagged lrrk2, expressed gfp, such as lrrk2, transfection, cell, foreign dna, primary cell, large sized protein, lysosomal, 250kd protein, low transfection efficiency, electroporation

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000463

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Transfection of foreign DNA and Immunofluorescent (IF) microscopy are powerful tools used in cellular and molecular biology to monitor subcellular localisation of proteins. Mouse embryonic fibroblasts (MEFs) are primary cells notorious for their low transfection efficiency by mostly available chemical methods. This efficiency become even lower if the aim is to express a large sized protein, such as LRRK2, a 250kD protein. Here, we described a method where we used electroporation to transiently transfect GFP-tagged LRRK2 into MEFs. We also used IF microscopy to visualise the subcellular localisation of the transiently expressed GFP-LRRK2. Furthermore, we investigated the colocalization of GFP-LRRK2 with a lysosomal localised TMEM192-3xHA.

Attachments



855-2211.docx

422KB

Materials

Materials

Cell lines

- Mouse Embryonic Fibroblast VPS35 WT (stably expressing TMEM192-3xHA)
- Mouse Embryonic Fibroblast VPS35 D620N (stably expressing TMEM192-3xHA)

Plasmids

- GFP-LRRK2 (DU13363). Plasmid available at MRCPPU depository at MRCPPUreagents@dundee.ac.uk

Antibodies

Table 1: List of primary antibodies

	Antibody	Company	Cat. number	Host species
	GFP	Abcam	AB13970	Chicken
	HA	Roche	47877600	Rat

Table 2: List of fluorophore-conjugated secondary antibodies

	A	B	C	D	E
	Antibody	Conjugated Fluorophore	Company	Cat. number	Host Species
	anti-Chicken	Alexa 488	Invitrogen	A11039	Goat
	anti-Rat	Alexa 594	Invitrogen	21209	Donkey

Media and Reagents

▪ Growth Media:

A	B
Dulbecco's Modified Eagle's Medium (DMEM) (GIBCO 11960-085)	
Foetal Bovine Serum (FBS) (Sigma F7524 Batch# BCBW6817)	10%

A	B
L-Glutamine (GIBCO 25030024)	1%
Penicillin-Streptomycin (GIBCO 15140122)	1%
NEAA (GIBCO 11140-035)	1X
Sodium Pyruvate (GIBCO 11360-039)	1mM

- Transfection media: Growth media without Penicillin-Streptomycin (Pen-Step)
- Dulbecco's phosphate-buffered saline (PBS) (GIBCO 14190169)
- Bovine Serum Albumin, BSA (Roche, 10735094001)
- Sodium Azide (Sigma, S2002).
- Hoechst 33342 solution (Thermo, 62249)
- VECTASHIELD antifading Mounting media (Vector Laboratories, H1000)

Equipment

- NEPA21 Super Electroporator with electroporator chamber (SONIDEL)
- 2mm gap Cuvette with individual pipette (SONIDEL, EC-002S)
- Incubator with FPI-sensor system and display controller MB1 (BINDER GmbH. Model: CB150. Power Output: 1.40kW, 230V, 6.1 Amp). This incubator has CO₂ and O₂ control.
- Leica TCS SP8 MP Multiphoton Microscope.
- Super Premium microscope slides (Frosted on one side) (VWR, 631-0114)
- Borosilicate Glass square coverslips (VWR, 631-0125)
- DeNovix® CellDrop Brightfield cell counter

Consumables

- 6-well tissue culture Petri Dishes (ThermoFisher. Catalog# 140675).
- Standard 1ml and 200 µl Pipette tips (Greiner bio-one. Catalog# 686271 and 685261 respectively).

☒ DMEM, high glucose, no glutamine **Thermo Fisher Catalog #11960085**

☒ L-Glutamine (200mM) **Thermo Fisher Scientific Catalog #25030024**

☒ Penicillin-Streptomycin (10,000 U/mL) **Thermo Fisher Scientific Catalog #15140122**

☒ MEM Non-Essential Amino Acids Solution (100X) **Thermo Fisher Catalog #11140035**

☒ DPBS no calcium no magnesium **Gibco - Thermo Fisher Scientific Catalog #14190169**



⊗ Bovine Serum Albumin Fraction V **Merck MilliporeSigma (Sigma-Aldrich) Catalog #10735094001**

⊗ Sodium azide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S2002**

⊗ Hoechst 33342 Solution (20 mM) **Thermo Fisher Catalog #62249**

⊗ VECTASHIELD® Antifade Mounting Medium **Vector Laboratories Catalog #H-1000**

⊗ Nunc® Cell-Culture Treated Multidishes, 6 well **Thermo Fisher Catalog #140675**

⊗ PIPETTE TIPS 100- 1000 µL BLUE SUITABLE FOR EPPENDORF STERILE 60 PIECES PER RACK **greiner bio-one Catalog #686271**

⊗ PIPETTE TIP 10 - 100 µL SUITABLE FOR EPPENDORF 96 PIECES / ST RACK **greiner bio-one Catalog #685261**

⊗ Anti-GFP antibody (ab13970) **Abcam Catalog #ab13970**







⊗ Goat anti-Chicken IgY (H L) Secondary Antibody, Alexa Fluor 488 **Thermo Fisher Scientific Catalog #A11039**

Troubleshooting



Transfection of cells with GFP-LRRK2 plasmid by electroporation










18h

- 1 Place coverslips in 6 well plate (one coverslip per well) and add  2 mL of media. Place the plate in an incubator. 
- 2 Pellet cells from 100% confluent 10cm plate and resuspend in  1 mL media.
- 3 Count cells and resuspend in media in a way that there are 30000-40000 cells per  10 μ L .
- 4 Power on Electroporator and plug in the electroporator chamber to the output socket.
- 5 Set the following Poring Pulse parameters:
 - Voltage: 200V
 - Length: 5ms
 - Interval: 50ms
 - Number of cycles: 2
 - Decay rate: 10%
 - Polarity: +
- 6 Set Transfer Pulse parameters:
 - Voltage: 20V
 - Length: 50ms
 - Interval: 50ms
 - Number of cycles: 5
 - Decay rate: 40%
 - Polarity: +/-
- 7 Add  3-4 μ g of plasmid into a 1.5ml Eppendorf tube. 

Note





Ensure this is not more than 10 μ l. If plasmid is too concentrated, add required amount and top-up to 10 μ l with media.





















- 8 Add  90 μL of resuspended cells into the tube. In total, cell number should be approximately 300000-400000. Pipette gently up and down 3 times to mix cell and plasmid. 
- 9 Transfer cell-plasmid mixture into 2mm gap Cuvette and close the cap. Add gently to the side to ensure there are no bubbles. If there are bubbles, tap the cuvette gently on the side until the bubbles move to the top.
- 10 Insert the cuvette into the electroporator chamber and close the lid.
- 11 Briefly measure electrical impedance by pressing the " Ω " button and check reading. This should be between 0.03 and 0.055. 
 - If reading is too low, add  1-2 μL of cells.
 - If too high, add  1-2 μL of media.
- 12 Press start button to begin electroporation process.
- 13 On completion, use Pasteur pipette (provided with the cuvette) to transfer cells into the wells containing a coverslip (from Step 1).
- 14 Rotate plate gently to spread the cells within the well.
- 15 Return plate into the incubator and incubate for at least  18:00:00 at  37 $^{\circ}\text{C}$. 
- 16 Replace media with full Growth media (i.e., media with Pen-Strip) and incubate for another 12-18 hours. 

Preparing cells for Immunofluorescence imaging

3h 25m

- 17 Remove media and wash cells. 
- 17.1 Remove media and wash cells with  3 mL PBS +0.2% BSA+0.02% sodium azide for  00:05:00 . (1/3) 



- 17.2 Wash cells with  3 mL PBS +0.2% BSA+0.02% sodium azide for  00:05:00 . (2/3) 5m
- 17.3 Wash cells with  3 mL PBS +0.2% BSA+0.02% sodium azide for  00:05:00 . (3/3) 5m
- 18 Fixed cells in 4% w/v PFA. Add  3 mL of dissolved PFA and incubate at  Room temperature for  00:10:00 . 10m  
- 19 Permeabilise cells with 1% NP-40 (v/v in PBS +0.2%BSA+0.02% sodium azide).
- 20 Block with 3% BSA (w/v in PBS) for  00:30:00 . 30m
- 21 Prepare a combination of primary antibodies (Table 1) as shown below. Antibodies are diluted in PBS +0.2% BSA+0.02% sodium azide.
- Rat anti-HA (1:1000) and Chicken anti-GFP (1:1000)
- 22 Incubate cells at  Room temperature with diluted primary antibodies for  01:00:00 . Do this in a humid chamber on a piece of Parafilm. Put a  60 μ L drop of diluted antibodies on the parafilm. Carefully place coverslip on the droplet, with the side containing attached cells, facing inward, making contact with the droplet. 1h 
- 23 Wash cells, 3 times, with  3 mL PBS +0.2%BSA+0.02% sodium azide. 
- 24 Prepare a combination of Secondary antibodies as described below (see Table 2 for more information about the secondary antibodies). Antibodies are diluted in PBS +0.2% BSA+0.02% sodium azide.
- anti-Rat Alexa 594 (1:500) and anti-Chicken Alexa 488 (1:500).
- 25 Add  0.5 μ L Hoechst 33342 solution for nuclear staining. 



- 26 Incubate cells at Room temperature with diluted secondary antibodies for 01:00:00 . Do this in a humid chamber on a piece of Parafilm. Put a 60 μL drop of diluted antibodies on the parafilm. Carefully place coverslip on the droplet, with the side containing attached cells, facing inward, making contact with the droplet.
- 27 Wash cells, 3 times, with 3 mL PBS +0.2%BSA+0.02% sodium azide.
- 28 Rinse cells by dipping briefly in MilliQ water and gently dry on Kleenex wipes.
- 29 Label microscope glass slides (preferably the one with frosted side) according to the primary antibody used. Take note of the emission wavelength of the probe on the secondary antibodies.
- 30 Add a drop of VECTASHIELD antifading Mounting media.
- 31 Mount cover slip (containing cells) on the glass slide, ensuring that the side containing the cells is facing inward, making contact with the oil. Allow to dry for 00:30:00 , ensuring slides are prevented from direct light.
- 32 Slides can be stored at 4 $^{\circ}\text{C}$ or viewed immediately on a confocal microscope.

1h



30m





Figure 1: Immunofluorescence images of mouse embryonic fibroblasts (MEFs) expressing GFP-LRRK2 and TMEM192-3xHA. MEFs VPS35 wildtype and D620N mutants stably expressing TMEM192-3xHA and transiently expressing GFP-LRRK2 were co-immunostained with anti-HA and anti-GFP antibodies. Scale bar is 2 μ m.