

Dec 18, 2019

## TRANSFECTION OF i<sup>3</sup>NEURONS (Support Protocol 3)

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.5wug7ew](https://dx.doi.org/10.17504/protocols.io.5wug7ew)

iPSCs

Michael S. Fernandopulle<sup>1</sup>, Ryan Prestil<sup>1</sup>, Christopher Grunseich<sup>1</sup>, Chao Wang<sup>2</sup>, Li Gan<sup>2</sup>, Michael E. Ward<sup>1</sup>

<sup>1</sup>National Institute of Neurological Disorders and Stroke, National Institutes of Health, Bethesda, Maryland;

<sup>2</sup>Gladstone Institute of Neurological Disease, Gladstone Institutes, San Francisco, California

Neurodegeneration Method Development Community

Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



Julia Rossmanith

protocols.io

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5wug7ew](https://dx.doi.org/10.17504/protocols.io.5wug7ew)

External link: <https://doi.org/10.1002/cpcb.51>

**Protocol Citation:** Michael S. Fernandopulle, Ryan Prestil, Christopher Grunseich, Chao Wang, Li Gan, Michael E. Ward 2019. TRANSFECTION OF i<sup>3</sup>NEURONS (Support Protocol 3) . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5wug7ew>

### Manuscript citation:

Fernandopulle, M. S., Prestil, R., Grunseich, C., Wang, C., Gan, L., & Ward, M. E. (2018). Transcription-factor mediated differentiation of human iPSCs into neurons. Current Protocols in Cell Biology, e51. doi:<https://doi.org/10.1002/cpcb.51>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 26, 2019

**Last Modified:** December 18, 2019

Protocol Integer ID: 26292

**Keywords:** i3LMN, i3Neurons, iPSC, iPSC-derived neurons, transcription factor-mediated differentiation

## Abstract

Transient protein expression can easily be studied in i<sup>3</sup>Neurons using lipid-based transfection. This protocol is identical to that in iPSCs (see **Basic Protocol 2**). i<sup>3</sup>Neurons are modestly transfectable, with 5 % to 10 % of cells showing fluorescent protein expression after 24 hr. We have found that refreshing neuronal medium 1 to 2 hr after transfection both allows successful DNA entry into cells and largely prevents cytotoxicity resulting from the transfection reagent. Unlike iPSCs, i<sup>3</sup>Neurons show increased protein expression/accumulation over time, with greater fluorescence 48 to 72 hr after transfection than at 24 hr. Transient transfections also show more durable expression in i<sup>3</sup>Neurons than iPSCs, likely because episomes are not diluted by cell division. i<sup>3</sup>Neurons can be transfected in suspension (i.e., re-plating after day 3 of differentiation) or as an adherent culture, although better results are observed in adherent cultures. They are also amenable to serial transfections (i.e., re-transfecting with the same construct 24 hr apart) if higher-percentage transfections are desired.

## Attachments



[fernandopulle2018.pd...](#)

1.7MB

## Guidelines

This protocol is identical to that in iPSCs (see **Basic Protocol 2**). i<sup>3</sup>Neurons are modestly transfectable, with 5 % to 10 % of cells showing fluorescent protein expression after 24 hr. We have found that refreshing neuronal medium 1 to 2 hr after transfection both allows successful DNA entry into cells and largely prevents cytotoxicity resulting from the transfection reagent. Unlike iPSCs, i<sup>3</sup>Neurons show increased protein expression/accumulation over time, with greater fluorescence 48 to 72 hr after transfection than at 24 hr. Transient transfections also show more durable expression in i<sup>3</sup>Neurons than iPSCs, likely because episomes are not diluted by cell division. i<sup>3</sup>Neurons can be transfected in suspension (i.e., re-plating after day 3 of differentiation) or as an adherent culture, although better results are observed in adherent cultures. They are also amenable to serial transfections (i.e., re-transfecting with the same construct 24 hr apart) if higher-percentage transfections are desired.

## Safety warnings

⚠ Please see SDS (Safety Data Sheet) for hazards and safety warnings.

