

Feb 13, 2019 Version 1

# Transfection of constructs in diplomonads to block NHEJ pathway using the inhibitors of KU70/80 proteins. V.1

DOI

[dx.doi.org/10.17504/protocols.io.x3zfqp6](https://dx.doi.org/10.17504/protocols.io.x3zfqp6)

Binnypreet Kaur<sup>1, 2</sup>, Drahomíra Faktorová<sup>1, 2</sup> and Julius Lukeš<sup>3, 2, 3</sup>

<sup>1</sup>Institute of Parasitology, Czech Academy of Sciences;

<sup>2</sup>Faculty of Sciences, University of South Bohemia, České Budějovice, Czech Republic;

<sup>3</sup>Institute of Parasitology, Czech Academy of Sciences<sup>3</sup>;

<sup>4</sup>Canadian Institute for Advanced Research, Toronto, Canada



**Binnypreet Kaur**

1 Institute of Parasitology, Biology Centre, Czech Academy o...

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.x3zfqp6](https://dx.doi.org/10.17504/protocols.io.x3zfqp6)

**Protocol Citation:** Binnypreet Kaur<sup>1, 2</sup>, Drahomíra Faktorová<sup>1, 2</sup> and Julius Lukeš<sup>1, 2, 3</sup> 2019. Transfection of constructs in diplomonads to block NHEJ pathway using the inhibitors of KU70/80 proteins.. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x3zfqp6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development

**We are still developing and optimizing this protocol**

**Created:** February 13, 2019

**Last Modified:** February 13, 2019

**Protocol Integer ID:** 20313

**Keywords:** Drahomíra Faktorová



## Abstract

The process of NHEJ is completed when the DNA ligase complex ligates the broken ends of the DNA. Ku70-Ku80 are the DNA-dependent protein kinase heterodimeric regulatory factor that forms a link between two broken DNA ends and structurally support, align and protect them from further degradation. W7 inhibitor (W7(N-(6-aminohexyl)-5-chloro-1-naphthalenesulfonamide) inhibits the cofactor of inositol hexakisphosphate which effects the regulatory factor of Ku heterodimers.

## Materials

Nucleofector™ 2b Device , Human T Cell Nucleofector™ Kit, Cell counter, Normal growth medium, Tissue culture plates, Microcentrifuge tubes

- 1 Step 0: Pre-incubation of Diplonema cells for 4 hours with W7 (5 mg/ml) before electroporation. EC value of W7 was calculated by Alamar blue assay
- 2 Count the cells and plan to 1 nucleofection with  $5 \times 10^7$  cells for each construct.
- 3 Harvest the cells by centrifugation at 1300xg for 5 min at room temperature in Swing Bucket Rotor
- 4 Resuspend the cell pellet in 100ul of AMAXA Human T- cell solution at 4C (from refrigerator combine 81.8ul of Human T-cell nucleofector solution + 18.2ul Supplement).
- 5 Add 5-10ug of (PCR) (linearized DNA) into the cuvette (resuspend in 10ul of H<sub>2</sub>O).
- 6 : Put everything into the cuvette, close the cap and place in the electroporator, cuvette should only fit in one direction, but metal sides should face towards you.
- 7 Press for the Program X-001 to electroporate.
- 8 Transfer the entire cell- DNA transfectant into 10ml of Diplonema media (chloramphenicol) with the supplied disposable micropipettes
- 9 Immediately observe the flasks containing transfectants under a microscope and place them in the incubator at 27C + shaker and let them grow for 6-8h.
- 10 After 6-8 hours took out the flasks and make controls which should be without antibiotics - put 1ml of the electroporated cells in the first well of 24 well plates.
- 11 In the rest - 9ml -transfectants add Puromycin antibiotics in the flasks.
- 12 Add 1.5ml of electroporated cells in the first row (6 wells) of 24 well plates. Put 900ul of the media (+ selectable drug) in the rest of the wells. Take 100ul from first well and transfer to the next to make 10x serial dilution and so on...up to the lowest concentration  $5 \times 10^2$  cells/ml and let them grow until selection of clones is done.



- 13 Result: Unfortunately, targeting to the planned position (N-terminal tagging of alpha-tubulin with mCherry under puromycin<sup>R</sup> selection) did not work in any of the obtained clones.