


Aug 01, 2024

Transfections for gene delivery and genome editing in *S. rosetta* (VERSION 4)

DOI

dx.doi.org/10.17504/protocols.io.j8nlk86o5l5r/v1



David Booth¹

¹University of California, San Francisco

Protist Research to Opti...



David Booth

University of California, San Francisco

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk86o5l5r/v1>

Protocol Citation: David Booth 2024. Transfections for gene delivery and genome editing in *S. rosetta* (VERSION 4). protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlk86o5l5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 01, 2024

Last Modified: August 01, 2024

Protocol Integer ID: 104401

Keywords: cas9 genome editing, execution of crispr, genome editing, crispr, transfections for gene delivery, gene delivery, purified cas9 ribonucleoprotein, guide rna, genome, rna, dna oligonucleotide, dna oligonucleotides as template, cas9 ribonucleoprotein from streptomyces pyogene, cas9 cleavage site, spcas9 rnp, dna, transfection, mutation, gene

Abstract

This protocol details the preparation and execution of CRISPR/Cas9 genome editing in *S. rosetta*. The protocol builds on a [method to transfect macromolecules into *S. rosetta*](#) for delivering a purified Cas9 ribonucleoprotein from *Streptomyces pyogenes* (*SpCas9* RNP) into *S. rosetta*. Upon cleaving the *S. rosetta* genome at locations specified by the guide RNA (gRNA) of the *SpCas9* RNP, *S. rosetta* can use DNA oligonucleotides as templates to repair the double-stranded break. Those repair templates can encode foreign sequences and mutations for editing the *S. rosetta* genome, so long as DNA oligonucleotides have >30 bases of sequence that is homologous to both sides of the Cas9 cleavage site.

Guidelines


Perform cell culturing and transfection procedure inside of a biosafety cabinet to maintain sterility.

Protocol materials

 Papain from papaya latex **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3125-100MG**

Troubleshooting



Before start

Please consult the attached file of media recipes for artificial seawater, high nutrient media, and low nutrient media.  20240731_MediaTable.xlsx 20KB





Culture Cells

1 Seed a large culture of *S. rosetta*.

- 1.1 Two days prior to transfection, inoculate  80 mL of media (15% Red Algae + 2% Peptone-Yeast-Glycerol) with a culture of *S. rosetta* feeding on *E. pacifica* (ATCC PRA-390) to a final concentration of *S. rosetta* of $[M] 1.5 \times 10^4$ Mass Percent. Also add  400 μ L of $[M] 10$ Mass Percent of *E. pacifica* to the flask to provide ample food for growth.

Note

Prepare the *E. pacifica* solution by resuspending a frozen, 10 mg pellet of *E. pacifica* in 1 ml of artificial seawater.


- 1.2 Grow the culture for  40:00:00 in a 300 cm² flask at  22 °C .

1d 16h

Prepare cargo for transgene expression from plasmids

2 Prepare expression plasmids for transfection

Note

Plasmids can be prepared in advance and stored at  -20 °C

- 2.1 Use standard DNA plasmid preparation protocols to purify plasmids from *E. coli* that lack adenine and cytidine methyl transferases (dam-/dcm-).
- 2.2 If needed, concentrate DNA plasmids to a final concentration of $[M] 5 \mu\text{g}/\mu\text{L}$ using a standard ethanol precipitation.

**Note**

If you have purified your plasmid with magnetic beads, be sure to thoroughly remove the magnetic beads by centrifuging them at maximum speed for 5 min. Only take the clarified supernatant.










Note**DNA precipitation with ethanol and a salt solution**

1. Add [M] 3 Molarity (M) Potassium Acetate $\text{pH } 5.4$ to purified plasmid to achieve a final concentration of at least [M] 300 millimolar (mM) Potassium Acetate .
2. Measure the volume of DNA and potassium acetate and then multiply this volume by 2.5x to calculate the volume of 100% ethanol to add to the DNA/Potassium Acetate mixture.
3. After adding the ethanol, thoroughly homogenize the solution. At this point, you may see a white precipitate of DNA.
4. Centrifuge the solution at $\text{21000} \times g$, 4°C , 00:30:00 .
5. Remove the supernatant and wash the pellet with $200 \mu\text{L}$ of [M] 70 % volume Ethanol, ice-cold.
6. Centrifuge the pellet at $\text{21000} \times g$, 4°C , 00:05:00 .
7. Remove the supernatant.
8. Resuspend the pellet in a minimal volume [M] 10 millimolar (mM) HEPES-KOH $\text{pH } 7.5$
9. Dilute the sample ten-fold in water and read the concentration with a spectrophotometer. $1 \mu\text{L} =$ [M] 0.05 $\mu\text{g}/\mu\text{L}$ of dsDNA
10. Adjust the concentration of the plasmid achieve a final concentration of [M] 5 $\mu\text{g}/\mu\text{L}$.

3 Prepare a mixture of carrier molecules to add to the purified plasmid.**Note**

The mixture can be prepared ahead of time and stored at -20°C . Note that the final mixture is very viscous and should be thoroughly homogenized by pipetting up and down before usage.








- 3.1 For one transfection, add  2 μL of  20 $\mu\text{g}/\mu\text{L}$ pUC19 plasmid that has been resuspended in  10 millimolar (mM) Tris-HCl  8.0
- 3.2 For one transfection, add  1 μL of  250 millimolar (mM) Adenosine Triphosphate-KOH  7.5 .
- 3.3 For one transfection, add  1 μL of  100 Mass Percent Heparin, Sodium Salt .
- 3.4 Thoroughly mix the viscous solution by pipetting up and down.

4 Combine expression plasmids and the carrier mixture for each transfection.

Note

This step can be performed ahead of time, but it is usually done the same day as the transfection.

- 4.1 For one transfection, place  4 μL of the carrier mixture ( [go to step #3](#)) in the bottom of a 1.5 ml conical tube.
- 4.2 Add  1 μL of  5 $\mu\text{g}/\mu\text{L}$ Expression Plasmid ( [go to step #2](#)) to the carrier mixture and slowly pipette up and down to thoroughly mix the solution. This solution is called the "Plasmid Delivery Mix."











Note






The expression of different transgenes may require higher or lower concentrations of plasmid. Performing a titration of plasmid for your particular assay may be necessary. Typically, a range of 1 to 10 μg of plasmid is recommended for one transfection.

Prepare cargo for genome editing with Cas9 RNP

5 **Prepare a guide RNA (gRNA) that binds to *SpCas9* and targets DNA by annealing CRISPR RNA (crRNA) with the trans-activating CRISPR RNA (tracrRNA) .**








- 5.1 Resuspend crRNA in duplex buffer ([M] 30 millimolar (mM) HEPES-KOH  7.5 ; [M] 100 millimolar (mM) Potassium Acetate) to a final concentration of [M] 200 micromolar (μM) .
- 5.2 Resuspend tracrRNA in duplex buffer to a final concentration of [M] 200 micromolar (μM) .
- 5.3 Mix equal volumes of crRNA ( go to step #5.1) and tracrRNA ( go to step #5.2) to have a final gRNA concentration of [M] 100 micromolar (μM) (gRNA is the annealed complex of crRNA and tracrRNA).
- 5.4 Incubate the gRNA solution at  95 °C in an aluminum block for  00:05:00 .
- 5.5 Place the aluminum block at  Room temperature to slowly cool the gRNA to  25 °C .
- 5.6 Store the gRNA at  -20 °C .

6 **Prepare DNA oligonucleotides that serve as repair templates after *SpCas9* cleavage.**

- 6.1 Dissolve oligonucleotides to a final concentration of [M] 250 micromolar (μM) in 10 mM HEPES-KOH, pH 7.5.
- 6.2 Incubate the dissolved oligonucleotides at  55 °C for  01:00:00 .
- 6.3 Store oligonucleotides at  -20 °C .
- 6.4 Before starting nucleofections, ensure that the oligonucleotides are fully dissolved by incubating them at  55 °C for  01:00:00 , which concurs while the *SpCas9*/gRNA complex assembles.





7 Assemble *SpCas9* with the gRNA to form the *SpCas9* RNP.

- 7.1 For one transfection, place  2 μL of  20 micromolar (μM) *SpCas9* in the bottom of a 0.2 ml PCR tube.
- 7.2 Add  2 μL of  100 micromolar (μM) gRNA ( go to step #5) by slowly pipetting up and down with *SpCas9* to gently mix the gRNA together. This solution is called the "*SpCas9* ribonucleoprotein (RNP)."
- 7.3 Incubate the *SpCas9* RNP at  Room temperature for  01:00:00 (roughly the time to complete the preparation of *S. rosetta* for priming, see below).

Prepare transfection Reagents

8 Prepare SF Buffer (Lonza) for transfections.

- 8.1 Add all of buffer B (smaller volume that may also be called supplement 1) to buffer A (larger volume). 
- 8.2  0 μL Store on ice until ready for use. The combined buffer can also be stored at 4°C for up to 3 months.

Note

The combined buffer can be stored at 4°C for up to 3 months.

Note

Because the Lonza kits are quite expensive, we recommend aliquoting large volumes of the SF components (900 μl aliquots for buffer A and 200 μl aliquots for buffer B) to prevent SF buffer from quickly spoiling after buffers A and B have been combined.

9 Prepare the priming buffer.



- 9.1 Dilute papain to a final concentration of [M] 100 micromolar (μM) in dilution buffer ([M] 50 millimolar (mM) HEPES-KOH pH 7.5 , [M] 200 millimolar (mM) Sodium Chloride , [M] 20 % volume Glycerol and [M] 10 millimolar (mM) Cysteine) from a stock solution of [M] 1 millimolar (mM) Papain (the recommended Papain is already at this concentration), and incubate at Room temperature just before priming cells for transfection.



Papain from papaya latex **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3125-100MG**

Note

The dilution buffer should be sterile filtered through a 0.22 μm filter. This buffer may also be prepared ahead of time and stored in a -80°C freezer until just before its use.

- 9.2 Prepare the remaining components of the priming buffer ([M] 40 millimolar (mM) HEPES-KOH pH 7.5 , [M] 34 millimolar (mM) Lithium Citrate , [M] 15 Mass / % volume PEG 8000 and [M] 50 millimolar (mM) Cysteine). DO NOT combine the papain and priming buffer until just before adding the priming buffer to cells.




Note

The priming buffer without papain should be sterile filtered through a 0.22 μm filter, which can also be made ahead of time and stored at -80°C until it is used. Just be sure that the priming before is warmed to room temperature prior to use.

Wash Cells


- 10 Prepare *S. rosetta* for transfection by washing away feeder bacteria.



10.1 Homogenize  80 mL culture of *S. rosetta* feeding on *E. pacifica* ( go to step #1) by vigorously shaking the flask and then splitting the culture into  40 mL aliquots in 50 ml conical tubes. Vigorously shake the tubes for 30 sec to further break up bacterial clumps and loosely associated *S. rosetta* cells

10.2 Centrifuge the cells at  2400 x g, 4°C, 00:03:00 in a swinging bucket rotor. 

10.3 Use a serological pipette to gently remove all but 2 ml of the supernatant from the cell pellet. With a fine tip transfer pipette, gently remove the remaining liquid near the pellet. Note that the supernatant may probably be cloudy with *E. pacifica* bacteria.

10.4 Resuspended the two cell pellets in a total volume of  25 mL Cell Wash Buffer , combine into one conical tube, and homogenize the cells by vigorously shaking the tube for 30 sec.

Note

Cell Wash Buffer:

420 mM NaCl



50 mM MgCl₂

30 mM Na₂SO₄

10 mM KCl

titrate to pH 8.0 with ~2.4 mM NaHCO₃




sterile filtered through 0.22 µm filter


10.5 For a second time, centrifuge the resuspended cells at  2400 x g, 4°C, 00:03:00 in a swinging bucket rotor. 

10.6 Remove the supernatant as before ().

10.7 Resuspend the cell pellet in  100 µL of Cell Wash Buffer . This resuspension is called the "washed cells."

11 **Prepare**  100 µL **aliquots of** [M] 5*10⁷ Mass Percent .

11.1 Dilute  2 μL of "washed cells" ( go to step #10.7) into  196 μL Artificial Seawater .

11.2 Fix the diluted cells with  2 μL of [M] 37.5 Mass / % volume Formaldehyde and homogenize by vortexing.

11.3 Pipette the fixed cells into a fixed chamber slide and determine the cell concentration.

Note


Remember that concentration of diluted and fixed cells is a 100-fold dilution from the "washed cells." Be sure to factor that dilution into your concentration.

Note

Cells can be counted on a hemacytometer (Neubauer with brightlines) or with an automated cell counter. We recommend a Luna-FL automated cell counter.

Equipment
























LUNA-FL	NAME
Dual Fluorescence Cell Counter	TYPE
Logos Biosystems	BRAND
L20001	SKU
https://logosbio.com/automated-cell-counters/fluorescence/luna-fl	LINK

11.4 After determining the cell concentration, dilute the "washed cells" with Cell Wash Buffer to final concentration of [M] 5×10^7 Mass Percent and split into  100 μL aliquots.

**Note**

One aliquot provides enough cells for 12 nucleofections.

Prime Cells**12 Prime cells for nucleofection by degrading the glycocalyx that surrounds *S. rosetta*.**




- 12.1 **Spin the**  100 μL **aliquots of washed cells** () at  800 x g and  22 $^{\circ}\text{C}$ for  00:02:00 .
- 12.2 Gently remove the supernatant from the cell pellet with a gel-loading pipette tip.
- 12.3 Combine the priming buffer components () to make a final priming buffer ()
 40 millimolar (mM) HEPES-KOH  7.5 ,  34 millimolar (mM) Lithium Citrate ,
 15 Mass / % volume PEG 8000 ,  50 millimolar (mM) Cysteine , and
 1.5 micromolar (μM) Papain)
- 12.4 Resuspend each cell pellet in  100 μL of priming buffer.
- 12.5 Incubate cells for  00:45:00 at  Room temperature . 45m 
- 12.6 Add  10 μL of  50 Mass Percent Bovine Serum Albumin to each aliquot of primed cells for quenching proteolysis from the priming buffer.
- 12.7 Centrifuge cells at  1200 x g, 22 $^{\circ}\text{C}$,  00:04:00 .
- 12.8 Gently remove the supernatant from the cell pellet with a gel-loading pipette tip.
- 12.9 Resuspended each cell pell in  25 μL of SF Buffer ( go to step #8). This suspension of cells is called the "primed cells."




- 12.10 Store the "primed cells" on ice while preparing nucleofection reactions.



Transfect Cells



13 Deliver cargo via nucleofection.

- 13.1 For one transfection, add  16 μL of ice-cold SF Buffer to the Plasmid Delivery Mix () or the *SpCas9* RNP ().

Note


For reactions that use two different gRNAs, assemble each *SpCas9* RNP separately then combine each *SpCas9* RNP at this step. After the *SpCas9* RNPs have been combined, add  16 μL of ice-cold SF Buffer

- 13.2 For genome editing, add  2 μL of the repair oligonucleotide template to the *SpCas9* RNP and SF Buffer ( [go to step #6](#)).

- 13.3 Add  2 μL of "primed cells" (from  [go to step #12.10](#)) to the tube with delivery cargo and SF Buffer. This solution, which is called the "nucleofection mix."

- 13.4 Transfer the entire nucleofection mix into one well of a 96-well nucleofection plate.

Note

At this point, prepare for the recovery step, by transferring the recovery buffer into a convenient vessel and setting the pipette to  100 μL .

- 13.5 Pulse the nucleofection plate with either the **CU 154 (harsher)** or the **CT 151 (milder)** pulse.



Equipment

4D-Nucleofector Core Unit

Control system for performing nucleofection

Lonza

AAF-1002B

https://bioscience.lonza.com/lonza_bs/US/en/Transfection/p/000000000000203684/4D-Nucleofector-Core-Unit

Equipment

96-well Shuttle Device

Add-on for Nucelofector 4d device to perform plate-based nucleofections

Lonza






AAM-1001S

https://bioscience.lonza.com/lonza_bs/US/en/Transfection/p/000000000000191639/96-well-Shuttle-Device


Rest and Recover Cells

- 14 **Allow membranes to reseal by resting cells in recovery buffer before growing cells again in media.**




- 14.1 Immediately after transfection, add  100 μL of ice-cold recovery buffer ( 10 millimolar (mM) HEPES-KOH  7.5 ,  900 millimolar (mM) Sorbitol , and  8 Mass / % volume PEG 8000) to each transfection and gently mixed by firmly tapping the side of the plate.


Note

The recovery buffer should be made ahead of time, sterile filtered through a 0.22 μm membrane, and stored at  4 $^{\circ}\text{C}$.

- 14.2 Allow cells to rest in recovery buffer for  00:05:00 .



- 14.3 Gently mix the well in the nucleofection plate by pipetting up and down before transferring the entire volume in nucleofection well (the nucleofection mix plus the recovery buffer) into  2 mL 10% Red Algae Media in one well of a 6- well cell culture plate.

Note

Cells may also be transferred into  1 mL 10% Red Algae Media in one well of a 12- well cell culture plate.

- 14.4 Incubate at  27 $^{\circ}\text{C}$ for  00:30:00

15 **Add *E. pacifica* food and grow transfected cells.**

- 15.1 Add  10 μL of  10 Mass Percent of *E. pacifica* to the wells in the 6-well plate. For a 12-well plate, only add 5 μL .



Note

Prepare the *E. pacifica* solution by resuspending a frozen, 10 mg pellet of *E. pacifica* in 1 ml of artificial seawater.

15.2 Incubate the cell culture plate at  27 °C for downstream experiments.

