

Sep 20, 2019

Transformation

DOI

<https://dx.doi.org/10.17504/protocols.io.7idhka6>

宏亮¹

¹Northeast Forest University

2019 iGEM NEFU_China

Tech. support email: shengyiyanwork@gmail.com



宏亮 董

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7idhka6>

Protocol Citation: 宏亮 2019. Transformation. protocols.io <https://dx.doi.org/10.17504/protocols.io.7idhka6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: September 20, 2019

Last Modified: September 20, 2019

Protocol Integer ID: 27941

Keywords: Transformation, transformation

Materials

DH5α chemically competent

LB medium

Safety warnings

 Please wear gloves for the experiment.



- 1 Put DH5 α chemically competent bacterium on ice for melting about 15min.
- 2 Add 50 μ l DH5 α chemically competent bacterium into 1.5 ml tubes, mix with ligation product or digested product thoroughly, incubate the bacterium on ice for 30 min.
- 3 Heat shock the bacterium for exactly 45 s at 42 °C water bath.
- 4 Place back on ice for 2 min.
- 5 Add 500 μ l LB medium to each tube, shake the tubes at 37 °C, 200 rpm for 45 min.
- 6 Centrifuge the bacterium at 3000r for 3 min.
- 7 Discard 300 μ l supernatant and resuspend the sediment.
- 8 Plate 100 μ l for a reaction product.
- 9 Incubate plates upside down overnight at 37°C.