

Oct 18, 2019

Transformation and Preparation of Chemically Competent *Bacillus subtilis* cells

DOI

dx.doi.org/10.17504/protocols.io.8fvhtn6

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8fvhtn6>

Protocol Citation: Igem Dusseldorf 2019. Transformation and Preparation of Chemically Competent *Bacillus subtilis* cells. protocols.io <https://dx.doi.org/10.17504/protocols.io.8fvhtn6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 18, 2019

Last Modified: October 18, 2019

Protocol Integer ID: 28885

Keywords: preparation of chemically competent bacillus subtilis cell, chemically competent bacillus subtilis cell, cell

Troubleshooting



Preparation of Chemically Competent *Bacillus subtilis* cells

1

For 20 ml or 50 ml PARIS-Medium mix:

final concentration	Stock	for 50 ml	for 20 ml
100 mM potassium phosphate buffer pH 7.0			
60 mM K ₂ HPO ₄	1 M (3.4 g in 20 ml for stock)	3 ml	1.2 ml
40 mM KH ₂ PO ₄	0.5 M (1.36 g in 20 ml for stock)	4 ml	1.6 ml
3 mM trisodium citrate	0.5 M (0.735 g in 5 ml)	300 μ l	120 μ l
20 mM Potassium-L-glutamate	1 M (4.06 g in 20 ml)	1 ml	400 μ l
21 mM MgSO ₄	1 M (1.23 g in 5 ml)	1050 μ l (7x)	420 μ l (7x)
1 % Glucose	50 % (10g in 20 ml)	1 ml	400 μ l



20 mg/ml L-Tryptophan	5 mg/ml (25 mg in 5 ml)	200 μ l	80
0.1 % Caseinhydrolysate (DIFCO!)	10 % (1 g in 10 ml)	500 μ l	200 μ l

- 2 Inoculation of precultures: Spread *B. subtilis* from Glycerol stock on LB plate and incubate over night for 37°C or over the weekend at 30°C.
- 3 Resuspend single clones in 5 ml Paris Medium (test tube) on incubator roller or plate washer with 1 ml Paris Medium (saves one day work)
- 4 Inoculation of 2.5 ml Paris Medium (test tube) from pre-culture (wash or liquid-overnight culture) to OD₅₈₀=0.2
- 5 4h incubation on incubator roller at 37°C
- 6 Centrifuge 1 ml pellet cells in reaction tubes for 1 min at maximum speed. Remove the supernatant.
- 7 Resuspend the pellet in 1 ml Paris medium with 10 % (v/v) glycerol.
- 8 Store 100 μ l aliquots at -80 °C

Transformation of *Bacillus subtilis*

- 9 Thaw aliquots at 37 °C, add 900 μ l Paris medium and 500-1000 ng plasmid DNA (test tube)
- 10 Incubate 6 h in the incubator roller at 37°C



- 11 Pellet cells, remove 800 μ l, resuspend and plate the rest (for normal transformation on LB+antibiotic)