

Oct 16, 2019

## Transformation of Chemically Competent Cells

DOI

[dx.doi.org/10.17504/protocols.io.8achsaw](https://dx.doi.org/10.17504/protocols.io.8achsaw)

NUS iGEM<sup>1</sup>

<sup>1</sup>National University of Singapore



**NUS iGEM**

National University of Singapore

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.8achsaw>

**Protocol Citation:** NUS iGEM 2019. Transformation of Chemically Competent Cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8achsaw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** October 15, 2019

**Last Modified:** September 11, 2023

**Protocol Integer ID:** 28708

**Keywords:** cell

## Troubleshooting



1 Add  1  $\mu\text{L}$  of miniprep plasmid or  5  $\mu\text{L}$  of Gibson product into competent cell

### Protocol

NAME

**Hi-Fi Gibson Assembly**

CREATED BY

NUS iGEM

Preview

1.1 Calculate for volumes of respective fragments to assemble based on their length and concentration

1.2 Add calculated volume of each fragment (maximum volume:  5  $\mu\text{L}$  )

1.3 Add in  5  $\mu\text{L}$  of 2x Hi-Fi DNA Assembly Master Mix into the PCR tube



NEBuilder HiFi DNA Assembly Master Mix - 10 rxns **New England**  
Biolabs Catalog #E2621S


1.4 Vortex to mix

1.5 Spin down PCR tube

1.6 Incubate samples at  50 °C for  00:45:00

2 Incubate cells in  42 °C water bath for  00:00:45



3 Incubate cells on ice for  00:02:00

4 Add  1 mL of LB media into cell sample


### Protocol

NAME

## Preparation of LB Media


CREATED BY  
NUS iGEM

Preview

4.1 Weigh  25 g of Luria Broth Base powder.



Luria Broth Base (Miller's LB Broth Base)&trade;, powder **Thermo**  
**Fisher Catalog #12795027**

4.2 Add the powder into  1 L of water.





Water refers to sterilized deionized water

4.3 Autoclave entire bottle of LB media.

5 Incubate cells at  37 °C for  01:00:00

6 Spin down the cells at  6000 x g for  00:01:00



- 7 Remove  700  $\mu\text{L}$  of the supernatant
- 8 Re-suspend pellet in remaining media
- 9 Transfer  100  $\mu\text{L}$  of the culture onto the agar plate containing appropriate antibiotic

#### Protocol

NAME


### Preparation of LB Agar

CREATED BY

NUS iGEM

Preview

- 9.1 Weigh  32 g of LB agar powder.

 LB Agar, powder (Lennox L agar) **Thermo Fisher Catalog #22700025**

- 9.2 Add the powder into  1 L of deionized water.

 Water refers to sterilized deionized water

- 9.3 Autoclave entire bottle of LB media.

- 10 Spread the cells evenly

- 11 Incubate at  37  $^{\circ}\text{C}$  overnight

