



Dec 01, 2025

Transformation of Tn7 insertion elements across strains of *Vibrio fischeri*

 [PLOS One](#)

 Peer-reviewed method

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg31k11l25/v1>

Andrew G. Cecere^{1,2}, Chris A. Muriel-Mundo^{1,2}, Derek J. Fisher^{3,4}, Tim I. Miyashiro^{1,2,5}

¹Department of Biochemistry and Molecular Biology, Pennsylvania State University, University Park, PA, USA;

²The One Health Microbiome Center, Huck Institutes of the Life Sciences, Pennsylvania State University, University Park, PA, USA;

³Multidisciplinary Biomedical & Biological Sciences, Southern Illinois University Carbondale, Carbondale, IL, USA;

⁴School of Biological Sciences, Southern Illinois University Carbondale, Carbondale, IL, USA;

⁵Correspondence: tim14@psu.edu

PLOS ONE Lab Protocols

Tech. support email: plosone@plos.org



Tim I Miyashiro

Pennsylvania State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg31k11l25/v1>

External link: <https://doi.org/10.1371/journal.pone.0338632>

Protocol Citation: Andrew G. Cecere, Chris A. Muriel-Mundo, Derek J. Fisher, Tim I. Miyashiro 2025. Transformation of Tn7 insertion elements across strains of *Vibrio fischeri*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg31k11l25/v1>

Manuscript citation:

Cecere AG, Muriel-Mundo C, Fisher DJ, Miyashiro TI (2025) Transformation of Tn7 insertion elements across strains of *Vibrio fischeri*. PLOS One 20(12). doi: [10.1371/journal.pone.0338632](https://doi.org/10.1371/journal.pone.0338632)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 08, 2025

Last Modified: December 01, 2025

Protocol Integer ID: 226768

Keywords: strains of *vibrio fischeri*, genetic content at the tn7 insertion, transformation of tn7 insertion element, tn7 insertion, tn7 insertion element, *vibrio fischeri*, strain, strain es114, dependent natural transformation, genetic content

Funders Acknowledgements:

National Institute of General Medical Sciences

Grant ID: R35 GM152259

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

This protocol details how to use pLosTfoX-dependent natural transformation to transform non-canonical strains of *Vibrio fischeri* with genetic content at the Tn7 insertion of strain ES114.

Materials

Materials needed:

1. Culture tubes, e.g., 14-mL Round-Bottom Polystyrene Test Tubes (Falcon)
2. Defined minimal medium (DMM) with N-acetylglucosamine (GlcNAc) – [50 mM MgSO₄, 10 mM CaCl₂, 300 mM NaCl, 10 mM KCl, 0.01 mM FeSO₄, 0.33 mM K₂HPO₄, 50 mM Tris-HCl (pH 7.5), 0.2% GlcNAc]. Store at 4°C, use within 24 hours.
3. Chloramphenicol (cam)
4. 125-mL baffled Erlenmeyer flask
5. LBS medium [1% (w/v) tryptone, 0.5% (w/v) yeast extract, 2% (w/v) NaCl, 50 mM Tris-HCl (pH 7.5)], with 1.5% w/v agar for solid medium
6. Shaking incubator at 28°C (New Brunswick)
7. Spectrophotometer and cuvettes (Eppendorf)
8. Genomic DNA (extracted according to manufacturer's instructions, Biosearch Technologies Inc.)
9. LBS + 5.0 µg/mL erythromycin (erm) plates

Preparation of *Vibrio fischeri* Cultures

- 1 Initiate a starter culture by inoculating a Falcon tube containing 2 mL DMM-GlcNAc + 2.5 µg/mL cam with an isolated colony of a *V. fischeri* strain harboring pLosTfoX. Incubate starter culture overnight (~16 h) at 28°C shaking at 200 rpm.
- 0.2 Transfer 20 mL DMM + 2.5 µg/mL cam into a 125-mL baffled Erlenmeyer flask and incubate at 28°C to prewarm media.
- 0.3 Measure the OD₆₀₀ of the starter culture. (Note: typically OD₆₀₀ = 1.0–1.5)
- 0.4 Inoculate flask with 1 mL starter culture. Incubate at 28°C shaking at 200 rpm.

Preparation of natural transformation reactions

- 0.5 When the turbidity of culture is OD₆₀₀ = 0.25, transfer 0.5 mL of the culture to two Falcon tubes. To one tube (experimental), add 0.5 µg genomic DNA substrate. To the other tube (no DNA control), add equivalent volume of vehicle. Gently vortex both tubes and incubate at room temperature statically overnight.

Selection of transformants

- 0.6 Add 0.5 mL LBS to each culture tube. Incubate shaking at 28°C.
- 0.7 After 90 minutes, plate 0.2 mL onto LBS-erm. To assess total viable cell count, perform 10-fold serial dilutions to 10⁻⁷ using 10-µL volumes into 90 µL LBS. Spot 10-µL volumes of 10⁻²–10⁻⁷ onto LBS plate with pre-drawn grids. Incubate plates at 28°C overnight.
- 0.8 For each sample, calculate the concentration of viable cells and transformants per mL. Calculate the transformation efficiency by dividing concentration of transformants by the concentration of DNA.
- 0.9 Streak purify ermR colonies onto LBS-erm for further isolation.