

Oct 27, 2025

## Transformation Protocol ATMT Trichoderma atroviride V2

DOI

<https://dx.doi.org/10.17504/protocols.io.261gekk87g47/v1>

José anuel Villalobos Escobedo<sup>1</sup>

<sup>1</sup>Tecnologico de Monterrey

José anuel Villalobos Escobedo: First author

ATMT



José Manuel Villalobos-Escobedo

Instituto Tecnológico y de Estudios Superiores de Monterrey ...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.261gekk87g47/v1>

**Protocol Citation:** José anuel Villalobos Escobedo 2025. Transformation Protocol ATMT Trichoderma atroviride V2. protocols.io <https://dx.doi.org/10.17504/protocols.io.261gekk87g47/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 26, 2025

**Last Modified:** October 27, 2025

**Protocol Integer ID:** 230778

**Keywords:** transformation protocol atmt trichoderma atroviride v2, trichoderma atroviride, agrobacterium tumefaciens, protocol for the transformation, conidia, barseq technique, transformation

**Funders Acknowledgements:**

U.S. Department of Energy

Grant ID: DE-AC02-05CH11231

## Abstract

This is the protocol for the transformation of conidia by *Agrobacterium tumefaciens* to perform the BarSeq technique in *Trichoderma atroviride*, which is reported by Villalobos-Escobedo et al., 2025.

This is the protocol for the transformation of conidia by *Agrobacterium tumefaciens* to perform the BarSeq technique in *Trichoderma atroviride*, which is reported by Villalobos-Escobedo et al., 2025.

2w 4d

## 1 DAY 1

### 1.1

Inoculate 3 plates of *Trichoderma* @ Temperature 27 °C Grow them for 4 days

Sterilize flasks and bottles

5 bottles of 500 ml

2 100 mL

2 250 ml

### 1.2 Prepare LB + Kan (500 mL)

50 ug/mL is the working concentration for Kanamycin; the amount added depends on the stock concentration

Ex: 50 mg/mL stock concentration

1 L of LB would require 1 mL of 50 mg/ml Kan

-Usually, we put 1 ml in 500 ml of LB

## 2 DAY 3

### 2.1 Make liquid agro induction media (ABI)

For 100 transformations, 5 mL of ABI is needed, water is autoclaved, and all others are filtered and sterilized

200 mL to wash, 200 mL to resuspend (grow)

MIX AND FILTER INSIDE HOOD

-- CAS amino acids (10% CAS amino acids) (Heat up and stir)

-- Glucose (36% Glucose) (Heat up and stir)

-- Thiamine (10% Thiamine)



## Sterilized bottles

For AB Buffer, check pH, then filter. The pH of the AB needs to be 5.5.

|  | A                             | B   | C          |
|--|-------------------------------|-----|------------|
|  | ATMT.                         |     |            |
|  | ABI plates(25ml x# of plates) |     |            |
|  |                               | 1L  | liquid ABI |
|  | ddH <sub>2</sub> O(ml)        | 437 | 0.2        |
|  | Bacto Agar(g)                 | 20  | 87.4       |
|  |                               |     |            |
|  | AB salts (ml)                 | 50  | 10         |
|  | 500mM phosphate solution (ml) | 2.4 | 0.48       |
|  | 2x AB buffer (ml)             | 500 | 100        |
|  | 10% thiamine (ul)             | 100 | 20         |
|  | 36% glucose stock (ml)        | 5.5 | 1.1        |
|  | 10% CASamino acids (ml)       | 5   | 1          |
|  | acetosyring one (ml)          | 1   | 0.2        |

## Preparation of buffers (ABI)

- 2.2 Prepare: 1 mL 3,5-dimethoxy-4-hydroxy-acetophenone (acetosyringone) dissolved in dimethyl sulfoxide (DMSO). (make this fresh every time)



-- Acetosyringone is filtered through syringe and 0.22  $\mu$ m filter top in DMSO (0.098 g/mL)

- 2.3
1. Autoclave 437ml ddH<sub>2</sub>O with 20g Bacto Agar
  2. Add MIXED induction liquid to agar.
  3. Mix together inside hood for a few minutes
  4. Pour plates

Pour solid agro induction media plates and leave on the bench to dry for 2 days.

- 2.4 To prepare medium for Agrobacterium LB + Kan.

- 2.5 Inoculate *A. tumefaciens* strain harboring the plasmid for Trichoderma transformation onto LB + Kan.

250  $\mu$ L in 50 mL

- 2.6 Shake at 250 rpm at 30C overnight  
(Place around 7 pm- 8 pm)

- 3 Day 4

- 3.1 Measure the OD<sub>600</sub> of the agro cell culture. (It will likely be around OD<sub>600</sub> 1.)  
(Around 9 am )

- 3.2 To obtain an accurate measurement, dilute the sample 1:10 in H<sub>2</sub>O in the cuvette.  
(100  $\mu$ L agro cell culture, 900  $\mu$ L H<sub>2</sub>O)

- 3.3 Back dilute the agro culture to OD<sub>600</sub> 0.5 in LB+Kan with twice the number of ml of media to the number of transformations you plan to do (i.e., if you are doing 100 transformations, back dilute into 200ml LB).

100 mL from stock and 100 mL new LB+Kan

- 3.4 Shake at 250rpm at 30C for 2h

- 3.5 Measure the OD<sub>600</sub> of the agro cell culture after 2 hrs



3.6 Centrifuge the agrobacterium culture at 4000 rpm, 00:15:00



3.7 Pour supernatant into waste container.

3.8 Wash cells in the Agrobacterium induction media.

-- The number of mL should correspond to twice the number of transformations that you are doing

ex: 20 ml for 10 transformations

New flasks are used

3.9 Centrifuge the agro cells at 4000 rpm, 00:15:00



3.10 Pour supernatant into waste container.

3.11 Resuspend cells in *Agrobacterium* induction media:

The number of mL should correspond to twice the number of transformations that you are doing

For 100 transformations, 200 ml of induction media, and each transformation has  $4 \times 10^7$  conidia.

It's  $2 \times 10^7$  per ml—concentration of *Agrobacterium*.

**We can use tubes with caps for roller.**

3.12 Put cultures on shaker @ 25 °C and rpm 250, for 24 hrs

To prepare materials for the next day:

- Filter top

- 0.45 um Filter

- Vacuum

4 DAY 5

6d



4.1 (using plates of *T. atroviride* that are 4 days old)

NOTE: Make sure the big square plates are at room temperature.

Approximately 23h after putting agro cells into induction media, collect *Trichoderma* conidia from plates into 5-10ml of H<sub>2</sub>O.

4.2 Shake and vortex vigorously to suspend conidia.

4.3 Measure conidial concentration

Count four corner squares, and middle square

Avg spores: Sum of both sides / 10 squares

Math: (Avg. spores)(1×10<sup>4</sup>)(25)(100)

4.4 Aliquot 2e<sup>7</sup> conidia into 2ml Eppendorf tubes

4.5 Spin cells at 10000 rpm, 00:05:00  
or 4000 rpm for 20 min




4.6 Pipet off supernatant

5m



4.7 Incubate at room temperature for at least 5 min

4.8 Put a 0.2 um sterile PES bottle-top filter on a bottle and set up the vacuum.



Transformation in the hood.

Equipment:

Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane

NAME

0.2 µm Bottle Top Filter or equivalent

TYPE

Nalgene

BRAND

09-741-07

SKU

<https://www.fishersci.com/shop/products/nalgene-rapid-flow-sterile-disposable-bottle-top-filters-pes-membrane/0974107>

LINK


500 mL, 0.2 µm PES Bottle Top Filter

SPECIFICATIONS



%% Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane"  
title="Bottletopfilter.Jpg">

- 4.9 Place a 0.45  $\mu\text{m}$  HAWP filter onto the bottle top filter using sterile forceps:



Equipment  
NAME



HAWP MF-Millipore Membrane Filter, 0.45 µm pore size

TYPE

Membrane filter

BRAND

Millipore

SKU

HAWP04700

LINK


[https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?](https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?lang=en&region=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMuPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD_BwE)

[lang=en&region=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMuPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD\\_BwE](https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?lang=en&region=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMuPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD_BwE)

SPECIFICATIONS

0.45 um 47 mm

- 4.10 Pipet the 2ml of Trichoderma/agrobacterium cell mixture onto 0.45um HAWP filter  
  
-- Being careful to only get cell mixture on the HAWP filter.
- 4.11 Wait for the liquid induction media (ABI) to filter through the bottle top filter
- 4.12 Using sterile forceps, carefully transfer the filter with cells on it to the solid induction media petri dish, making sure to keep the side with the cells on it facing up!
- 4.13 Incubate induction media plates with filters at 21 °C in the dark for 6 days.



Dried filters ready for incubation

- 4.14 Materials for next part: - 50 mL or 500 mL tubes with sterile beads and sterile water. - Layered PDA plates. - MM plates for quantifying. - Sterile water.

Arrow

Materials for next part:

- 50 mL or 500 mL tubes with sterile beads and sterile water





- Layered PDA plates
- MM plates for quantifying
- Sterile water

5 Day 11

5.1 Label 2 Falcon tubes for each HAWP filter  
-- 25 - 50 filters fit into 500 mL tubes

5.2 Add 4ml sterile distilled H<sub>2</sub>O to a 50ml conical tube  
-- 500 mL conical tube

5.3 Using sterile forceps, carefully transfer the HAWP filter with cells on it to the 50 mL conical  
-- Use beads and VORTEX until it has all fallen off

5.4 Vortex vigorously to get the rest of the cells off of the HAWP filter

5.5 Transfer 2 mL of cell suspension to a sterile 2 mL Eppendorf tube.

5.6 Spin at 10000 rpm, 00:05:00 .



5.7 Pipet off supernatant into a waste bottle.



5.8 Transfer the rest of the cell suspension to the 2 mL Eppendorf tube and Spin at 10000 rpm, 00:05:00 .



5.9 Pipet off supernatant into a waste bottle



5.10 Resuspend cell mixture into 100  $\mu$ L sterile distilled H<sub>2</sub>O

5.11 To spread in large plates (150 mm) PDA.



- 1.- PDA with Hygromycin B and Carbenicillin (4 ml Hygromycin B and 1 ml Carbenicillin per L): dispense 50 ml.
- 2.- PDA with Carbenicillin (1 ml per Liter): dispense 13.5 ml.
- 3.- Conidia transformants: 1 ml
- 4.- Low-melt PDA with Carbenicillin (1 ml per L). 13.5 ml. 37 C. We can mix with conidia.
- 5.- PDA layer Carbenicillin and Hygromycin B. 13.5 ml, kept at 55 C.

5.12 1 plate with filter goes to 3 plates of top/bottom media

5.13 Plate cells on  
5mL Top: MM + Carb  
25 mL bottom: MM + Carb + Hyg

Note: You don't need to dilute the control plate, but for the experimental ones:

1. 100%:
2. 1:25 **4 ul of cells in 96 ul of water**

Note: Make up media with half the amount of water, autoclave the water separately, then combine.

5.14 Place at 25C in the light for 7 days (checking efficiency with Triton)

## Final Harvest (Day 18)

6



Washing of plates with water

Materials needed:

- 50% glycerol
- Sterile water
- Glass spreaders
- 1 ml Pipettes
- 25 ml pipettors and pipet tip
- Tips 100-1000 ul
- 15 mL Falcon tubes

Dispense 10 mL of sterile water onto the plate

- 7 Use the sterile glass spreader to detach as many as possible.
- 8 Transfer into a 50 mL Falcon tubes.
- 8.1 Add another 10 mL of sterile water to plate and transfer into tube.



9 Spin down 50 mL Falcon tubes at 4,000 RPM for 20 minutes

20m



10 Resuspend pellets in 10 ml of 50% glycerol

10.1 Store in 1.5 mL tubes that are no more than half full, stored at an angle to avoid rupturing the tube.

--Store to -80 °C in Cryotubes.

## Acknowledgements

This material by m-CAFEs Microbial Community Analysis & Functional Evaluation in Soils (m-CAFEs@lbl.gov), a Science Focus Area led by Lawrence Berkeley National Laboratory, is based upon work supported by the U.S. Department of Energy, Office of Science, and Office of Biological & Environmental Research under contract number DE-AC02-05CH11231. This work was supported by a University of California Berkeley Innovative Genomics Institute grant to N.L.G. and a grant from the National Institute of General Medical Sciences of the National Institutes of Health under award number R35GM150926 to L.B.H. We thank Alfredo Herrera-Estrella (Center for Research and Advanced Studies, Irapuato, Guanajuato, Mexico) for his help in designing the conidia experiments. The authors also wish to thank María Belén Mercado-Esquivias for her support in writing the protocol in Protocols.io and for her assistance with the ATMT experiments.