

Sep 28, 2021 Version 2

# Transformation Protocol for BL21(DE3) Competent Cells (C2527I) V.2

Forked from a private protocol



DOI

[dx.doi.org/10.17504/protocols.io.bgtxjwpn](https://dx.doi.org/10.17504/protocols.io.bgtxjwpn)

New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)



New England Biolabs

New England Biolabs

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bgtxjwpn](https://dx.doi.org/10.17504/protocols.io.bgtxjwpn)

External link: <https://www.neb.com/protocols/0001/01/01/transformation-protocol-for-bl21-de3-competent-cells-c2527>

**Protocol Citation:** New England Biolabs 2021. Transformation Protocol for BL21(DE3) Competent Cells (C2527I). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bgtxjwpn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** May 25, 2020

**Last Modified:** September 28, 2021

**Protocol Integer ID:** 37463

**Keywords:** T7 Expression Strain, comp cells, transforming for BL21(DE3)

## Abstract

This transformation protocol is for the C2527I cells. (For the C2527H protocol, see [here](#).)

## Guidelines

### Transformation Protocol Variables

**Thawing:** Cells are best thawed on ice and DNA added as soon as the last bit of ice in the tube disappears. Cells can also be thawed by hand, but warming above 0°C will decrease the transformation efficiency.

**Incubation of DNA with Cells on Ice:** For maximum transformation efficiency, cells and DNA should be incubated together on ice for 30 minutes. Expect a 2-fold loss in transformation efficiency for every 10 minutes you shorten this step.

**Heat Shock:** Both the temperature and the timing of the heat shock step are important and specific to the transformation volume and vessel. Using the transformation tube provided, 10 seconds at 42°C is optimal.

**Outgrowth:** Outgrowth at 37°C for 1 hour is best for cell recovery and for expression of antibiotic resistance. Expect a 2-fold loss in transformation efficiency for every 15 minutes you shorten this step. SOC gives 2-fold higher transformation efficiency than LB medium; and incubation with shaking or rotating the tube gives 2-fold higher transformation efficiency than incubation without shaking.

**Plating:** Selection plates can be used warm or cold, wet or dry without significantly affecting the transformation efficiency. However, warm, dry plates are easier to spread and allow for the most rapid colony formation.

Chemically competent *E. coli* cells suitable for transformation and protein expression.

### Highlights

- Transformation efficiency:  $1\text{--}5 \times 10^7$  cfu/μg pUC19 DNA
- T7 Expression Strain
- Deficient in proteases Lon and OmpT
- Resistant to phage T1 (*fhuA2*)
- B Strain
- Free of animal products


### Genotype

*fhuA2 [lon] ompT gal (λ DE3) [dcm] ΔhsdS*


*λ DE3* = *λ sBamHI ΔEcoRI-B int::(lacI::PlacUV5::T7 gene1) i21 Δnin5*

## Materials

### MATERIALS

 BL21(DE3) Competent E.coli - 6×0.2 ml **New England Biolabs Catalog #C2527I**





































## Safety warnings

 Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

## Before start





Perform steps 2–9 in the tube provided.



- 1 Thaw a tube of BL21(DE3) Competent *E. coli* cells  On ice until the last ice crystals disappear. 
- 2 Mix gently and carefully pipette  50 µL cells into a transformation tube  . 
- 3 Add  1 µL –  5 µL containing  1 pg –  100 ng plasmid DNA to the cell mixture. 
- 4 Carefully flick the tube **4–5 times** to mix cells and DNA. **Do not vortex.** 
- 5 Place the mixture  On ice for  00:30:00 . Do not mix. 
- 6 Heat shock at exactly  42 °C for exactly  00:00:10 . Do not mix. 
- 7 Place  On ice for  00:05:00 . Do not mix. 
- 8 Pipette  950 µL of  Room temperature SOC into the mixture. 
- 9 Place at  37 °C for  01:00:00 . Shake vigorously (  250 rpm ) or rotate. 
- 10 Warm selection plates to  37 °C .
- 11 Mix the cells thoroughly by flicking the tube and inverting. 
- 12 Perform several 10-fold serial dilutions in SOC. 
- 13 Spread  50 µL –  100 µL of each dilution onto a selection plate and incubate  Overnight at  37 °C .  



### Note

Alternatively, incubate at  30 °C for  24:00:00 –  36:00:00 or at  25 °C for  48:00:00 .