



Jun 16, 2025

Transformation Protocol for Targeted Homologous Recombination in *Auxenochlorella protothecoides*

[The Plant Cell](#)


DOI

<https://dx.doi.org/10.17504/protocols.io.x54v922mql3e/v1>Marco A Dueñas¹, Yang-Tsun Lin², Jeffrey Moseley², Sabeeha S. Merchant²¹UC Berkeley; ²University of California, Berkeley

Merchant Lab UC Berkel...

**Marco A. Dueñas**

University of California, Berkeley



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.x54v922mql3e/v1>External link: <https://doi.org/10.1093/plcell/koaf259>



Protocol Citation: Marco A Dueñas, Yang-Tsun Lin, Jeffrey Moseley, Sabeeha S. Merchant 2025. Transformation Protocol for Targeted Homologous Recombination in *Auxenochlorella protothecoides*. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v922mqj3e/v1>

Manuscript citation:

Craig RJ, Dueñas MA, Camacho DJ, Gallaher SD, Avendaño-Monsalve MC, Lin Y, Blaby-Haas CE, Moseley JL, Merchant SS (2025) Targeted genetic manipulation and yeast-like evolutionary genomics in the green alga *Auxenochlorella*. *The Plant Cell* 37(11). doi: [10.1093/plcell/koaf259](https://doi.org/10.1093/plcell/koaf259)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 29, 2024

Last Modified: June 16, 2025

Protocol Integer ID: 100870

Keywords: transformation protocol for targeted homologous recombination, *auxenochlorella protothecoide*, targeted homologous recombination, gene cassettes into the genome, targeted gene cassette, lithium acetate transformation method, *auxenochlorella*, transformant strain, transformation protocol, genome

Funders Acknowledgements:

Department of Energy (DOE) Office of Science, Biological and Environmental Research


Grant ID: DE-SC0023027

NIH T32 Genetics Dissection of Cells and Organisms Grant

Grant ID: 1T32GM132022-01

Abstract

Detailed instructions for the introduction of targeted gene cassettes into the genome of *Auxenochlorella Protothecoides* using a lithium acetate transformation method. Following this protocol should yield transformant strains able to be grown on selectable media. This protocol is modified and adapted from Moseley et al. (1) and Dueñas et al. (2).





Materials

Auxenochlorella protothecoides cells

1.5 mL Graduated Flat Cap Tube (Fisherbrand, #02-681-320)

nuclease-free H₂O

High-Fidelity® restriction enzyme (New England Biolabs)

rCutSmart™ Buffer (New England Biolabs, #B6004S)

Purple Gel Loading Dye, 6X (New England Biolabs, #B7024S)

50 mL Polypropylene Conical Tube (Falcon, #352098)

15 mL Polypropylene Conical Tube (Falcon, #05-539-4)

Polyethylene glycol 4,000 (Alfa Aesar, # A16151)

Tris Hydrochloride (Fisher #BP153-500)

Lithium Acetate dihydrate, 98% (Acros Organics, #6109-17-4)

EDTA (Fisher, #S311-500)

phenol-chloroform-isoamyl alcohol (Fisher #BP17521-400)

Isopropanol (Sigma-Aldrich, #W292907)

Tris-HCl (Sigma-Aldrich, #T5941)

Ethanol (Sigma-Aldrich, #E7023)

Glucose (referred to as Dextrose Anhydrous, Fisher #BP350500)


sterile Milli-Q H₂O

Thiamine-HCL (Macron Fine Chemicals™, # 0000238814)

ApM1 Media (Refer to)

Safety warnings



phenol-chloroform-isoamyl alcohol  (Fisher #BP17521-400) is corrosive. Steps involving this chemical should be done in a hood and with proper PPE and glove wear.

Preparation of DNA for transformation

10m

- 1 Before you start, check your designed construct to ensure your gene cassette can be properly digested. Ensure you have restriction sites for digestion at each end of your targeting flanks and no internal restriction sites that will be digested within the cassette.
- 2 Digest  100 µg of plasmid construct in  750 µL volume to linearize DNA (usually overnight). Volume of plasmid solution will be dependent on your initial plasmid concentration. It is recommended to perform a maxiprep of your plasmid to ensure adequate plasmid stock. Use the following formula to determine what concentrations to use for your overnight mix:

Plasmid volume to use (Y) = $100000 / (X \text{ ng/}\mu\text{L of plasmid})$

Overnight digestion mix

Y µL Plasmid solution

 75 µL rCutSmart™ Buffer (New England Biolabs, #B6004S)

 10 µL of each High-Fidelity® restriction enzyme (New England Biolabs)

Z µL of nuclease-free H₂O (For adjustment to a final volume of 750 µL)

750 µL Total


Note

If calculated plasmid volume exceeds total, do not add any water and proceed with the maximum volume of plasmid for the  750 µL reaction. Note that this may decrease transformation efficiency.

- 3 Run a small aliquot on a gel. Make a mix consisting  3 µL of digestion mix,  2 µL of 6X loading dye (New England Biolabs, #B7024S), and  10 µL nuclease-free H₂O. Ensure complete digestion and proper cassette size.

Note

If digestion appears incomplete, allow restriction digest to proceed for additional time before proceeding,



Example of a linearized DNA for transformation. The top bands represents the gene cassette of interest

- 4 Extract digest with  750 μL 25:24:1 phenol-chloroform-isoamyl alcohol  08 (Fisher #BP17521-400) to kill enzyme. Proceed with centrifugation at  10,000 rcf, Room temperature for  00:10:00 . 10m
- 5 Transfer aqueous phase (top layer) to a new  1.5 mL tube (Fisherbrand, #02-681-320). Precipitate DNA with  525 μL isopropanol (Sigma-Aldrich, #W292907). Incubate at room temperature for  01:00:00 or overnight at  4 $^{\circ}\text{C}$. 1h
- 6 Air-dry the DNA with cap open for 30-60 minutes in sterile hood. After this step, a small, white pellet should be visible. Dissolve DNA in  50 μL in nuclease-free H_2O or sterile elution buffer ( 10 millimolar (mM) Tris-HCl (Sigma-Aldrich, #T5941),  7.5). DNA should be between  1-2 $\mu\text{g}/\mu\text{L}$.
- 7 Spin down at max speed in microfuge for  00:30:00 . Discard the supernatant and wash the DNA pellet 2 times with  500 μL  70 % (v/v) Ethanol (Sigma-Aldrich, #E7023) . Do the wash steps in the sterile hood to avoid contamination. Centrifuge at max speed for  00:05:00 between each wash and remove all traces of Ethanol with a pipette. 35m
- 8 Air-dry the DNA with cap open for 30-60 minutes in sterile hood. Dissolve DNA in $\sim 50 \mu\text{L}$ of sterile elution buffer ( 10 millimolar (mM) Tris-HCl,  7.5). NanoDrop to determine final concentration, which should be between 1-2 $\mu\text{g}/\mu\text{L}$.



Note

A DNA concentration of 1-2 µg/µL will be adequate for about 3-5 sets transformations. If final concentrations are lower than this amount, transformation can still proceed but ensure to add a higher amount of DNA in step 15.

Making a Transformation Starter Culture


- 9 Start a culture of wild-type UTEX 250 cells (or A. pro strain of interest) in  50 mL ApM1 media  26-28 °C in a dark incubator, shaking at 200 rpm for  36:00:00 to  48:00:00. This can be adjusted as needed but aim for an OD 750 of about 2.0-3.5 ($\sim 3 \times 10^7$ cells/mL) for the transformation. This can be done from an aliquot as well.

3d 12h

Note

For THIC selection ONLY, the ApM1 start culture should be grown without thiamine added to the media so that the cell culture can deplete residual amounts. To ensure adequate cell concentration and proper thiamine depletion, start the culture with a full loop of cells and allow to grow for 2 days. If this is not done, final transformant cell plates will be overrun and selection will fail.

For information on how to make ApM1 media, plates, glucose stocks, and thiamine stocks, refer to *Auxenochlorella* ApM1/ApRM1 media recipes and stocks V.2 (dx.doi.org/10.17504/protocols.io.q26g71pq8gwz/v2)



Example of a UTEX 250 transformation starter culture grown in ApM1 2% glucose, 1x thiamine. Pictured in the left flask is a culture just after initial inoculation. On the right is the cells culture after 2 days of growth, with a final OD 750 of 3.5.

Preparing Solutions for Transformation

10 Prepare the following stock solutions:

Dissolve 90.02 g of Lithium Acetate dihydrate, 98% (Acros Organics, #6109-17-4) to 1 L MQ water to make a 1 Mass Percent **Lithium Acetate (LiAC) solution**.

Dissolve 1211.4 grams of Tris-HCl (Tris Hydrochloride, Fisher #BP153-500) and 372.24 grams of EDTA (Fisher, #S311-500) in 1L of MQ water to make a **1000x TE Buffer**.

Dissolve 500 g of Polyethylene glycol 4,000 (PEG-4000, Alfa Aesar, # A16151) to 500 mL MQ water to make a 50% (w/V) **PEG-4000 solution**.

11 Before starting the transformation, make a .1M LiAc/1X TE solution. For a single round of transformation, the following stock is recommended:

LiAc/1xTE Mix (10 mL)

- 1 mL 1 Molarity (M) Lithium Acetate
- 1 mL 1X TE (1m L from 1000x TE stock)
- 8 mL sterile Milli-Q H₂O

**Note**

The lithium acetate solution is critical in increasing permeability of the cell wall for DNA uptake.

- 12 Before starting the transformation, make a .40% PEG/.1M LiAc/1x TE solution. For a single round of transformation, the following stock is recommended:

40% PEG/.1M LiAc/1xTE (5 mL)

-  5 mL  1 Molarity (M) LiAc
-  5 mL 1x TE
-  4 mL  50 Mass / % volume PEG-4000

Note

For more than 5 transformations, a 10 mL 40% PEG/.1M LiAc/1xTE solution can be used.

Day 1 of Lithium Acetate Transformation**1h 45m**

- 13 Transfer the cells to a  50 mL Polypropylene Conical Tube (Falcon, #352098) and centrifuge for at  3750 rpm for  00:10:00 . 10m
- 14 Remove supernatant. Add  5 mL of LiAc/1xTE from step 11 to the cell pellet and mix thoroughly while keeping sterile. Repeat centrifugation at  3750 rpm for  00:10:00 . 10m
- 15 Remove supernatant and add  500 μ L of LiAc/1xTE and mix with a pipette. Transfer this mixture to a  1.7 mL Eppendorf tube (sterile) and allow it to shake for  01:00:00 . 1h

**Note**

Final amount within the tube will most likely exceed $500\ \mu\text{L}$ due to the volume of the cells and residual solution from the washing step. This should not effect the rest of the protocol and can be used for additional aliquots in step 15.

- 16 After $01:00:00$, aliquot $150\ \mu\text{L}$ into eppendorf tubes for each of your plasmid transformations plus a blank. Add DNA cassette to the tube based on the following formula:

$$\text{DNA to add} = 15 / (\text{NanoDrop concentration from step 8} / 1000).$$

This should ideally be within a range of $10\text{--}20\ \mu\text{L}$. If higher, use the suggested amount. If lower, use the suggested amount or add $10\ \mu\text{L}$ for simplicity.

In one aliquot, add $15\ \mu\text{L}$ MQ water. This will serve as a transformation blank.

Note

Including the blank will be CRITICAL to ensure the transformation worked properly. Do not skip this step.

- 17 Incubate in the shaker for $00:30:00$

- 18 Add $750\ \mu\text{L}$ of the 40% (v/v) PEG-400 / 0.1% Mass Percent LiAc / 1xTE from step 12. Incubate the tubes overnight in a shaker at $26\text{--}28\ ^\circ\text{C}$.

Note

The PEG solution is critical for DNA permeability, stability, aggregation, and eventual incorporation into the cells.



Overnight incubation of cells in 750 μ L of 40% PEG-400 / .1 M LiAc / 1xTE solutions.

Day 2 of Lithium Acetate Transformation

8h 0m 30s

- 19 For antibiotic selection, add cells to 10 mL ApM1, 2% glucose, 1X thiamine in a 25 mL flask and recover in the dark for ~6-8 hours at 24-26 °C , ~120 rpm shaking to allow expression of the antibiotic resistance gene. After this allotted time, transfer cells to 1 15 mL Polypropylene Conical Tube (Falcon, #05-539-4), and centrifuge at 3750 rpm for 00:05:00 and remove the supernatant.

5m

Note

This step is critical in antibiotic selection, skipping it will severely reduce viable colonies or none will appear entirely. For THIC selection, this step can be skipped entirely.

- 20 If step 19 was skipped, spin down cells in microfuge at 5000 rpm for 00:00:30 and remove supernatant. For all selections, resuspend cells in 250 μ L 1 Mass Percent sorbitol.

30s


**Note**

The sorbitol supports cell recovery as an osmotic stabilizer after membrane stress and aids in the adherence of cells to the surface of the plates.


- 21 Split each transformation between 2 selective plates, ~  180 μL per plate. Spread with sterile 4 mm glass beads (Fisher #11-312B) and dry in hood.
- 22 For heterotrophic growth, incubate in the dark at  26-28 °C Visible colonies should begin to appear usually after 5-7 days depending on the selection marker. If colonies do not appear after 14 days, the transformation may have failed. Redo and ensure all steps have been followed.


Note

Your colony numbers will vary depending on the method of selection and where your gene cassette is targeted. Generally, you should expect about 10-50+ colonies from a successful transformation plate. You should see close to no colonies (only a few escapes) appear for the blank selection plates.



Example of transformant colonies after 7 days post lithium acetate transformation using a THIC selection marker. Colonies were selected on agarose plates ApM1, 2% glucose media lacking any exogenous thiamine. The bottom right plate is a blank where wild type cells were not treated with DNA.





Example of transformant colonies after 7 days post lithium acetate transformation using a *neoR* selection marker. Colonies were selected on agar plates consisting ApM1, 2% glucose, 1xT Thiamine media with 100 ug/mL of the antibiotic G418. The leftmost plate is a blank where wild type cells were not treated with DNA.

Single Colony Purification of Transformant Colonies

- 23 Using a 24 well plate, add 1mL of sterile ApM1 media with 2% glucose. For antibiotic selection, include 1X thiamine and include the antibiotic for further selection. For THIC selection, do not add thiamine



24 Using either a sterile toothpick or pipette tip, pick 6-12 colonies and mix into a well. Repeat this for 6-12 colonies per strain. We also recommend using a well for a negative control (WT cells) and a positive control (If applicable) to ensure selection media is working as intended.

25 Allow well cultures to grow for a period of 3-5 days.

26 Perform dilutions in either a separate 24 well plate or 1.5 mL Eppendorf tubes. Fill each with either sterile MQ water or ApM1 media. For each well plate culture, you will need a set of 2 dilutions

27 For the dilution 1, add  10 μ L of your initial selection culture into the first  1 mL well/tube to make a 1/100 dilution.

Note

Avoid mixing for all dilutions steps is CRITICAL to ensure accurate dilutions. This can be done continuously with a P1000 pipette for well plates or capping and shaking with an eppendorf tube.

28 Take dilution 1 from the previous step and further dilute to make dilution 2 . We recommend a final dilution ranging from $\sim 1/1000$ to $\sim 1/5000$. See the following amounts for the approximate final dilution:

1/1000: Add 100 μ L from dilution 1 into dilution 2

1/2500: Add 40 μ L from dilution 1 into dilution 2

1/5000: Add 20 μ L from dilution 1 into dilution 2

Note

A higher dilution will ensure more spaced but less colonies. Less dilution will ensure more colonies but tighter spacing for colony picking. If no colonies appear, decrease the dilution amount and repeat single colony purification.

29 For each dilution, add 200 μ L of dilution 2 onto a selection plate and spread using 4 mm glass beads (Fisher #11-312B). Allow to dry for about 30 minutes to an hour, and then wrap and store in the dark in 26-28C. Wait 4-14 days for colonies to appear.



Note

Colony abundance and time of appearance will be dependent on the dilution used. For a lower dilution, expect colonies to appear around 5-10 days post incubation. For a higher dilution, allow up to 14 days for colonies to appear. If no colonies appear after 21 days, redo the single colony purification at a lower dilution.

Protocol references

1. Moseley et al. Production of lipids and terpenoids in *Auxenochlorella protothecoides* 12 (2021). US Patent US12037630B2, filed November 5th, 2021, granted July 16th, 2024.
2. M.A. Dueñas, R.J. Craig, S.D. Gallaher, J.L. Moseley, & S.S. Merchant, Leaky ribosomal scanning enables tunable translation of bicistronic ORFs in green algae, *Proc. Natl. Acad. Sci. U.S.A.* 122 (9) e2417695122, <https://doi.org/10.1073/pnas.2417695122> (2025).

Acknowledgements

This work was supported by the Department of Energy (DOE) Office of Science, Biological and Environmental Research program under award no. DE-SC0023027 (to J.L.M. and S.S.M.). M.A.D. was supported, in part, by the NIH T32 Genetics Dissection of Cells and Organisms Grant (1T32GM132022-01) and the Newton Graduate Fellowship in Synthetic Biology (QB3-Berkeley).