

Sep 22, 2019

Transformation using Electroporation

DOI

dx.doi.org/10.17504/protocols.io.7jzhkp6

Alba Balletbó¹

¹Wageningen University

iGEM Wageningen 2019



Alba Balletbó

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.7jzhkp6

Protocol Citation: Alba Balletbó 2019. Transformation using Electroporation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.7jzhkp6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 22, 2019

Last Modified: September 22, 2019

Protocol Integer ID: 27993

Materials

DNA solution

SOC Medium (See SOC Medium, Protocols)

Antibiotic (If applicable)



- 1 Pre-chill the electroporation cuvette on ice, thaw the electrocompetent cells on ice as well and mix them by flicking gently.
- 2 Transfer the amount of cells required for the electroporation cuvette in question (e.g. 50 μ L). Add 1 μ L of DNA solution.
- 3 Insert the correct specifications (Voltage, Resistance, etc.) for the used transformation strain into the pulse generator. Place the electroporation cuvette and apply the pulse.
- 4 Immediately add up the volume to 1mL with SOC medium (e.g. 950 μ L) and incubate the culture for 1 hour at 37°C.
- 5 Spin the culture down at 3000rpm for 3 minutes, discard 900 μ L and resuspend the pellet in the remaining 100 μ L.
- 6 Plate accordingly on agar plates with the appropriate antibiotic, if applicable.
- 7 Incubate the plates overnight at 37°C.