

Mar 05, 2024

Transforming pJC8 into HB101

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj3kkwlk5/v1>

Bonnie Evans^{1,2}

¹MRC Laboratory of Medical Sciences; ²Imperial College London



Bonnie Evans

MRC Laboratory of Medical Sciences, Imperial College London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj3kkwlk5/v1>

Protocol Citation: Bonnie Evans 2024. Transforming pJC8 into HB101. protocols.io

<https://dx.doi.org/10.17504/protocols.io.n2bvj3kkwlk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 05, 2024



Last Modified: March 05, 2024

Protocol Integer ID: 96163

Keywords: pjc8 into hb101, transforming pjc8, pjc8, hb101, empty cosmid, control strain

Abstract

Taken from *Mix and Go E. coli* Transformation Kit (see attachment).

Transforming *E. coli* HB101 with pJC8-empty cosmid to use as a control strain.

Attachments



t3001_t3002_mix_go_...

870KB

Materials

Mix and Go E. coli Transformation Kit (Zymo Research, T3001/T3002)

Contents of SOB medium:

20g Bacto Tryptone

5g Yeast Extract

0.58g Sodium Chloride

0.19g Potassium Chloride

1L sterile H₂O



Before starting

- 1 Prepare tetracycline LB agar plates

Protocol

NAME

Making tetracycline LB agar plates

CREATED BY

Bonnie Evans

Preview

- 2 Warm plates in 37 °C incubator

Note

Chilled plates will decrease *Mix & Go* cell transformation efficiency.

- 3 Prepare SOC medium

- 3.1 Get SOB medium and 20 % glucose from media kitchen

- 3.2 Add 2 mL 20 % glucose to 100 ml SOB medium

Transformation

- 4 Isolate pJC8 cosmid from *E. coli* DH5α



Protocol

NAME

Qiagen QIAprep Spin Miniprep for cosmids from metagenomic library

CREATED BY

Bonnie Evans

[Preview](#)

- 5 Thaw 100 uL aliquot of *E. coli* HB101 Mix & Go competent cells on ice
- 6 Add 5 uL pJC8 cosmid
- 7 Mix by tapping the tube and shaking downwards once
- 8 Place on ice for 5 minutes
- 9 Add 400 uL SOC medium
- 10 Incubate at 37 °C with shaking at 200 rpm for 1 hour

Note

An outgrowth in SOC medium is required for efficient transformation when selecting with tetracycline.

- 11 Pipette all (505 uL) of the transformation mix onto a tetracycline LB agar plate and spread using plating beads.
- 12 Dry under laminar flow for a few minutes



13 Incubate plate at 37 °C overnight

Expected result

If the transformation is successful, you will have individual colonies.