



Feb 11, 2025

Transmission electron microscopy

 [Nature Cell Biology](#)

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l291opv1y/v1>

Ching-Chieh Chou¹, Judith Frydman¹

¹Stanford University



Ching-Chieh Chou

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l291opv1y/v1>

External link: <https://doi.org/10.1038/s41556-025-01623-y>

Protocol Citation: Ching-Chieh Chou, Judith Frydman 2025. Transmission electron microscopy. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l291opv1y/v1>

**Manuscript citation:**

Chou, CC., Vest, R., Prado, M.A. *et al.* Proteostasis and lysosomal repair deficits in transdifferentiated neurons of Alzheimer's disease. *Nat Cell Biol* **27**, 619–632 (2025). <https://doi.org/10.1038/s41556-025-01623-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 11, 2025

Last Modified: February 11, 2025

Protocol Integer ID: 119919

Keywords: ASAPCRN, transmission electron microscopy processing, transmission electron microscopy, transmission electron microscopy analysis, human cell

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Aligning Science Across Parkinson's

Grant ID: ASAP-024268

Abstract

Processing human cells for transmission electron microscopy analysis.

- 1 Cells are grown on Ibidi dishes: μ -Dish 35 mm, high Grid-50 Glass Bottom is a 35 mm.
- 2 Cells are fixed in Karnovsky's fixative: 2% Glutaraldehyde (EMS Cat# 16000) and 4% PFA (EMS Cat# 15700) in 0.1 M Sodium Cacodylate (EMS Cat# 12300) pH 7.4 for 1 hr at room temperature. After fixation samples are placed at 4°C.
- 3 Cells are then post-fixed in cold 1% Osmium tetroxide (EMS Cat# 19100) in water and allowed to warm for 2 hr in a hood, washed 3X with ultra-filtered water, then en bloc stained 2 hr in 1% Uranyl Acetate at room temperature.
- 4 Samples are dehydrated in a series of ethanol washes for 10 min each at room temperature beginning at 30%, 50%, 70%, 95%, changed to 100% 2X, then Propylene Oxide (PO) for 10 min.
- 5 Samples are infiltrated with EMBED-812 resin (EMS Cat#14120) mixed 1:1, and 2:1 with PO for 2 hr each.
- 6 Samples are placed into EMBED-812 for 2 hr opened then placed into flat molds w/labels and fresh resin and placed into 65°C oven overnight.
- 7 Cells of interest are located using the grid pattern and cut out with a gem saw and remounted on pre-labeled resin blocks with fresh resin and polymerized overnight again.
- 8 Once full polymerized the glass coverslip is etched away using hydrofluoric acid for 20 min.
- 9 Using the finder grid pattern left behind the block faces are trimmed down allowing for serial sectioning of the cells of interest.
- 10 Sections are taken around 90 nm, picked up on formvar/Carbon coated slot Cu grids, stained for 40seconds in 3.5% Uranyl Acetate in 50% Acetone followed by staining in 0.2% Lead Citrate for 6 min.
- 11 Observed in the JEOL JEM-1400 120kV and photos are taken using a Gatan Orius 2k X 2k digital camera.