

Aug 21, 2025

Transmission electron microscopy of fibroblasts

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l6zxxzgqe/v1>

Ashley Seifert^{1,2}, Ajoy Aloysius^{1,2}

¹University of Kentucky; ²Department of Biology



Ashley Seifert

University of Kentucky

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l6zxxzgqe/v1>

Protocol Citation: Ashley Seifert, Ajoy Aloysius 2025. Transmission electron microscopy of fibroblasts. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6zxxzgqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 15, 2025



Last Modified: August 21, 2025

Protocol Integer ID: 224800

Keywords: ASAPCRN, transmission electron microscopy of fibroblast, cultured cells by transmission electron microscopy, transmission electron microscopy, imaging cultured cell, seifert lab, fibroblast

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

NIH

Grant ID: R01AR070313

Abstract

This protocol describes the procedure for preparing and imaging cultured cells by transmission electron microscopy used in the Seifert Lab.

Protocol materials

- ⊗ 4% Paraformaldehyde solution
- ⊗ Glutaraldehyde **Electron Microscopy Sciences Catalog #16220**
- ⊗ Sorensens Phosphate Buffer **Electron Microscopy Sciences Catalog #11600**
- ⊗ Osmium Tetroxide **Electron Microscopy Sciences Catalog #19134**
- ⊗ 100% Ethanol
- ⊗ Uranyl acetate **Electron Microscopy Sciences Catalog #22400**
- ⊗ Lead Citrate (Reynolds) **Delta Microscopies Catalog #11300**



- 1 Fix cells by replacing the growth medium with a fixative solution containing
 4% Paraformaldehyde solution and 3.5% EM-grade
 Glutaraldehyde **Electron Microscopy Sciences Catalog #16220** in 0.1 M
 Sorensens Phosphate Buffer **Electron Microscopy Sciences Catalog #11600**
at  4 °C
- 2 Incubate for  00:05:00 at  Room temperature
- 3 Replace solution with fresh fixative solution
- 4 Incubate for  01:00:00 at  4 °C
- 5 Wash cells thoroughly with 0.1 M
 Sorensens Phosphate Buffer **Electron Microscopy Sciences Catalog #11600**
- 6 Post-fix in 1%  Osmium Tetroxide **Electron Microscopy Sciences Catalog #19134**
for  01:00:00 at  4 °C
- 7 Dehydrate sample in a graded ethanol series from 50% through 100% for  00:05:00
each
- 8 Wash twice with  100% Ethanol at  4 °C
- 9 Infiltrate with epoxy resin directly in the culture dish
- 10 Transfer into beam capsules and embed overnight
- 11 Section polymerized blocks at  80-90 nm using an ultramicrotome



- 12 Mount sections on 300-mesh copper EM grids
- 13 Stain with  Uranyl acetate **Electron Microscopy Sciences Catalog #22400** and  Lead Citrate (Reynolds) **Delta Microscopies Catalog #11300** to enhance contrast
- 14 Image on a transmission electron microscope