



Jun 19, 2023

TSS competent cells and transformation

DOI

dx.doi.org/10.17504/protocols.io.5jyl8jzbrg2w/v1



David Ruano Gallego¹

¹Universidad Autónoma de Madrid



David Ruano Gallego

Universidad Autónoma de Madrid

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8jzbrg2w/v1>

Protocol Citation: David Ruano Gallego 2023. TSS competent cells and transformation . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl8jzbrg2w/v1>

Manuscript citation:



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: April 11, 2023

Last Modified: June 19, 2023

Protocol Integer ID: 80304

Keywords: efficient transformation of plasmid, tss competent cell, plasmid, competent cell, efficient transformation



Abstract

This protocol allows the efficient transformation of plasmid into E. coli and related strains.



Materials

Transformation Storage Solution (TSS) buffer

- PEG3350 10%. Sigma. Ref: 202444
- DMSO 5%.  DMSO (dimethyl sulfoxide) Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8418
- Glycerol 10%
 Glycerol - for molecular biology, ≥99% Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516
- MgCl₂ 20 mM.
 1 M Magnesium Chloride (MgCl₂) Merck MilliporeSigma (Sigma-Aldrich) Catalog #M8266
- LB medium 2x.

Adjust to pH 6.1 with 6 M HCl.

Autoclave.

Add previously sterilized MnCl₂ 4 M to a final concentration of 140 mM.

Add water to adjust the modified LB medium to 1x.

Prepare aliquots of 1 ml and store at -20°C.

5x KCM solution

- KCl 0.5 M.  Potassium Chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9541
- CaCl₂ 150 mM.  Calcium chloride dihydrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #C7902
- MgCl₂ 250 mM.  Magnesium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #M8266

Sterilize at 121°C for 20 min. Store at 4°C.

Troubleshooting

Before start

















Prepare solution TSS and KCM.

When prepared, TSS has a cloudy appearance that disappears once autoclaved. Also, over time, TSS can seem to be contaminated without being so, since the polyethylene glycol tends to precipitate.



Preparing competent cells

2d

- 1 Grow o/n a 5 mL culture of the strain to be transformed.
 5 mL  37 °C  Overnight
- 2 Dilute 1:100 of the preinoculum in 50 mL prewarmed LB ($OD_{600} = 0,05$).
 50 mL  37 °C
- 3 Grow at 37°C with agitation (170 rpm) until it reaches $OD_{600} = 0,5$.
 170 rpm, 37°C
- 4 Centrifuge the culture at 4000 g and 4°C for 10 min.
 4000 x g, 4°C, 00:10:00
- 5 Discard supernatant and wash with 1 mL of chilled TSS.
 1 mL  On ice
- 6 Centrifuge at 4000 g and 4°C for 3 min.
 4000 x g, 4°C, 00:03:00
- 7 Discard supernatant and resuspend in 1 mL of TSS.
 1 mL  On ice
- 8 Incubate 10 min on ice.
 On ice  00:10:00
- 9 Store 30 µL aliquots at -80°C.
 30 µL  -80 °C

12h

1h 30m



10m

3m

10m

Transformation protocol

1h 40m

- 10 Mix 5 µL of 5x KCM solution with DNA (around 1-5 µL) and H₂O to a total of a 25 µL mixture.  25 µL
- 11 Mix the mixture with 25 µL of competent cells.
 25 µL
- 12 Incubate 30 min on ice.

30m



🌡️ On ice ⌚ 00:30:00

13 Thermal shock: incubate at 42°C for 90 s.

1m 30s

🌡️ 42 °C ⌚ 00:01:30

14 Incubate 2 min on ice.

2m

🌡️ On ice ⌚ 00:02:00

15 Recovery: add 250 µl of LB and incubate at 37°C for 1 h.

1h

🌡️ 37 °C ⌚ 01:00:00

16 Plate 100 µl of the solution on LB agar plates supplemented with the indicated antibiotic and incubate o/n at 37°C.

12h

🌡️ 37 °C ⌚ Overnight

Protocol references

<https://pubmed.ncbi.nlm.nih.gov/2648393/>

<https://pubmed.ncbi.nlm.nih.gov/36317996/>