

May 09, 2019

U Cinn - NEFA Concentration

DOI

dx.doi.org/10.17504/protocols.io.xmnfk5e



Patrick Tso¹, Dana Lee¹

¹University of Cincinnati

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xmnfk5e>

External link: <https://mmpc.org/shared/document.aspx?id=196&docType=Protocol>

Protocol Citation: Patrick Tso, Dana Lee 2019. U Cinn - NEFA Concentration. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xmnfk5e>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 31, 2019

Last Modified: May 09, 2019

Protocol Integer ID: 19854

Keywords: non-esterified fatty acids, NEFA-HR enzymatic colorimetric method, nefa concentration summary, hr enzymatic colorimetric method assay, enzymatic colorimetric method assay, nefa, serum, assay

Abstract

Summary:

Quantitative determinations of non-esterified fatty acids in plasma/serum/lymph will be made using the NEFA-HR enzymatic colorimetric method assay.



Materials

MATERIALS

⊗ HR Series NEFA-HR(2) Color Reagent A FUJIFILM Wako Pure Chemical Corporation Catalog #999-34691

⊗ HR Series NEFA-HR(2) Solvent A FUJIFILM Wako Pure Chemical Corporation Catalog #995-34791

⊗ HR Series NEFA-HR(2) Color Reagent B FUJIFILM Wako Pure Chemical Corporation Catalog #991-34891

⊗ HR Series NEFA-HR(2) Solvent B FUJIFILM Wako Pure Chemical Corporation Catalog #993-35191

⊗ NEFA Standard Solution FUJIFILM Wako Pure Chemical Corporation Catalog #276-76491

Reagent Preparation:

Working Color Reagent Solutions A:

Reagents and Materials:

Color Reagent A

Solvent A

Procedure:

Reconstitute **Color Reagent A** with a portion of **Solvent A** and then transfer entire contents into **Solvent A** bottle, rinsing **Color Reagent** vial several times.

Working Color Reagent Solutions B:

Reagents and Materials:

Color Reagent B

Solvent B

Procedure:

Reconstitute **Color Reagent B** with a portion of **Solvent B** and then transfer entire contents into **Solvent B** bottle, rinsing **Color Reagent** vial several times.

Note:

FUJIFILM Wako RRID:SCR_013651

Troubleshooting

- 1 Prepare working Color Reagent Solutions A and B.
 - A. Reconstitute **Color Reagent A** with a portion of **Solvent A** and then transfer entire contents into **Solvent A** bottle, rinsing Color Reagent vial several times.
 - B. Reconstitute **Color Reagent B** with a portion of **Solvent B** and then transfer entire contents into **Solvent B** bottle, rinsing Color Reagent vial several times.
- 2 Locate working Standard (1mmol/L or 1 mEq/L).
THIS ASSAY DOES NOT REQUIRE A SERIAL DILUTION
- 3 Using a 96 well flat bottom plate, into separate wells, pipette 5µL of deionized water, 1mMstandard, or sample to be assayed.
- 4 Add 200µL of **Color Reagent Solution A** to all wells.
- 5 Mix well and Incubate plate for 5 minutes at 37°C.
- 6 Measure the absorbance of each well at 550nm (sub:660nm). This measurement (Abs1) will serve as the sample blank.
- 7 Add 100µL of **Color Reagent Solution B** to all wells.
- 8 Mix well and Incubate plate for 5 minutes at 37°C.
- 9 Measure the absorbance of each well at 550nm (sub:660nm). This will be your Abs2 value.
- 10 Obtain the final absorbance (Sample_{abs}) by subtracting the first reading (step 5) from the second reading (step 8). *
- 11 Plot the absorbance vs. concentration to construct the calibration curve. A linear calculation model should be used.
- 12 To calculate sample concentration by calculation use the following formula:



Sample Conc. = (Sample Absorbance/Standard Absorbance) * Standard Concentration

*The sample blank absorbance (Abs1) from the first measurement (step 5) should be multiplied by a Factor (F) in order to correct for changes in volume, as follows:

$$F = (\text{Sample vol} + \text{R1 vol}) / (\text{Sample vol} + \text{R1 vol} + \text{R2 vol})$$

For this assay: $F = (5+200) / (5+200+100) = 0.67$

Therefore: $\text{Sample}_{\text{abs}} = \text{Abs2} - (\text{Abs1} * 0.67)$

Specimen: Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

Assay Linearity: 4.0 mEq/L

Reagent Stability: 7 days at 2-8°C

Stability of Final Reaction: 60 minutes