

May 09, 2019

U Cinn - Phospholipids Assay

DOI

dx.doi.org/10.17504/protocols.io.xmvfk66



Patrick Tso¹, Dana Lee¹

¹University of Cincinnati

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xmvfk66>

External link: <https://mmpc.org/shared/document.aspx?id=198&docType=Protocol>

Protocol Citation: Patrick Tso, Dana Lee 2019. U Cinn - Phospholipids Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xmvfk66>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 31, 2019

Last Modified: May 09, 2019

Protocol Integer ID: 19861

Keywords: phospholipids, determinations of phospholipid, phospholipids assay summary, phospholipids in the sample, phospholipid, blue pigment, amount of phospholipid, absorbance of the blue color, enzymatic assay, assay, serum, absorbance,


Abstract

Summary:

Determinations of phospholipids in plasma/serum/lymph will be made using the Wako Phospholipids C enzymatic assay. In this assay, phospholipids in the sample are hydrolyzed ultimately producing a blue pigment. The amount of phospholipids in the sample is determined by measuring the absorbance of the blue color.

Materials

MATERIALS

 Phospholipid C assay kits standard included in kit **FUJIFILM Wako Pure Chemical Corporation Catalog #433-36201**

Working Reagent:

Reagents and Materials:

Color Reagent

Buffer

Procedure:


Reconstitute one vial of **Color Reagent** with a portion of **Buffer** then transferring entire contents to **Buffer** bottle, rinsing **Color** vial several times.

Note:

FUJIFILM Wako [RRID:SCR_013651](#)

Troubleshooting

- 1 Prepare working standards by making a serial dilution of the stock 300mg/dl standard.
(See diagram below) *Stock standard included in kit.



- 2 Prepare **Working Reagent** by reconstituting one vial of **Color Reagent** with a portion of **Buffer** then transferring entire contents to **Buffer** bottle, rinsing **Color Reagent** vial several times.
- 3 Using a 96 well flat bottom plate, into separate wells, pipette 2µL of deionized water, standard, or sample to be assayed.
- 4 Add 300µL of **Working Reagent** to all wells.
- 5 Incubate plate for 5 minutes at 37°C.
- 6 Determine the absorbance (abs) of the standards and of each unknown at 600nm.
- 7 Calculate values of unknowns from the standard curve.

Specimen: Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

Assay Linearity: 1000 mg/dl

Reagent Stability: 7 days at 2-8°C

Stability of Final Reaction: 60 minutes