

May 09, 2019

## U Cinn - Triglyceride Assay

DOI

[dx.doi.org/10.17504/protocols.io.xm4fk8w](https://dx.doi.org/10.17504/protocols.io.xm4fk8w)



Patrick Tso<sup>1</sup>, Dana Lee<sup>1</sup>

<sup>1</sup>University of Cincinnati

Mouse Metabolic Phenotyping Centers  
Tech. support email: [info@mmpc.org](mailto:info@mmpc.org)



Lili Liang

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.xm4fk8w](https://dx.doi.org/10.17504/protocols.io.xm4fk8w)

External link: <https://mmpc.org/shared/document.aspx?id=200&docType=Protocol>

**Protocol Citation:** Patrick Tso, Dana Lee 2019. U Cinn - Triglyceride Assay. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.xm4fk8w>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** January 31, 2019

**Last Modified:** May 09, 2019

**Protocol Integer ID:** 19868

**Keywords:** Triglyceride Assay,



## Abstract

### Summary:

Determinations of triglycerides in plasma/serum/lymph will be made using a Randox Triglycerides colorimetric kit. The triglycerides are determined after enzymatic hydrolysis with lipases.

## Materials

### MATERIALS

 Triglyceride Assay Kit **Randox Catalog #TR213**

### Working Reagent:

*Reagents and Materials:*

**Enzyme Reagent R1b**

**Buffer R1a**

*Procedure:*

Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.

### Note:

**Randox Life Sciences, RRID:SCR\_005525**

- 1 Prepare **Working Standards** by making a serial dilution of the stock 200mg/dl standard. (See diagram below) \*Stock standard is included in kit



- 2 Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.
- 3 Using a 96 well flat bottom plate, into separate wells, pipette 2µL of deionized water, standard, or sample to be assayed.
- 4 Add 200µL of reconstituted **Working Reagent** to all wells.
- 5 Incubate plate for 5 minutes at 37°C
- 6 Determine the absorbance (abs) of the standards and of each unknown at 500nm.
- 7 Calculate values of unknowns from the standard curve.

**Specimen:** Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

**Assay Linearity:** 1172 mg/dl

**Working Reagent Stability After Reconstitution:** 21 days stored at 2-8°C. PROTECT FROM LIGHT

**Stability of Final Reaction:** 60 minutes