

May 09, 2019

U Cinn - Triglyceride Assay

DOI

dx.doi.org/10.17504/protocols.io.xm4fk8w



Patrick Tso¹, Dana Lee¹

¹University of Cincinnati

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xm4fk8w>

External link: <https://mmpc.org/shared/document.aspx?id=200&docType=Protocol>

Protocol Citation: Patrick Tso, Dana Lee 2019. U Cinn - Triglyceride Assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xm4fk8w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 31, 2019

Last Modified: May 09, 2019

Protocol Integer ID: 19868

Keywords: Triglyceride Assay, determinations of triglyceride, randox triglycerides colorimetric kit, triglyceride assay summary, randox triglyceride, triglyceride, lipase, enzymatic hydrolysis with lipase, serum, assay, colorimetric kit,

Abstract

Summary:

Determinations of triglycerides in plasma/serum/lymph will be made using a Randox Triglycerides colorimetric kit. The triglycerides are determined after enzymatic hydrolysis with lipases.

Materials

MATERIALS

 Triglyceride Assay Kit **Randox Catalog #TR213**

Working Reagent:

Reagents and Materials:

Enzyme Reagent R1b

Buffer R1a

Procedure:

Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.

Note:

Randox Life Sciences, [RRID:SCR_005525](#)

Troubleshooting

- 1 Prepare **Working Standards** by making a serial dilution of the stock 200mg/dl standard.
(See diagram below) *Stock standard is included in kit



- 2 Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.
- 3 Using a 96 well flat bottom plate, into separate wells, pipette 2µL of deionized water, standard, or sample to be assayed.
- 4 Add 200µL of reconstituted **Working Reagent** to all wells.
- 5 Incubate plate for 5 minutes at 37°C
- 6 Determine the absorbance (abs) of the standards and of each unknown at 500nm.
- 7 Calculate values of unknowns from the standard curve.

Specimen: Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

Assay Linearity: 1172 mg/dl

Working Reagent Stability After Reconstitution: 21 days stored at 2-8°C. PROTECT FROM LIGHT

Stability of Final Reaction: 60 minutes