

May 10, 2019

U Mass - Glucose Tolerance Test with insulin secretion

DOI

dx.doi.org/10.17504/protocols.io.xxbfpin



Jason Kim¹

¹University of Massachusetts

Mouse Metabolic Phenotyping Centers
Tech. support email: info@mmpc.org



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.xxbfpin

External link: <https://mmpc.org/shared/document.aspx?id=142&docType=Protocol>

Protocol Citation: Jason Kim 2019. U Mass - Glucose Tolerance Test with insulin secretion. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xxbfpin>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2019

Last Modified: May 10, 2019

Protocol Integer ID: 20163

Keywords: Glucose tolerance test, insulin secretion, obesity, Mice



Abstract

Summary:

Glucose tolerance test with insulin secretion measures systemic clearance of glucose and systemic appearance of insulin following an intraperitoneal bolus injection of 20% dextrose. This experiment measures insulin sensitivity and insulin secretion (pancreatic β -cell function) in awake mice assuming that there are no alterations in systemic insulin clearance. Insulin sensitivity and pancreatic β -cell function are altered in obese mice.

Materials

MATERIALS

⊗ 20% Dextrose injection USP **Pfizer (Hospira) Catalog #NDC0409-7935-19**

⊗ Insulin Ultrasensitive ELIZA **Alpco Catalog #80-INSMSU-E01**

Note:

Hospira, RRID:SCR_003985

Insulin Ultrasensitive ELIZA, Cite this, **(Alpco Diagnostics Cat# 80-INSMSU-E01, RRID:AB_2792981)**

- 1 Mice may be fasted overnight (~15 hours) or for 5 hours prior to the start of experiment.
- 2 Collect plasma sample (20 μ l) before the start of experiment (basal-0 min) to measure basal glucose and insulin levels.
- 3 Administer intraperitoneal injection of 20% dextrose (1 or 2 g/kg body weight) using an insulin syringe.
- 4 Collect plasma samples (20 μ l) at 10, 20, 30, 60, 90, and 120 min following injection to measure circulating glucose and insulin concentrations.
- 5 For data analysis, plasma glucose levels vs. time after injection are plotted, and area-under-curve may be calculated to estimate insulin sensitivity. Also, plasma insulin levels vs. time after injection are plotted, and area-under-curve may be calculated to estimate insulin secretion (pancreatic β -cell function).